# Explorer™ Balances

## User Guide

Scan the QR code to download  
the EU Declaration of Conformity

English

Español

Français

Deutsch

Italiano

Português

Svenska

Nederlands

Dansk

Polski

Čeština

Magyar

## 1. INTRODUCTION

Thank you for choosing our product. This user guide is designed to assist you in the installation and safe operation of your balance. We strongly recommend that you read this guide thoroughly before proceeding with the installation and use of the device. For comprehensive instructions and additional information, please visit our website at [www.ohaus.com/products](http://www.ohaus.com/products).

## 2. ASSEMBLY

### 2.1 For 0.01mg, 0.1mg, and 1mg Models:

Unpack the balance and install the weighing pan

### 2.2 For 0.01g, and 0.1g Models:

Unpack the balance and install the weighing pan

## 3. INSTALLATION

### 3.1 Selecting a Location

Avoid excessive vibrations, heat sources, air current, or rapid temperature changes. Allow sufficient space.

### 3.2 Connecting Power

Connect the AC adapter or the power cord to the instrument's power input connector, and then connect the AC plug to a suitable electrical outlet.

**CAUTION: FOR USE WITH CSA CERTIFIED (OR EQUIVALENT APPROVED) POWER SOURCE, WHICH MUST HAVE A LIMITED CURRENT OUTPUT.**

Note: Allow equipment to warm up for 60 minutes for optimal weighing performance.

### 3.3 Turn on the Balance

Short press the power button to turn on the balance.

### 3.4 Leveling the Balance

Navigate to QUICK SETUP > DIGITAL LEVELING  
Adjust the balance feet according to arrow directions

### 3.5 Initial Calibration

To achieve optimal performance, it is recommended to verify the calibration of the balance before using it for the first time. Navigate to Main Menu > CALIBRATION > INTERNAL CALIBRATION to calibrate the balance.

### 3.6 Connecting the Interface

Enhanced communication options include 2 USB Host (Type A) x 2, USB Device (Type C) x 1, Ethernet Lan port (RJ45) x 1, RS232 x 1, and Optional Wi-Fi, Bluetooth Dongle.


## 4. OPERATION

The Explorer balance features a 7-inch display with vibrant colors and a glass panel that is easy to clean and protects against chemical spills.

### 4.1 Overview of Display

### 4.2 Main Screen

### 4.3 Using the Balance

Before using any application, be sure the balance has been leveled and calibrated.

- If required, place an empty container on the pan and press Tare.
- Add sample to the pan or container. The display shows the weight of the sample.

## 4.4 Weighing Application Modes

Basic Weighing	The Basic Weighing application is used to determine the weight of items in the selected unit of measure.
Parts Counting	Parts Counting displays the number of pieces or parts based on an average piece weight.
Check Counting	Checks if the current sample pieces are within tolerances (e.g. an over and under limit).
Percent Weighing	Percent Weighing displays the current weight as a percentage of a reference weight.
Check Weighing	Check Weighing checks if the current weight is within tolerances (e.g. an over and under limit).
Dynamic Weighing	Dynamic Weighing is used to weigh objects that are not stable, such as animals.
Totalization	Totalization is used to sum the weights of multiple samples and report the statistical data for the series of samples.
Formulation	Formulation is used to combine various elements in proportionate amounts.
Differential	Differential is used to calculate the difference in weights of multiple samples taken at different times.
Density Determination	Density Determination is used to determine the density of a solid or a liquid.
Peak Hold	Peak Hold captures the maximum weight in a series of weighing.

## 5. MAINTENANCE

### 5.1 Calibration

Always use calibrated and certified weights for checking the balance. These weights should be traceable to national or international standards. Allow the balance to warm up for the recommended period (usually 60 minutes) before performing the calibration check. This helps stabilize the balance and ensures accurate readings.

### 5.2 Cleaning

Clean the balance thoroughly before calibration to remove any dust or debris that might affect the readings.

**WARNING:** Electric Shock Hazard. Disconnect the equipment from the power supply before cleaning. Make sure that no liquid enters the interior of the instrument.

**ATTENTION:** Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents. The housing may be cleaned with a cloth dampened with a mild detergent if necessary.

### 5.3 Service Information

If the troubleshooting section in the user manual doesn't solve your problem, contact an Authorized Ohaus Service Agent. Visit [www.ohaus.com](http://www.ohaus.com) to locate the nearest Ohaus office. An Ohaus Product Service Specialist will be ready to help you.

## 6. TECHNICAL DATA

### 6.1 Specifications

#### Ambient Conditions

- Indoor use only
- Altitude: 2000 m
- Specified Temperature range: 10°C to 30°C
  - EXR 0.1mg NTEP models: 15°C to 25°C
  - Operability is assured at ambient temperatures between 5°C and 40°C
- Humidity: Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 40°C.
- Electrical supply: 12VDC, 1.5A. (For use with certified or approved power supply, which must have a SELV and limited energy output)
- Mains supply voltage fluctuations: up to ±10% of the nominal voltage
- Installation category: II
- Pollution degree: 2

#### Materials

- Bottom Housing: Die-cast Aluminum, Painted
- Top Housing: Die-cast Aluminum, Painted
- Terminal: Glass, Die-cast Aluminum
- Weighing Platforms: Plastic (PC)
- Weighing Pan:
  - Zink Alloy (0.01mg models)
  - 316 SST (0.1mg, 1mg, 0.01g, 0.1g models)
- In-use Cover: Plastic (PET)
- Draft Shield; Glass, Aluminum, Plastic

## 7. SAFETY PRECAUTIONS

### 7.1 Intended Use

This instrument is intended for use in laboratories, pharmacies, schools, businesses and light industry. It must only be used for measuring the parameters described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from OHAUS, is considered as not intended. This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the instrument is not used according to these operating instructions, the intended protection provided by the instrument may be impaired.

#### Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

WARNING	For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.
CAUTION	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
ATTENTION	For important information about the product. May lead to equipment damage if not avoided.
NOTE	For useful information about the product.

#### Warning Symbols

General hazard

Explosion hazard

Electrical shock hazard

### 7.2 Safety Notes

- CAUTION: Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.
- Before connecting power, verify that the AC adapter's input voltage range and plug type are compatible with the local AC mains power supply.
- Do not position the equipment such that it is difficult to reach the power connection.
- Make sure that the power cord does not pose a potential obstacle or tripping hazard.
- Operate the equipment only under ambient conditions specified in these instructions.
- This equipment is for indoor use only.
- Do not operate the equipment in wet, hazardous or unstable environments.
- Do not allow liquids to enter the equipment.
- Do not load the equipment above its rated capacity.
- Do not drop loads on the platform.
- Do not place the equipment upside down on the platform.
- Use only approved accessories and peripherals.
- Disconnect the equipment from the power supply when cleaning.
- Service should only be performed by authorized personnel.

**WARNING:** Never work in an environment subject to explosion hazards! The housing of the instrument is not gas tight. (Explosion hazard due to spark formation, corrosion caused by the ingress of gases).

**WARNING:** Electrical shock hazards exist within the housing! The housing should only be opened by authorized and qualified personnel. Remove all power connections to the unit before opening.

**CAUTION:** The protection provided by the unit may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

**WARNING: ELECTRICAL SHOCK HAZARD**

Avoid touching two carbon brushes while the Ionizer is turned on.

## 8. COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC), 2014/35/EU (LVD) and 2014/31/EU (NAWI). The EU Declaration of Conformity is available online at <a href="http://www.ohaus.com/ce">www.ohaus.com/ce</a> .
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to <a href="http://www.ohaus.com/weee">www.ohaus.com/weee</a> .
	EN 61326-1
	CAN/CSA-C22.2 No. 61010-1 UL 61010-1

### Industry Canada Note

CAN ICES-003(A) / NMB-003(A)

ISO 9001 Registration

ISO 14001 Registration

The management system governing the production of this product is ISO 9001 certified.

## FCC Supplier Declaration of Conformity

Unintentional Radiator per 47CFR Part B  
Trade Name: OHAUS CORPORATION  
Model: Explorer™ EXR...

Party issuing Supplier's Declaration of Conformity:

Ohaus Instruments (Changzhou) Co., Ltd.  
Building C, No. 6 Zhengqiang Road, Xuejia Town, Xinbei District, Changzhou  
Jiangsu 213022,  
China  
Phone: +86 519 85287270

Responsible Party – U.S. Contact Information:

Ohaus Corporation  
8 Campus Drive, Suite 105  
Parsippany, NJ 07054  
United States  
Phone: +1 973 377 9000  
Web: [www.ohaus.com](http://www.ohaus.com)

FCC Compliance Statement:

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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P/N 31064681 C © 2025 Ohaus Corporation, all rights reserved.