

Elmasonic Accessories - Baskets

made of stainless steel and with Rilsan-coated handles



	10	15/20	30	40	60	80	100	120	100/120 Einlage
Order number	1004170	1004209	1004176	1004232	1004243	1004265	1004178	1004277	1076005
List price	36 €	40 €	46 €	52 €	65 €	120 €	140 €	140 €	62 €
Basket int. dim. W/D/H (mm)	175/70/30	110/100/50	195/105/50	190/105/75	255/115/75	455/106/75	255/200/70	250/190/110	285/195/30
Mesh size YxZ (mm)	6x1				9x1				
Loading weight max. (kg)	1,0			3,0	5,0		6,0	7,0	6,0
Suitable for									
• Elmasonic Easy	•	•	•	•	•		•	•	•
• Elmasonic Med			•		•		•	•	•
• Elmasonic Select			•	•	•	•	•	•	•
• Elmasonic S	•	•	•	•	•	•	•	•	•
• Elmasonic P			•		•				
• Elmadry TD			•				•	•	•

All prices ex works, excluding packaging, plus VAT.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.

Elmasonic Accessories - Baskets

- made of stainless steel (W Nr. 1.4301) and with Rilsan-coated handles
- a basket insert (order no. 1076005) is available for baskets 100 and 120 (Denta)
- the basket 200 (order no. 1108458) also serves as a holder for sterilisation baskets
 - DIN 1/1 short (W/D - 485/253 mm)
 - DIN 3/4 (W/D - 406/253 mm)



	100/120 Denta		150	180	200	300	500	900
Order number	1077166	1071801	1031243	1004275	1108458	1004271	1107675	1107641
List price	155 €	175 €	180 €	160 €	180 €	180 €	285 €	375 €
Basket int. dim. W/D/H (mm)	255/200/25	255/200/45	455/250/45	280/250/115	490/255/95	455/250/110	460/300/195	545/460/190
Mesh size YxZ (mm)	9x1				11	7x1	11	
Loading weight max. (kg)	6,0		7,0	8,0	4,0	10,0	20,0	30,0
Suitable for								
• Elmasonic Easy	•	•		•		•		
• Elmasonic Med	•	•		•	•	•	•	•
• Elmasonic Select	•	•	•	•		•	•	•
• Elmasonic S	•	•	•	•		•	•	•
• Elmasonic P	•	•		•		•		
• Elmadry TD	•	•				•		

All prices ex works, excluding packaging, plus VAT.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.


Elmasonic Accessories - Baskets

Dimension definition

Basket internal dimensions

Mesh size

Y = Mesh size; Z = Wire diameter



The drawings are merely included as a visualization aid and are not true to scale.

Technical specifications serve as orientation values and are subject to tolerances due to production and manufacturing. Proper operation can only be guaranteed when using original Elma devices.