

KERN OZO556



Individuality, variety and flexible working through our modular construction system â Stereo microscope heads



Category

Brand: Optics

Packaging & shipping

Delivery: 24 h
Dimensions packaging (WxDxH): 310 x 300 x 220 mm
Gross weight: 3 kg
Net weight: 2,451 kg
Shipping method: Parcel service

Construction

Eye distance (max): 76
Eye distance (min): 52
Optical system: Greenough
Rotation of tube 360°: yes
Tube: Binocular
Tube angle: 35°