



- Large range available, treated for cell culture
- Optically clear and uniform surface suitable for microscopy
- 60 and 40 mm dishes available with grids for cloning or determination of plating efficiency
- Nunclon™ Δ certified

**LITERATURE:**

Nunc Bulletin No: 2, 3

Page 293

**Dishes Nunclon™ Δ**  
Polystyrene. With lid. Sterile




| Cat. No.                          | 150318  | 153066  | 174926  | 171099* | 174888** | 150326  | 150288  | 169558  |
|-----------------------------------|---------|---------|---------|---------|----------|---------|---------|---------|
| Name                              | 35 x10  | 35 x 10 | 35 x 10 | 35 x 10 | 60 x 15  | 60 x 15 | 60 x 15 | 60x15   |
| Grid size, mm                     | -       | -       | 2 x 2   | -       | -        | -       | -       | 2 x 2   |
| Culture area, cm <sup>2</sup>     | 8.8     | 8.8     | 8.8     | 8.8     | 20.8     | 21.5    | 21.5    | 21.5    |
| Airvent                           | -       | +       | +       | +       | +        | -       | +       | +       |
| Max. External dimensions, mm      | 40 x 12 | 40 x 12 | 40 x 12 | 40 x 12 | 60 x 15  | 60 x 15 | 60 x 15 | 60 x 15 |
| Max. External diameter bottom, mm | 34.7    | 34.7    | 34.7    | 34.7    | 48       | 53.5    | 53.5    | 53,5    |
| Suggested working volume, ml      | 3       | 3       | 3       | 3       | 5        | 5       | 5       | 5       |
| Units per pack/case               | 10/500  | 10/500  | 20/500  | 20/500  | 20/500   | 10/400  | 10/400  | 10/400  |

\* For suspension culture - not Nunclon Delta treated

\*\*Made from Permanox™

**Dishes Nunclon™ Δ**  
Polystyrene. With lid. Sterile



| Cat. No.                          | 150350*  | 172958*  | 168381*  | 157150   | 166508         |
|-----------------------------------|----------|----------|----------|----------|----------------|
| Name                              | 100 x 15 | 100 x 20 | 150 x 20 | 150 x 20 | -              |
| Grid size, mm                     | -        | -        | -        | -        | -              |
| Culture area, cm <sup>2</sup>     | 56.7     | 56.7     | 145      | 150      | 500            |
| Airvent                           | +        | +        | +        | +        | -              |
| Max. External dimensions, mm      | 92 x 17  | 92 x 21  | 144 x 21 | 147 x 18 | 245 x 245 x 25 |
| Max. External diameter bottom, mm | 87.2     | 87.1     | 138.3    | 136      | -              |
| Suggested working volume, ml      | 12.5     | 12.5     | 35       | 35       | 135            |
| Units per pack/case               | 10/150   | 10/480   | 10/80    | 10/120   | 4/16           |

\* Packed in resealable bags