# pipetman®

4

## Air-Displacement Pipette STANDARD PIPETTING WITH REDUCED FORCES



#### Minimized Pipetting Forces

PIPETMAN<sup>®</sup> G offers one of the lowest pipetting forces commercially available. With high quality piston assembly system including superior piston grease, PIPETMAN G enables hours of comfortable pipetting.

#### Reliable Results

PIPETMAN G relies on the legendary quality, accuracy, and precision of PIPETMAN, with reduced pipetting forces.

### Built-to-Last and Convenient

Made of PVDF and stainless steel, PIPETMAN G offers premium durability. The PIPETMAN G clip ejector allows for convenient access to serviceable parts.

Durable and reliable, they are available either in single or multichannel models equipped with metal tip ejectors.





-

		PIPET	MAN® G Si	ingle Chann	el
	Model	PIPETMA DIAMON		Part Number	Volume
5 <sup>2-2</sup> 4 P2	P2G	D10 DL10	DF10ST DFL10ST	F144054M	0.2-2 μL
(PIO)	P10G	D10 DL10	DF10ST DFL10ST	F144055M	1-10 μL
(* <sup>20</sup> 4) P20	P20G	D200	DF30ST	F144056M	2-20 μL
51-100-4 P100	P100G	D200	DF100ST	F144057M	10-100 μL
\$ <sup>250</sup> 2+ P200	P200G	D200 D300	DF200ST DF300ST	F144058M	20-200 μL
PIDO	P1000G	D1000 D1200	DF1000ST DF1200ST	F144059M	100-1000 μL
P300	P5000G*	D5000		F144066	500-5000 μL
Piom.	P10mLG*	D10mL		F144067	1000-10000 μL
		PIPE	TMAN® G I	Multichanne	
(FB)	P8x10 P12x10	D10 DL10	DF10ST DFL10ST	F144068 F144069	1-10 μL
N-2044 P20	P8x20 P12x20	DL10 D200	DFL10ST DF30ST	F144070 F144071	2-20 μL
\$ <sup>200</sup> 4	P8x200 P12x200	D200 D300	DF200ST DF300ST	F144072 F144073	20-200 μL
9-500-4- P300	P8x300 P12x300	D200 D300	DF200ST DF300ST	F144074 F144075	30-300 μL

Models equipped with a tip ejector are only sold with a stainless steel tip ejector. Plastic tip ejector is available as an accessory. Refer to *PIPETMAN G User's Guide LT801122*. \*P5000G and P10mLG do not have tip ejectors.

Gilson maximum permissible errors are guaranteed only when PIPETMAN pipettes are used with the recommended PIPETMAN DIAMOND Tips.

#### Accessories

COMFORT HANDLE	
Model	Part Number
COMFORT HANDLE pack of 3, blue	F161486
COMFORT HANDLE pack of 3, white	F161488
COMFORT HANDLE pack of 50, blue	F161485
COMFORT HANDLE pack of 50, white	F161487
COMFORT HANDLE stickers, 120	F161499

DISPOSABLE FILTERS	
Model	Part Number
Bag of 10 Filters for P5000/P10mL	F161280
Bag of 100 Filters for P5000/P10mL	F161241

PLASTIC EJECTORS		
Model	Part Number	
P2/P10 with adapter	F107027	
P20	F107028	
P100	F107029	
P200	F107030	
P1000	F107031	

PIPETMAN G specifications are downloadable from <u>www.gilson.com</u>

## **4 GILSON**<sup>®</sup>

#### www.gilson.com/contactus

LT800591/C  $\mid$  @2022 Gilson S.A.S. All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.