







Compatibility Chart AHN myPlate® qPCR 96 Well Plate		AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
		Image						
		Standard Well (0.2 mL, height 21 mm)	■	■	■			
		Low Profile Well0 (0.1 mL, height 16 mm)				■	■	■

		qPCR 96 Well Plate						
Applied Biosystems (Life Technology)		AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
Device	Note / Holder No.							
7000	4350105			■				
7300	4367032			■				
7500	4367032			■				
7500Fast	4359652					■		
7900HT Fast Model(0.1 ml)						■		
7900HT Standard (0.2 ml)				■				
StepOne Plus						■		
ViiA7 Fast Model (0.1 ml)	4459846					■		
ViiA7 (0.2 ml)	4459845			■				
QuantStudio Fast Model (0.1 ml)	4459846					■		
QuantStudio (0.2 ml)	4459845			■				
QuantStudio 3/5 Fast (0.1 ml)						■		
QuantStudio 3/5 (0.2 ml)				■				
5700/7700	Sequencer			■	■			
3100/3130 (xL) (0.1 ml Fast 96 Well)	4367470 (Retainer)4315933 (Septa)					■		
310/3100/3130 (xL) (0.2 ml Standard 96 Well)	4317241 (Retainer)4315933 (Septa)			■	■			
3700/3730 (xL) (0.1 ml Fast 96 Well)	4367469 (Retainer)4315933 (Septa)					■		
3700/3730 (xL) (0.2 ml Standard 96 Well)	4410228 (Retainer)4315933 (Septa)			■	■			
3500(xL) (0.1 ml Fast 96 Well)	4409530 (Retainer)4412614 (Septa)					■		
3500(xL) (0.2 ml Standard 96 Well)	4410228 (Retainer)4412614 (Septa)			■	■			

		qPCR 96 Well Plate						
Roche		AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
Device	Note							
LightCycler 480 & 480II								■
LightCycler 96								■

Compatibility Chart AHN myPlate® qPCR 96 Well Plate		AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
		Image						
		Standard Well (0.2 mL, height 21 mm)	■	■	■			
		Low Profile Well0 (0.1 mL, height 16 mm)				■	■	■

		qPCR 96 Well Plate						
Bio-Rad		AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
Device	Note							
iCycler iQ			■		■			
iQ5			■		■			
MyiQ			■		■			
Chromo4			■		■			
Opticon I & II							■	■
CFX96							■	■

		qPCR 96 Well Plate						
Stratagene/Eppendorf/Biometra		AHN Cat. No.	3-047-01-0	3-047-03-0	3-047-02-0	3-046-02-0	3-046-01-0	3-046-03-9
Device	Note							
Mx4000/3000P/3005P			■		■			
AriaMx							■	■
MasterCycle ep Realplex			■		■			
Optical Thermocycler			■		■			