

PCR Tube Racks 96-Well

Perfect for pre- or post-PCR sampling

- Made of polypropylene
- Accommodates individual tubes or strips
- Removable hinged lid lies flat
- Alphanumeric identification for quick location and retrieval of tubes



Pack of 5

| Item No. | Color | Item No. | Color |
|----------|-----------|----------|-------|
| HS23461A | Assorted* | HS23461C | Blue |
| HS23461B | Natural | HS23461D | Green |

*Assorted colors (Blue, Green, Pink, Yellow, Orange)

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|------------------|
| 0.2 | 96 | 5.2 | 5 x 3.9 x 1.2 | 12.7 x 9.9 x 3.0 |



PCR Tube Racks Reversible

Removable hinged lid fits both sides

- Made of polypropylene
- PCR side of the rack holds 8 or 12 x 0.2 mL tube strips or individual tubes
- Reverse side holds 0.5 mL tubes
- Both sides of the rack have 12 wells for 1.5 mL tubes

Pack of 5

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in (cm) |
|----------|-------|-------------|-------|-----------------|----------------------|
| HS2344A | | 0.2 | 168 | 5.5 | 8.1 x 5.2 x 2.1 |
| | | 0.5 | 40 | 7.6 | (20.6 x 13.1 x 5.3) |
| | | 1.5/2.0 | 12 | 11.4 | |

*Assorted colors (Blue, Green, Pink, Yellow, Orange)



Engineering, OEM & Industrial Design

- We have in-house product engineering and R&D teams that are dedicated to delivering quality, innovative, efficient, and compact solutions that customers can rely on.
- We have a range of capabilities from private labeling to product development, design, packaging and manufacturing.

