


Inventory Systems

for Freezers and Refrigerators

2011



Type	Boxes	Grid Dividers	Vertical Racks	Horizontal Racks	Drawer Racks	Upright Bins
Polypropylene						
EPPI® Box 130 x 130 / 133 x 133	3–26	27	29	30–31	32–37	38
EPPI® Box 67 x 67	3		★	★	★	38
Polycarbonate						
Kryo Box 130 x 130 / 133 x 133	28		29	30–31	32–37	38
Fiberboard						
ALPHA Box 136 x 136	40	40	41	42–43	44–48	49
ALPHA Box long2 136 x 268	50	50	★	★	★	49
ALPHA Box short 136 x 94	51	51	★	★	★	49
BRAVO Box 133 x 133	52	52	53	54 + 30	55–57 32–33	38
BRAVO Box long2 133 x 257	58	58	★	★	★	38
BRAVO Box short 133 x 92	59	59	★	★	★	38
CellBox 122 x 122 / 122 x 237 148 x 148 / 148 x 287	60	60	61	61	★	38 + 49
DELTA Box 133 x 133	62	62	53	54 + 30	55–57 32–33	38
MIKE Box 142 x 142	63	63	★	★	★	38 + 49
NANU Box 76 x 76	64	64	64	64	★	38 + 49
Slide Boxes						
Obi Box	39		★	★	★	38 + 49
Mikrotiter Plates and DeepWells						
All brands 86 x 128			65	66–67	68–69	38 + 49
Marking System						
EPPI® Ident	11 + 39					

★ = on request

Polypropylene: -100 up to +121°C · Polycarbonate: -196 up to +121°C




EPi® 32 / 12x12 • height 32 mm fix.

12x12 conical holes for 144 PCR tubes 0.2 ml.
Individual or strips. With alphanumeric ID code.
Removable safety baseplate.

Cat. No.	Colour	Cat. No.	Colour
KEB32-CS144NA	natural	KEB32-CS144YE	yellow
KEB32-CS144NY	neon-yellow	KEB32-CS144BU	blue
KEB32-CS144NB	neon-blue	KEB32-CS144GR	green
KEB32-CS144NG	neon-green	KEB32-CS144RE	red
KEB32-CS144NP	neon-red/pink	KEB32-CS144VT	violet
KEB32-CS144NO	neon-orange	KEB32-CS144BL	black
KEB32-CS144WH	white	KEB32-CS144BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 37 / 10x10 • height 37 mm fix.

Inner height 35 mm. 10x10 divider holds for vials up to 11 mm diameter.

Cat. No.	Colour	Cat. No.	Colour
KEA37-C100NA	natural	KEA37-C100YE	yellow
KEA37-C100NY	neon-yellow	KEA37-C100BU	blue
KEA37-C100NB	neon-blue	KEA37-C100GR	green
KEA37-C100NG	neon-green	KEA37-C100RE	red
KEA37-C100NP	neon-red/pink	KEA37-C100VT	violet
KEA37-C100NO	neon-orange	KEA37-C100BL	black
KEA37-C100WH	white	KEA37-C100BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 37 / 1x1 • height 37 mm fix.

Inner height 35 mm. Without grid dividers. Useful for EPi® holes inserts or grid dividers (page 27).

Cat. No.	Colour	Cat. No.	Colour
KEA37-C1NA	natural	KEA37-C1YE	yellow
KEA37-C1NY	neon-yellow	KEA37-C1BU	blue
KEA37-C1NB	neon-blue	KEA37-C1GR	green
KEA37-C1NG	neon-green	KEA37-C1RE	red
KEA37-C1NP	neon-red/pink	KEA37-C1VT	violet
KEA37-C1NO	neon-orange	KEA37-C1BL	black
KEA37-C1WH	white	KEA37-C1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 45 Junior / 4x4 • height 45-60 mm variable.

16-place divider holds vials up to 14 mm diameter. Removable grid divider.
External dimensions 67x67 mm. Other colours on request.

Cat. No.	Grid	Colour
KEA45-VJ16NA	4x4	natural
KEA45-VJ16NB	4x4	neon-blue
KEA45-VJ16NG	4x4	neon-green
KEA45-VJ16NO	4x4	neon-orange
KEA45-VJ16NP	4x4	neon-red/pink
KEA45-VJ16NY	4x4	neon-yellow
KEA45-VJ16BB	4x4	black-black*
KEA45-VJ1NA	without	natural

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 34

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 34

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 34

Matching Racks


Vertical Racks
Please request!

Horizontal Racks
Please request!

Drawer Racks
Please request!

EPi® Box made of Polypropylene

130 x 130 mm




EPi® 45 / 9x9 • height 45-53 mm variable.

9x9 divider holds 81 vials up to 12.5 mm diameter. Without ID code. (See also page 6: KEA50)

Cat. No.	Colour
KEA45-V81NA	natural
KEA45-V81NY	neon-yellow
KEA45-V81NB	neon-blue
KEA45-V81NG	neon-green
KEA45-V81NP	neon-red/pink
KEA45-V81NO	neon-orange
KEA45-V81WH	white

Cat. No.	Colour
KEA45-V81YE	yellow
KEA45-V81BU	blue
KEA45-V81GR	green
KEA45-V81RE	red
KEA45-V81VT	violet
KEA45-V81BL	black
KEA45-V81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 45 / 9x9 • height 45-53 mm variable.

9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides. (See also page 6: KEB50)

Cat. No.	Colour
KEB45-V81NA	natural
KEB45-V81NY	neon-yellow
KEB45-V81NB	neon-blue
KEB45-V81NG	neon-green
KEB45-V81NP	neon-red/pink
KEB45-V81NO	neon-orange
KEB45-V81WH	white

Cat. No.	Colour
KEB45-V81YE	yellow
KEB45-V81BU	blue
KEB45-V81GR	green
KEB45-V81RE	red
KEB45-V81VT	violet
KEB45-V81BL	black
KEB45-V81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 45 / 9x9 • height 45-53 mm variable.

9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the bottom of the box (transparent). (See also page 6: KEC50)

Cat. No.	Colour
KEC45-V81NA	natural
KEC45-V81NY	neon-yellow
KEC45-V81NB	neon-blue
KEC45-V81NG	neon-green
KEC45-V81NP	neon-red/pink
KEC45-V81NO	neon-orange



EPi® 45 / 9x9 • height 45-53 mm variable.

9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the bottom of the box (transparent) and on the lid. (See also page 6: KEE50)

Cat. No.	Colour
KEE45-V81NA	natural
KEE45-V81NY	neon-yellow
KEE45-V81NB	neon-blue
KEE45-V81NG	neon-green
KEE45-V81NP	neon-red/pink
KEE45-V81NO	neon-orange

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

EPi® Box made of Polypropylene

130 x 130 mm



EPi® 45 / 9x9 • height 45-53 mm variable.

9x9 divider holds 81 vials up to 12.5 mm diameter.
With alphanumeric ID code on the lid.
(See also page 7: KEH50)

Cat. No.	Colour
KEH45-V81NA	natural
KEH45-V81NY	neon-yellow
KEH45-V81NB	neon-blue
KEH45-V81NG	neon-green
KEH45-V81NP	neon-red/pink
KEH45-V81NO	neon-orange
KEH45-V81WH	white

Cat. No.	Colour
KEH45-V81YE	yellow
KEH45-V81BU	blue
KEH45-V81GR	green
KEH45-V81RE	red
KEH45-V81VT	violet
KEH45-V81BL	black



EPi® 45 / 9x9 • height 45-53 mm variable.

9x9 divider holds 81 vials up to 12.5 mm diameter.
With alphanumeric ID code on all 4 sides as well as on the lid. (See also page 7: KEK50)

Cat. No.	Colour
KEK45-V81NA	natural
KEK45-V81NY	neon-yellow
KEK45-V81NB	neon-blue
KEK45-V81NG	neon-green
KEK45-V81NP	neon-red/pink
KEK45-V81NO	neon-orange
KEK45-V81WH	white

Cat. No.	Colour
KEK45-V81YE	yellow
KEK45-V81BU	blue
KEK45-V81GR	green
KEK45-V81RE	red
KEK45-V81VT	violet
KEK45-V81BL	black



EPi® 45 / 9x9 • height 45-53 mm variable.

9x9 divider holds 81 vials up to 12.5 mm diameter.
With numeric ID code (1-81) on the lid.
(See also page 7: KEM50)

Cat. No.	Colour
KEM45-V81NA	natural
KEM45-V81NY	neon-yellow
KEM45-V81NB	neon-blue
KEM45-V81NG	neon-green
KEM45-V81NP	neon-red/pink
KEM45-V81NO	neon-orange
KEM45-V81WH	white

Cat. No.	Colour
KEM45-V81YE	yellow
KEM45-V81BU	blue
KEM45-V81GR	green
KEM45-V81RE	red
KEM45-V81VT	violet
KEM45-V81BL	black



EPi® 45 / 9x9 • height 45-53 mm variable.

9x9 divider holds 81 vials up to 12.5 mm diameter.
With numeric ID code (1-81) on the bottom of the box (transparent) and on the lid. (See also page 7: KER50)

Cat. No.	Colour
KER45-V81NA	natural
KER45-V81NY	neon-yellow
KER45-V81NB	neon-blue
KER45-V81NG	neon-green
KER45-V81NP	neon-red/pink
KER45-V81NO	neon-orange

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks


Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

EPi® Box made of Polypropylene

130 x 130 mm




EPi® 50 / 9x9 • height 52 mm fix.

Inner height 50 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. Without ID code.

(See also page 4: KEA45)

Cat. No.	Colour	Cat. No.	Colour
KEA50-C81NA	natural	KEA50-C81YE	yellow
KEA50-C81NY	neon-yellow	KEA50-C81BU	blue
KEA50-C81NB	neon-blue	KEA50-C81GR	green
KEA50-C81NG	neon-green	KEA50-C81RE	red
KEA50-C81NP	neon-red/pink	KEA50-C81VT	violet
KEA50-C81NO	neon-orange	KEA50-C81BL	black
KEA50-C81WH	white	KEA50-C81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 50 / 9x9 • height 52 mm fix.

Inner height 50 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides. (See also page 4: KEB45)

Cat. No.	Colour	Cat. No.	Colour
KEB50-C81NA	natural	KEB50-C81YE	yellow
KEB50-C81NY	neon-yellow	KEB50-C81BU	blue
KEB50-C81NB	neon-blue	KEB50-C81GR	green
KEB50-C81NG	neon-green	KEB50-C81RE	red
KEB50-C81NP	neon-red/pink	KEB50-C81VT	violet
KEB50-C81NO	neon-orange	KEB50-C81BL	black
KEB50-C81WH	white	KEB50-C81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 50 / 9x9 • height 52 mm fix.

Inner height 50 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the bottom of the box (transparent). (See also page 4: KEC45)

Cat. No.	Colour
KEC50-C81NA	natural
KEC50-C81NY	neon-yellow
KEC50-C81NB	neon-blue
KEC50-C81NG	neon-green
KEC50-C81NP	neon-red/pink
KEC50-C81NO	neon-orange



EPi® 50 / 9x9 • height 52 mm fix.

Inner height 50 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the bottom of the box (transparent) and on the lid. (See also page 4: KEE45)

Cat. No.	Colour
KEE50-C81NA	natural
KEE50-C81NY	neon-yellow
KEE50-C81NB	neon-blue
KEE50-C81NG	neon-green
KEE50-C81NP	neon-red/pink
KEE50-C81NO	neon-orange

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

EPi® Box made of Polypropylene

130 x 130 mm




EPi® 50 / 9x9 • height 52 mm fix.

Inner height 50 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on the lid. (See also page 5: KEH45)

Cat. No.	Colour
KEH50-C81NA	natural
KEH50-C81NY	neon-yellow
KEH50-C81NB	neon-blue
KEH50-C81NG	neon-green
KEH50-C81NP	neon-red/pink
KEH50-C81NO	neon-orange
KEH50-C81WH	white

Cat. No.	Colour
KEH50-C81YE	yellow
KEH50-C81BU	blue
KEH50-C81GR	green
KEH50-C81RE	red
KEH50-C81VT	violet
KEH50-C81BL	black



EPi® 50 / 9x9 • height 52 mm fix.

Inner height 50 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the lid. (See also page 5: KEK45)

Cat. No.	Colour
KEK50-C81NA	natural
KEK50-C81NY	neon-yellow
KEK50-C81NB	neon-blue
KEK50-C81NG	neon-green
KEK50-C81NP	neon-red/pink
KEK50-C81NO	neon-orange
KEK50-C81WH	white

Cat. No.	Colour
KEK50-C81YE	yellow
KEK50-C81BU	blue
KEK50-C81GR	green
KEK50-C81RE	red
KEK50-C81VT	violet
KEK50-C81BL	black



EPi® 50 / 9x9 • height 52 mm fix.

Inner height 50 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With numeric ID code (1-81) on the lid. (See also page 5: KEM45)

Cat. No.	Colour
KEM50-C81NA	natural
KEM50-C81NY	neon-yellow
KEM50-C81NB	neon-blue
KEM50-C81NG	neon-green
KEM50-C81NP	neon-red/pink
KEM50-C81NO	neon-orange
KEM50-C81WH	white

Cat. No.	Colour
KEM50-C81YE	yellow
KEM50-C81BU	blue
KEM50-C81GR	green
KEM50-C81RE	red
KEM50-C81VT	violet
KEM50-C81BL	black



EPi® 50 / 9x9 • height 52 mm fix.

Inner height 50 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With numeric ID code (1-81) on the bottom of the box (transparent) and on the lid. (See also page 5: KER45)

Cat. No.	Colour
KER50-C81NA	natural
KER50-C81NY	neon-yellow
KER50-C81NB	neon-blue
KER50-C81NG	neon-green
KER50-C81NP	neon-red/pink
KER50-C81NO	neon-orange

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

EPi® Box made of Polypropylene

130 x 130 mm




EPi® 45 / 8x8 • height 45-53 mm variable.

8x8 holes for 64 reaction vessels 1.5 ml and 2.0 ml Safe-Lock® up to 11 mm diameter. With alphanumeric ID code. (See also page 9: KEB50)

Cat. No.	Colour
KEB45-V064NA	natural
KEB45-V064NY	neon-yellow
KEB45-V064NB	neon-blue
KEB45-V064NG	neon-green
KEB45-V064NP	neon-red/pink
KEB45-V064NO	neon-orange
KEB45-V064WH	white

Cat. No.	Colour
KEB45-V064YE	yellow
KEB45-V064BU	blue
KEB45-V064GR	green
KEB45-V064RE	red
KEB45-V064VT	violet
KEB45-V064BL	black
KEB45-V064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 45 / 7x7 • height 45-53 mm variable.

7x7 divider holds 49 vials up to 16.5 mm diameter. (See also page 9: KEA50)

Cat. No.	Colour
KEA45-V49NA	natural
KEA45-V49NY	neon-yellow
KEA45-V49NB	neon-blue
KEA45-V49NG	neon-green
KEA45-V49NP	neon-red/pink
KEA45-V49NO	neon-orange
KEA45-V49WH	white

Cat. No.	Colour
KEA45-V49YE	yellow
KEA45-V49BU	blue
KEA45-V49GR	green
KEA45-V49RE	red
KEA45-V49VT	violet
KEA45-V49BL	black
KEA45-V49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 45 / 7x7 • height 45-53 mm variable.

7x7 holes for 49 vials up to 16.5 mm diameter. With alphanumeric ID code. (See also page 9: KEB50)

Cat. No.	Colour
KEB45-V049NA	natural
KEB45-V049NY	neon-yellow
KEB45-V049NB	neon-blue
KEB45-V049NG	neon-green
KEB45-V049NP	neon-red/pink
KEB45-V049NO	neon-orange
KEB45-V049WH	white

Cat. No.	Colour
KEB45-V049YE	yellow
KEB45-V049BU	blue
KEB45-V049GR	green
KEB45-V049RE	red
KEB45-V049VT	violet
KEB45-V049BL	black
KEB45-V049BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 45 / 1x1 • height 45-53 mm variable.

Without grid divider. Useful for EPi® holes inserts or grid dividers (page 27). (See also page 9: KEA50)

Cat. No.	Colour
KEA45-V1NA	natural
KEA45-V1NY	neon-yellow
KEA45-V1NB	neon-blue
KEA45-V1NG	neon-green
KEA45-V1NP	neon-red/pink
KEA45-V1NO	neon-orange
KEA45-V1WH	white

Cat. No.	Colour
KEA45-V1YE	yellow
KEA45-V1BU	blue
KEA45-V1GR	green
KEA45-V1RE	red
KEA45-V1VT	violet
KEA45-V1BL	black
KEA45-V1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks


Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

EPi® Box made of Polypropylene

130 x 130 mm




EPi® 50 / 8x8 • height 52 mm fix.

Inner height 50 mm. 8x8 holes for 64 reaction vessels 1.5 ml and 2.0 ml Safe-Lock® up to 11 mm diameter. With alphanumeric ID code. (See also page 8: KEB45)

Cat. No.	Colour
KEB50-C064NA	natural
KEB50-C064NY	neon-yellow
KEB50-C064NB	neon-blue
KEB50-C064NG	neon-green
KEB50-C064NP	neon-red/pink
KEB50-C064NO	neon-orange
KEB50-C064WH	white

Cat. No.	Colour
KEB50-C064YE	yellow
KEB50-C064BU	blue
KEB50-C064GR	green
KEB50-C064RE	red
KEB50-C064VT	violet
KEB50-C064BL	black
KEB50-C064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 50 / 7x7 • height 52 mm fix.

Inner height 50 mm. 7x7 divider holds 49 vials up to 16.5 mm diameter. (See also page 8: KEA45)

Cat. No.	Colour
KEA50-C49NA	natural
KEA50-C49NY	neon-yellow
KEA50-C49NB	neon-blue
KEA50-C49NG	neon-green
KEA50-C49NP	neon-red/pink
KEA50-C49NO	neon-orange
KEA50-C49WH	white

Cat. No.	Colour
KEA50-C49YE	yellow
KEA50-C49BU	blue
KEA50-C49GR	green
KEA50-C49RE	red
KEA50-C49VT	violet
KEA50-C49BL	black
KEA50-C49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 50 / 7x7 • height 52 mm fix.

Inner height 50 mm. 7x7 holes for 49 vials up to 16.5 mm diameter. With alphanumeric ID code. (See also page 8: KEB45)

Cat. No.	Colour
KEB50-C049NA	natural
KEB50-C049NY	neon-yellow
KEB50-C049NB	neon-blue
KEB50-C049NG	neon-green
KEB50-C049NP	neon-red/pink
KEB50-C049NO	neon-orange
KEB50-C049WH	white

Cat. No.	Colour
KEB50-C049YE	yellow
KEB50-C049BU	blue
KEB50-C049GR	green
KEB50-C049RE	red
KEB50-C049VT	violet
KEB50-C049BL	black
KEB50-C049BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 50 / 1x1 • height 52 mm fix.

Inner height 50 mm. Without grid divider. Useful for EPi® holes inserts or grid dividers (page 27). (See also page 8: KEA45)

Cat. No.	Colour
KEA50-C1NA	natural
KEA50-C1NY	neon-yellow
KEA50-C1NB	neon-blue
KEA50-C1NG	neon-green
KEA50-C1NP	neon-red/pink
KEA50-C1NO	neon-orange
KEA50-C1WH	white

Cat. No.	Colour
KEA50-C1YE	yellow
KEA50-C1BU	blue
KEA50-C1GR	green
KEA50-C1RE	red
KEA50-C1VT	violet
KEA50-C1BL	black
KEA50-C1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

EPi® Box made of Polypropylene

130 x 130 mm



EPi® 70 / 9x9 • height 70-80 mm variable.

9x9 divider holds 81 vials up to 12.5 mm diameter.

Cat. No.	Colour
KEA70-V81NA	natural
KEA70-V81NY	neon-yellow
KEA70-V81NB	neon-blue
KEA70-V81NG	neon-green
KEA70-V81NP	neon-red/pink
KEA70-V81NO	neon-orange
KEA70-V81WH	white

Cat. No.	Colour
KEA70-V81YE	yellow
KEA70-V81BU	blue
KEA70-V81GR	green
KEA70-V81RE	red
KEA70-V81VT	violet
KEA70-V81BL	black
KEA70-V81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 70 / 8x8 • height 70-80 mm variable.

8x8 holes for 64 vials up to 11 mm diameter.
With alphanumeric ID code.

Cat. No.	Colour
KEB70-V064NA	natural
KEB70-V064NY	neon-yellow
KEB70-V064NB	neon-blue
KEB70-V064NG	neon-green
KEB70-V064NP	neon-red/pink
KEB70-V064NO	neon-orange
KEB70-V064WH	white

Cat. No.	Colour
KEB70-V064YE	yellow
KEB70-V064BU	blue
KEB70-V064GR	green
KEB70-V064RE	red
KEB70-V064VT	violet
KEB70-V064BL	black
KEB70-V064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 70 / 7x7 • height 70-80 mm variable.

7x7 divider holds 49 vials up to 16.5 mm diameter.

Cat. No.	Colour
KEA70-V49NA	natural
KEA70-V49NY	neon-yellow
KEA70-V49NB	neon-blue
KEA70-V49NG	neon-green
KEA70-V49NP	neon-red/pink
KEA70-V49NO	neon-orange
KEA70-V49WH	white

Cat. No.	Colour
KEA70-V49YE	yellow
KEA70-V49BU	blue
KEA70-V49GR	green
KEA70-V49RE	red
KEA70-V49VT	violet
KEA70-V49BL	black
KEA70-V49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 70 / 7x7 • height 70-80 mm variable.

7x7 holes for 49 vials up to 16.5 mm diameter.
With alphanumeric ID code.

Cat. No.	Colour
KEB70-V049NA	natural
KEB70-V049NY	neon-yellow
KEB70-V049NB	neon-blue
KEB70-V049NG	neon-green
KEB70-V049NP	neon-red/pink
KEB70-V049NO	neon-orange
KEB70-V049WH	white

Cat. No.	Colour
KEB70-V049YE	yellow
KEB70-V049BU	blue
KEB70-V049GR	green
KEB70-V049RE	red
KEB70-V049VT	violet
KEB70-V049BL	black
KEB70-V049BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks


Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

EPPi® Box made of Polypropylene

130 x 130 mm



EPPi® 70 / 1x1 • height 70-80 mm variable.

Without grid divider. Useful for EPPi® holes inserts or grid dividers (page 27).

Cat. No.	Colour	Cat. No.	Colour
KEA70-V1NA	natural	KEA70-V1YE	yellow
KEA70-V1NY	neon-yellow	KEA70-V1BU	blue
KEA70-V1NB	neon-blue	KEA70-V1GR	green
KEA70-V1NG	neon-green	KEA70-V1RE	red
KEA70-V1NP	neon-red/pink	KEA70-V1VT	violet
KEA70-V1NO	neon-orange	KEA70-V1BL	black
KEA70-V1WH	white	KEA70-V1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35



EPPi® Ident Identification Labels down to -90°C

Self-adhesive labels for marking boxes made of plastic and fibreboard as well as stainless steel racks, etc. Paper, white, printable.

Supplied on sheets in A4 format in packs of 50 sheets each.

Cat. No.	W x H	Quantity	Cat. No.	W x H	Quantity
LAL2813WH50	28x13 mm	50x154	LAL6006WH50	60x 6 mm	50x135
LAL5230WH50	52x30 mm	50x36	LAL4013WH50	40x13 mm	50x115
LAL7015WH50	70x15 mm	50x57	LAL4025WH50	40x25 mm	50x60
LAL7020WH50	70x20 mm	50x42	LAL4812WH50	48x12 mm	50x96
LAL7035WH50	70x35 mm	50x24	LAL3512WH50	35x12 mm	50x144
LAL6020WH50	60x20 mm	50x42	LAL010WH50	Ø10 mm	50x100



page 32+33



page 30+31




page 65



page 68

EPi® Box made of Polypropylene

130 x 130 mm




EPi® 75 / 9x9 • height 75 mm fix.

Inner height 73 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. Without ID code.

Cat. No.	Colour
KEA75-C81NA	natural
KEA75-C81NY	neon-yellow
KEA75-C81NB	neon-blue
KEA75-C81NG	neon-green
KEA75-C81NP	neon-red/pink
KEA75-C81NO	neon-orange
KEA75-C81WH	white

Cat. No.	Colour
KEA75-C81YE	yellow
KEA75-C81BU	blue
KEA75-C81GR	green
KEA75-C81RE	red
KEA75-C81VT	violet
KEA75-C81BL	black
KEA75-C81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 75 / 9x9 • height 75 mm fix.

Inner height 73 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides.

Cat. No.	Colour
KEB75-C81NA	natural
KEB75-C81NY	neon-yellow
KEB75-C81NB	neon-blue
KEB75-C81NG	neon-green
KEB75-C81NP	neon-red/pink
KEB75-C81NO	neon-orange
KEB75-C81WH	white

Cat. No.	Colour
KEB75-C81YE	yellow
KEB75-C81BU	blue
KEB75-C81GR	green
KEB75-C81RE	red
KEB75-C81VT	violet
KEB75-C81BL	black
KEB75-C81BB	black/black*


*black/black: Lower half and lid black for light-sensitive samples



EPi® 75 / 9x9 • height 75 mm fix.

Inner height 73 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the bottom of the box (transparent).

Cat. No.	Colour
KEC75-C81NA	natural
KEC75-C81NY	neon-yellow
KEC75-C81NB	neon-blue
KEC75-C81NG	neon-green
KEC75-C81NP	neon-red/pink
KEC75-C81NO	neon-orange



EPi® 75 / 9x9 • height 75 mm fix.

Inner height 73 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the bottom of the box (transparent) and on the lid.

Cat. No.	Colour
KEE75-C81NA	natural
KEE75-C81NY	neon-yellow
KEE75-C81NB	neon-blue
KEE75-C81NG	neon-green
KEE75-C81NP	neon-red/pink
KEE75-C81NO	neon-orange

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31



Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35



EPi® 75 / 9x9 • height 75 mm fix.

Inner height 73 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on the lid.

Cat. No.	Colour
KEH75-C81NA	natural
KEH75-C81NY	neon-yellow
KEH75-C81NB	neon-blue
KEH75-C81NG	neon-green
KEH75-C81NP	neon-red/pink
KEH75-C81NO	neon-orange
KEH75-C81WH	white

Cat. No.	Colour
KEH75-C81YE	yellow
KEH75-C81BU	blue
KEH75-C81GR	green
KEH75-C81RE	red
KEH75-C81VT	violet
KEH75-C81BL	black

EPi® 75 / 9x9 • height 75 mm fix.

Inner height 73 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the lid.

Cat. No.	Colour
KEK75-C81NA	natural
KEK75-C81NY	neon-yellow
KEK75-C81NB	neon-blue
KEK75-C81NG	neon-green
KEK75-C81NP	neon-red/pink
KEK75-C81NO	neon-orange
KEK75-C81WH	white

Cat. No.	Colour
KEK75-C81YE	yellow
KEK75-C81BU	blue
KEK75-C81GR	green
KEK75-C81RE	red
KEK75-C81VT	violet
KEK75-C81BL	black

EPi® 75 / 9x9 • height 75 mm fix.

Inner height 73 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With numeric ID code (1-81) on the lid.

Cat. No.	Colour
KEM75-C81NA	natural
KEM75-C81NY	neon-yellow
KEM75-C81NB	neon-blue
KEM75-C81NG	neon-green
KEM75-C81NP	neon-red/pink
KEM75-C81NO	neon-orange
KEM75-C81WH	white

Cat. No.	Colour
KEM75-C81YE	yellow
KEM75-C81BU	blue
KEM75-C81GR	green
KEM75-C81RE	red
KEM75-C81VT	violet
KEM75-C81BL	black

EPi® 75 / 9x9 • height 75 mm fix.

Inner height 73 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With numeric ID code (1-81) on the bottom of the box (transparent) and on the lid.

Cat. No.	Colour
KER75-C81NA	natural
KER75-C81NY	neon-yellow
KER75-C81NB	neon-blue
KER75-C81NG	neon-green
KER75-C81NP	neon-red/pink
KER75-C81NO	neon-orange

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

EPi® Box made of Polypropylene

130 x 130 mm




EPi® 75 / 8x8 • height 75 mm fix.

Inner height 73 mm. 8x8 holes for 64 vials up to 11 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB75-C064NA	natural
KEB75-C064NY	neon-yellow
KEB75-C064NB	neon-blue
KEB75-C064NG	neon-green
KEB75-C064NP	neon-red/pink
KEB75-C064NO	neon-orange
KEB75-C064WH	white

Cat. No.	Colour
KEB75-C064YE	yellow
KEB75-C064BU	blue
KEB75-C064GR	green
KEB75-C064RE	red
KEB75-C064VT	violet
KEB75-C064BL	black
KEB75-C064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 75 / 7x7 • height 75 mm fix.

Inner height 73 mm. 7x7 divider holds 49 vials up to 16.5 mm diameter.

Cat. No.	Colour
KEA75-C49NA	natural
KEA75-C49NY	neon-yellow
KEA75-C49NB	neon-blue
KEA75-C49NG	neon-green
KEA75-C49NP	neon-red/pink
KEA75-C49NO	neon-orange
KEA75-C49WH	white

Cat. No.	Colour
KEA75-C49YE	yellow
KEA75-C49BU	blue
KEA75-C49GR	green
KEA75-C49RE	red
KEA75-C49VT	violet
KEA75-C49BL	black
KEA75-C49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 75 / 7x7 • height 75 mm fix.

Inner height 73 mm. 7x7 holes for 49 vials up to 16.5 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB75-C049NA	natural
KEB75-C049NY	neon-yellow
KEB75-C049NB	neon-blue
KEB75-C049NG	neon-green
KEB75-C049NP	neon-red/pink
KEB75-C049NO	neon-orange
KEB75-C049WH	white

Cat. No.	Colour
KEB75-C049YE	yellow
KEB75-C049BU	blue
KEB75-C049GR	green
KEB75-C049RE	red
KEB75-C049VT	violet
KEB75-C049BL	black
KEB75-C049BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 75 / 1x1 • height 75 mm fix.

Inner height 73 mm. Without grid divider. Useful for EPi® holes inserts or grid dividers (page 27).

Cat. No.	Colour
KEA75-C1NA	natural
KEA75-C1NY	neon-yellow
KEA75-C1NB	neon-blue
KEA75-C1NG	neon-green
KEA75-C1NP	neon-red/pink
KEA75-C1NO	neon-orange
KEA75-C1WH	white

Cat. No.	Colour
KEA75-C1YE	yellow
KEA75-C1BU	blue
KEA75-C1GR	green
KEA75-C1RE	red
KEA75-C1VT	violet
KEA75-C1BL	black
KEA75-C1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31


Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35




EPi® 80 / 9x9 • height 80 mm fix.

Inner height 77 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter.

Cat. No.	Colour
KEA80-C81NA	natural
KEA80-C81NY	neon-yellow
KEA80-C81NB	neon-blue
KEA80-C81NG	neon-green
KEA80-C81NP	neon-red/pink
KEA80-C81NO	neon-orange
KEA80-C81WH	white

Cat. No.	Colour
KEA80-C81YE	yellow
KEA80-C81BU	blue
KEA80-C81GR	green
KEA80-C81RE	red
KEA80-C81VT	violet
KEA80-C81BL	black
KEA80-C81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 80 / 8x8 • height 80 mm fix.

Inner height 77 mm. 8x8 holes for 64 vials up to 11 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB80-C064NA	natural
KEB80-C064NY	neon-yellow
KEB80-C064NB	neon-blue
KEB80-C064NG	neon-green
KEB80-C064NP	neon-red/pink
KEB80-C064NO	neon-orange
KEB80-C064WH	white

Cat. No.	Colour
KEB80-C064YE	yellow
KEB80-C064BU	blue
KEB80-C064GR	green
KEB80-C064RE	red
KEB80-C064VT	violet
KEB80-C064BL	black
KEB80-C064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 80 / 7x7 • height 80 mm fix.

Inner height 77 mm. 7x7 divider holds 49 vials up to 16.5 mm diameter.

Cat. No.	Colour
KEA80-C49NA	natural
KEA80-C49NY	neon-yellow
KEA80-C49NB	neon-blue
KEA80-C49NG	neon-green
KEA80-C49NP	neon-red/pink
KEA80-C49NO	neon-orange
KEA80-C49WH	white

Cat. No.	Colour
KEA80-C49YE	yellow
KEA80-C49BU	blue
KEA80-C49GR	green
KEA80-C49RE	red
KEA80-C49VT	violet
KEA80-C49BL	black
KEA80-C49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 80 / 7x7 • height 80 mm fix.

Inner height 77 mm. 7x7 holes for 49 vials up to 16.5 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB80-C049NA	natural
KEB80-C049NY	neon-yellow
KEB80-C049NB	neon-blue
KEB80-C049NG	neon-green
KEB80-C049NP	neon-red/pink
KEB80-C049NO	neon-orange
KEB80-C049WH	white

Cat. No.	Colour
KEB80-C049YE	yellow
KEB80-C049BU	blue
KEB80-C049GR	green
KEB80-C049RE	red
KEB80-C049VT	violet
KEB80-C049BL	black
KEB80-C049BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

EPi® Box made of Polypropylene

130 x 130 mm




EPi® 80 / 1x1 • height 80 mm fix.

Inner height 77 mm. Without grid divider. Useful for EPi® holes inserts or grid dividers (page 27).

Cat. No.	Colour
KEA80-C1NA	natural
KEA80-C1NY	neon-yellow
KEA80-C1NB	neon-blue
KEA80-C1NG	neon-green
KEA80-C1NP	neon-red/pink
KEA80-C1NO	neon-orange
KEA80-C1WH	white

Cat. No.	Colour
KEA80-C1YE	yellow
KEA80-C1BU	blue
KEA80-C1GR	green
KEA80-C1RE	red
KEA80-C1VT	violet
KEA80-C1BL	black
KEA80-C1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 95 / 9x9 • height 95 mm fix.

Inner height 93 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. Such as for 5 ml cryo vials.

Cat. No.	Colour
KEA95-C81NA	natural
KEA95-C81NY	neon-yellow
KEA95-C81NB	neon-blue
KEA95-C81NG	neon-green
KEA95-C81NP	neon-red/pink
KEA95-C81NO	neon-orange
KEA95-C81WH	white

Cat. No.	Colour
KEA95-C81YE	yellow
KEA95-C81BU	blue
KEA95-C81GR	green
KEA95-C81RE	red
KEA95-C81VT	violet
KEA95-C81BL	black
KEA95-C81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 95 / 9x9 • height 95 mm fix.

Inner height 93 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides. Such as for 5 ml cryo vials.

Cat. No.	Colour
KEB95-C81NA	natural
KEB95-C81NY	neon-yellow
KEB95-C81NB	neon-blue
KEB95-C81NG	neon-green
KEB95-C81NP	neon-red/pink
KEB95-C81NO	neon-orange
KEB95-C81WH	white

Cat. No.	Colour
KEB95-C81YE	yellow
KEB95-C81BU	blue
KEB95-C81GR	green
KEB95-C81RE	red
KEB95-C81VT	violet
KEB95-C81BL	black
KEB95-C81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 95 / 9x9 • height 95 mm fix.

Inner height 93 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the bottom of the box (transparent). Such as for 5 ml cryo vials.

Cat. No.	Colour
KEC95-C81NA	natural
KEC95-C81NY	neon-yellow
KEC95-C81NB	neon-blue
KEC95-C81NG	neon-green
KEC95-C81NP	neon-red/pink
KEC95-C81NO	neon-orange

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 35

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31


Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31


Drawer Racks
page 36



EPi® 95 / 9x9 • height 95 mm fix.

Inner height 93 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the bottom of the box (transparent) and on the lid. Such as for 5 ml cryo vials.

Cat. No.	Colour
KEE95-C81NA	natural
KEE95-C81NY	neon-yellow
KEE95-C81NB	neon-blue
KEE95-C81NG	neon-green
KEE95-C81NP	neon-red/pink
KEE95-C81NO	neon-orange




EPi® 95 / 9x9 • height 95 mm fix.

Inner height 93 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on the lid. Such as for 5 ml cryo vials.

Cat. No.	Colour
KEH95-C81NA	natural
KEH95-C81NY	neon-yellow
KEH95-C81NB	neon-blue
KEH95-C81NG	neon-green
KEH95-C81NP	neon-red/pink
KEH95-C81NO	neon-orange
KEH95-C81WH	white

Cat. No.	Colour
KEH95-C81YE	yellow
KEH95-C81BU	blue
KEH95-C81GR	green
KEH95-C81RE	red
KEH95-C81VT	violet
KEH95-C81BL	black




EPi® 95 / 9x9 • height 95 mm fix.

Inner height 93 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With alphanumeric ID code on all 4 sides as well as on the lid. Such as for 5 ml cryo vials.

Cat. No.	Colour
KEK95-C81NA	natural
KEK95-C81NY	neon-yellow
KEK95-C81NB	neon-blue
KEK95-C81NG	neon-green
KEK95-C81NP	neon-red/pink
KEK95-C81NO	neon-orange
KEK95-C81WH	white

Cat. No.	Colour
KEK95-C81YE	yellow
KEK95-C81BU	blue
KEK95-C81GR	green
KEK95-C81RE	red
KEK95-C81VT	violet
KEK95-C81BL	black



EPi® 95 / 9x9 • height 95 mm fix.

Inner height 93 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With numeric ID code (1-81) on the lid. Such as for 5 ml cryo vials.

Cat. No.	Colour
KEM95-C81NA	natural
KEM95-C81NY	neon-yellow
KEM95-C81NB	neon-blue
KEM95-C81NG	neon-green
KEM95-C81NP	neon-red/pink
KEM95-C81NO	neon-orange
KEM95-C81WH	white

Cat. No.	Colour
KEM95-C81YE	yellow
KEM95-C81BU	blue
KEM95-C81GR	green
KEM95-C81RE	red
KEM95-C81VT	violet
KEM95-C81BL	black

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

EPPI® Box made of Polypropylene

130 x 130 mm



EPPI® 95 / 9x9 • height 95 mm fix.

Inner height 93 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter. With numeric ID code (1-81) on the bottom of the box (transparent) and on the lid. Such as for 5 ml cryo vials.

Cat. No.	Colour
KER95-C81NA	natural
KER95-C81NY	neon-yellow
KER95-C81NB	neon-blue
KER95-C81NG	neon-green
KER95-C81NP	neon-red/pink
KER95-C81NO	neon-orange



EPPI® 95 / 8x8 • height 95 mm fix.

Inner height 93 mm. 8x8 holes for 64 vials up to 11 mm diameter. With alphanumeric ID code.

Cat. No.	Colour	Cat. No.	Colour
KEB95-C064NA	natural	KEB95-C064YE	yellow
KEB95-C064NY	neon-yellow	KEB95-C064BU	blue
KEB95-C064NB	neon-blue	KEB95-C064GR	green
KEB95-C064NG	neon-green	KEB95-C064RE	red
KEB95-C064NP	neon-red/pink	KEB95-C064VT	violet
KEB95-C064NO	neon-orange	KEB95-C064BL	black
KEB95-C064WH	white	KEB95-C064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPPI® 95 / 7x7 • height 95 mm fix.

Inner height 93 mm. 7x7 divider holds 49 vials up to 16.5 mm diameter.

Cat. No.	Colour	Cat. No.	Colour
KEA95-C49NA	natural	KEA95-C49YE	yellow
KEA95-C49NY	neon-yellow	KEA95-C49BU	blue
KEA95-C49NB	neon-blue	KEA95-C49GR	green
KEA95-C49NG	neon-green	KEA95-C49RE	red
KEA95-C49NP	neon-red/pink	KEA95-C49VT	violet
KEA95-C49NO	neon-orange	KEA95-C49BL	black
KEA95-C49WH	white	KEA95-C49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPPI® 95 / 1x1 • height 95 mm fix.

Inner height 93 mm. Without grid divider. Useful for EPPI® holes inserts or grid dividers (page 27).

Cat. No.	Colour	Cat. No.	Colour
KEA95-C1NA	natural	KEA95-C1YE	yellow
KEA95-C1NY	neon-yellow	KEA95-C1BU	blue
KEA95-C1NB	neon-blue	KEA95-C1GR	green
KEA95-C1NG	neon-green	KEA95-C1RE	red
KEA95-C1NP	neon-red/pink	KEA95-C1VT	violet
KEA95-C1NO	neon-orange	KEA95-C1BL	black
KEA95-C1WH	white	KEA95-C1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36



EPi® 96 / 9x9 • height 96-106 mm variable.

9x9 divider holds 81 vials up to 12.5 mm diameter.
(See also page 16: KEA95)

Cat. No.	Colour
KEA96-V81NA	natural
KEA96-V81NY	neon-yellow
KEA96-V81NB	neon-blue
KEA96-V81NG	neon-green
KEA96-V81NP	neon-red/pink
KEA96-V81NO	neon-orange
KEA96-V81WH	white

Cat. No.	Colour
KEA96-V81YE	yellow
KEA96-V81BU	blue
KEA96-V81GR	green
KEA96-V81RE	red
KEA96-V81VT	violet
KEA96-V81BL	black
KEA96-V81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 96 / 8x8 • height 96-106 mm variable.

8x8 holes for 81 vials up to 11 mm diameter. With alphanumeric ID code. (See also page 18: KEB95)

Cat. No.	Colour
KEB96-V064NA	natural
KEB96-V064NY	neon-yellow
KEB96-V064NB	neon-blue
KEB96-V064NG	neon-green
KEB96-V064NP	neon-red/pink
KEB96-V064NO	neon-orange
KEB96-V064WH	white

Cat. No.	Colour
KEB96-V064YE	yellow
KEB96-V064BU	blue
KEB96-V064GR	green
KEB96-V064RE	red
KEB96-V064VT	violet
KEB96-V064BL	black
KEB96-V064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 96 / 7x7 • height 96-106 mm variable.

7x7 divider holds 49 vials up to 16.5 mm diameter.
(See also page 18: KEA95)

Cat. No.	Colour
KEA96-V49NA	natural
KEA96-V49NY	neon-yellow
KEA96-V49NB	neon-blue
KEA96-V49NG	neon-green
KEA96-V49NP	neon-red/pink
KEA96-V49NO	neon-orange
KEA96-V49WH	white

Cat. No.	Colour
KEA96-V49YE	yellow
KEA96-V49BU	blue
KEA96-V49GR	green
KEA96-V49RE	red
KEA96-V49VT	violet
KEA96-V49BL	black
KEA96-V49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 96 / 7x7 • height 96-106 mm variable.

7x7 holes for 49 vials up to 16.5 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB96-V049NA	natural
KEB96-V049NY	neon-yellow
KEB96-V049NB	neon-blue
KEB96-V049NG	neon-green
KEB96-V049NP	neon-red/pink
KEB96-V049NO	neon-orange
KEB96-V049WH	white

Cat. No.	Colour
KEB96-V049YE	yellow
KEB96-V049BU	blue
KEB96-V049GR	green
KEB96-V049RE	red
KEB96-V049VT	violet
KEB96-V049BL	black
KEB96-V049BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

EPi® Box made of Polypropylene

130 x 130 mm




EPi® 96 / 5x5 • height 96-106 mm variable.

5x5 holes for 25 culture vials up to 17 mm diameter.

Cat. No.	Colour
KEA96-V025NA	natural
KEA96-V025NY	neon-yellow
KEA96-V025NB	neon-blue
KEA96-V025NG	neon-green
KEA96-V025NP	neon-red/pink
KEA96-V025NO	neon-orange
KEA96-V025WH	white

Cat. No.	Colour
KEA96-V025YE	yellow
KEA96-V025BU	blue
KEA96-V025GR	green
KEA96-V025RE	red
KEA96-V025VT	violet
KEA96-V025BL	black
KEA96-V025BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 96 / 10 L • height 96-106 mm variable.

For 10 culture vials up to 30 mm diameter.

Cat. No.	Colour
KEA96-V010NA	natural
KEA96-V010NY	neon-yellow
KEA96-V010NB	neon-blue
KEA96-V010NG	neon-green
KEA96-V010NP	neon-red/pink
KEA96-V010NO	neon-orange
KEA96-V010WH	white

Cat. No.	Colour
KEA96-V010YE	yellow
KEA96-V010BU	blue
KEA96-V010GR	green
KEA96-V010RE	red
KEA96-V010VT	violet
KEA96-V010BL	black
KEA96-V010BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 96 / 1x1 • height 96-106 mm variable.

Without grid divider. Useful for EPi® holes inserts or grid dividers (page 27).

Cat. No.	Colour
KEA96-V1NA	natural
KEA96-V1NY	neon-yellow
KEA96-V1NB	neon-blue
KEA96-V1NG	neon-green
KEA96-V1NP	neon-red/pink
KEA96-V1NO	neon-orange
KEA96-V1WH	white

Cat. No.	Colour
KEA96-V1YE	yellow
KEA96-V1BU	blue
KEA96-V1GR	green
KEA96-V1RE	red
KEA96-V1VT	violet
KEA96-V1BL	black
KEA96-V1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 102 / 9x9 • height 102 mm fix.

Inner height 100 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter.

Cat. No.	Colour
KEA102-C81NA	natural
KEA102-C81NY	neon-yellow
KEA102-C81NB	neon-blue
KEA102-C81NG	neon-green
KEA102-C81NP	neon-red/pink
KEA102-C81NO	neon-orange
KEA102-C81WH	white

Cat. No.	Colour
KEA102-C81YE	yellow
KEA102-C81BU	blue
KEA102-C81GR	green
KEA102-C81RE	red
KEA102-C81VT	violet
KEA102-C81BL	black
KEA102-C81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36



EPPI® 102 / 8x8 • height 102 mm fix.

Inner height 100 mm. 8x8 holes for 64 vials up to 11 mm diameter. With alphanumeric ID code.

Cat. No.	Colour	Cat. No.	Colour
KEB102-C064NA	natural	KEB102-C064YE	yellow
KEB102-C064NY	neon-yellow	KEB102-C064BU	blue
KEB102-C064NB	neon-blue	KEB102-C064GR	green
KEB102-C064NG	neon-green	KEB102-C064RE	red
KEB102-C064NP	neon-red/pink	KEB102-C064VT	violet
KEB102-C064NO	neon-orange	KEB102-C064BL	black
KEB102-C064WH	white	KEB102-C064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPPI® 102 / 7x7 • height 102 mm fix.

Inner height 100 mm. 7x7 divider holds 49 vials up to 16.5 mm diameter.

Cat. No.	Colour	Cat. No.	Colour
KEA102-C49NA	natural	KEA102-C49YE	yellow
KEA102-C49NY	neon-yellow	KEA102-C49BU	blue
KEA102-C49NB	neon-blue	KEA102-C49GR	green
KEA102-C49NG	neon-green	KEA102-C49RE	red
KEA102-C49NP	neon-red/pink	KEA102-C49VT	violet
KEA102-C49NO	neon-orange	KEA102-C49BL	black
KEA102-C49WH	white	KEA102-C49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPPI® 102 / 7x7 • height 102 mm fix.

Inner height 100 mm. 7x7 holes for 49 vials up to 16.5 mm diameter. With alphanumeric ID code.

Cat. No.	Colour	Cat. No.	Colour
KEB102-C049NA	natural	KEB102-C049YE	yellow
KEB102-C049NY	neon-yellow	KEB102-C049BU	blue
KEB102-C049NB	neon-blue	KEB102-C049GR	green
KEB102-C049NG	neon-green	KEB102-C049RE	red
KEB102-C049NP	neon-red/pink	KEB102-C049VT	violet
KEB102-C049NO	neon-orange	KEB102-C049BL	black
KEB102-C049WH	white	KEB102-C049BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPPI® 102 / 1x1 • height 102 mm fix.

Inner height 100 mm. Without grid divider. Useful for EPPI® holes inserts or grid dividers (page 27).

Cat. No.	Colour	Cat. No.	Colour
KEA102-C1NA	natural	KEA102-C1YE	yellow
KEA102-C1NY	neon-yellow	KEA102-C1BU	blue
KEA102-C1NB	neon-blue	KEA102-C1GR	green
KEA102-C1NG	neon-green	KEA102-C1RE	red
KEA102-C1NP	neon-red/pink	KEA102-C1VT	violet
KEA102-C1NO	neon-orange	KEA102-C1BL	black
KEA102-C1WH	white	KEA102-C1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36



EPi® 105 / 9x9 • height 105 mm fix.

Inner height 103 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter.

Cat. No.	Colour
KEA105-C81NA	natural
KEA105-C81NY	neon-yellow
KEA105-C81NB	neon-blue
KEA105-C81NG	neon-green
KEA105-C81NP	neon-red/pink
KEA105-C81NO	neon-orange
KEA105-C81WH	white

Cat. No.	Colour
KEA105-C81YE	yellow
KEA105-C81BU	blue
KEA105-C81GR	green
KEA105-C81RE	red
KEA105-C81VT	violet
KEA105-C81BL	black
KEA105-C81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 105 / 8x8 • height 105 mm fix.

Inner height 103 mm. 8x8 holes for 64 vials up to 11 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB105-C064NA	natural
KEB105-C064NY	neon-yellow
KEB105-C064NB	neon-blue
KEB105-C064NG	neon-green
KEB105-C064NP	neon-red/pink
KEB105-C064NO	neon-orange
KEB105-C064WH	white

Cat. No.	Colour
KEB105-C064YE	yellow
KEB105-C064BU	blue
KEB105-C064GR	green
KEB105-C064RE	red
KEB105-C064VT	violet
KEB105-C064BL	black
KEB105-C064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 105 / 7x7 • height 105 mm fix.

Inner height 103 mm. 7x7 divider holds 49 vials up to 16.5 mm diameter.

Cat. No.	Colour
KEA105-C49NA	natural
KEA105-C49NY	neon-yellow
KEA105-C49NB	neon-blue
KEA105-C49NG	neon-green
KEA105-C49NP	neon-red/pink
KEA105-C49NO	neon-orange
KEA105-C49WH	white

Cat. No.	Colour
KEA105-C49YE	yellow
KEA105-C49BU	blue
KEA105-C49GR	green
KEA105-C49RE	red
KEA105-C49VT	violet
KEA105-C49BL	black
KEA105-C49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 105 / 7x7 • height 105 mm fix.

Inner height 103 mm. 7x7 holes for 49 vials up to 16.5 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB105-C049NA	natural
KEB105-C049NY	neon-yellow
KEB105-C049NB	neon-blue
KEB105-C049NG	neon-green
KEB105-C049NP	neon-red/pink
KEB105-C049NO	neon-orange
KEB105-C049WH	white

Cat. No.	Colour
KEB105-C049YE	yellow
KEB105-C049BU	blue
KEB105-C049GR	green
KEB105-C049RE	red
KEB105-C049VT	violet
KEB105-C049BL	black
KEB105-C049BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31


Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37



EPi® 105 / 5x5 • height 105 mm fix.

Inner height 103 mm. For 25 vials up to 17 mm diameter. Such as for culture vials 14 ml.

Cat. No.	Colour
KEA105-C025NA	natural
KEA105-C025NY	neon-yellow
KEA105-C025NB	neon-blue
KEA105-C025NG	neon-green
KEA105-C025NP	neon-red/pink
KEA105-C025NO	neon-orange
KEA105-C025WH	white

Cat. No.	Colour
KEA105-C025YE	yellow
KEA105-C025BU	blue
KEA105-C025GR	green
KEA105-C025RE	red
KEA105-C025VT	violet
KEA105-C025BL	black
KEA105-C025BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 105 / 10 L • height 105 mm fix.

Inner height 103 mm. For 10 vials up to 30 mm diameter.

Cat. No.	Colour
KEA105-C010NA	natural
KEA105-C010NY	neon-yellow
KEA105-C010NB	neon-blue
KEA105-C010NG	neon-green
KEA105-C010NP	neon-red/pink
KEA105-C010NO	neon-orange
KEA105-C010WH	white

Cat. No.	Colour
KEA105-C010YE	yellow
KEA105-C010BU	blue
KEA105-C010GR	green
KEA105-C010RE	red
KEA105-C010VT	violet
KEA105-C010BL	black
KEA105-C010BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 105 / 1x1 • height 105 mm fix.

Inner height 103 mm. Without grid divider. Useful for EPi® holes inserts or grid dividers (page 27).

Cat. No.	Colour
KEA105-C1NA	natural
KEA105-C1NY	neon-yellow
KEA105-C1NB	neon-blue
KEA105-C1NG	neon-green
KEA105-C1NP	neon-red/pink
KEA105-C1NO	neon-orange
KEA105-C1WH	white

Cat. No.	Colour
KEA105-C1YE	yellow
KEA105-C1BU	blue
KEA105-C1GR	green
KEA105-C1RE	red
KEA105-C1VT	violet
KEA105-C1BL	black
KEA105-C1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 122 / 9x9 • height 122 mm fix.

Inner height 118 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter.

Cat. No.	Colour
KEA122-C81NA	natural
KEA122-C81NY	neon-yellow
KEA122-C81NB	neon-blue
KEA122-C81NG	neon-green
KEA122-C81NP	neon-red/pink
KEA122-C81NO	neon-orange
KEA122-C81WH	white

Cat. No.	Colour
KEA122-C81YE	yellow
KEA122-C81BU	blue
KEA122-C81GR	green
KEA122-C81RE	red
KEA122-C81VT	violet
KEA122-C81BL	black
KEA122-C81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37




EPPI® 122 / 8x8 • height 122 mm fix.

Inner height 118 mm. 8x8 holes for 64 vials up to 11 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB122-C064NA	natural
KEB122-C064NY	neon-yellow
KEB122-C064NB	neon-blue
KEB122-C064NG	neon-green
KEB122-C064NP	neon-red/pink
KEB122-C064NO	neon-orange
KEB122-C064WH	white

Cat. No.	Colour
KEB122-C064YE	yellow
KEB122-C064BU	blue
KEB122-C064GR	green
KEB122-C064RE	red
KEB122-C064VT	violet
KEB122-C064BL	black
KEB122-C064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPPI® 122 / 7x7 • height 122 mm fix.

Inner height 118 mm. 7x7 divider holds 49 vials up to 16.5 mm diameter.

Cat. No.	Colour
KEA122-C49NA	natural
KEA122-C49NY	neon-yellow
KEA122-C49NB	neon-blue
KEA122-C49NG	neon-green
KEA122-C49NP	neon-red/pink
KEA122-C49NO	neon-orange
KEA122-C49WH	white

Cat. No.	Colour
KEA122-C49YE	yellow
KEA122-C49BU	blue
KEA122-C49GR	green
KEA122-C49RE	red
KEA122-C49VT	violet
KEA122-C49BL	black
KEA122-C49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPPI® 122 / 7x7 • height 122 mm fix.

Inner height 118 mm. 7x7 holes for 49 vials up to 16.5 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB122-C049NA	natural
KEB122-C049NY	neon-yellow
KEB122-C049NB	neon-blue
KEB122-C049NG	neon-green
KEB122-C049NP	neon-red/pink
KEB122-C049NO	neon-orange
KEB122-C049WH	white

Cat. No.	Colour
KEB122-C049YE	yellow
KEB122-C049BU	blue
KEB122-C049GR	green
KEB122-C049RE	red
KEB122-C049VT	violet
KEB122-C049BL	black
KEB122-C049BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPPI® 122 / 1x1 • height 122 mm fix.

Inner height 118 mm. Without grid divider. Useful for EPPI® holes inserts or grid dividers (page 27).

Cat. No.	Colour
KEA122-C1NA	natural
KEA122-C1NY	neon-yellow
KEA122-C1NB	neon-blue
KEA122-C1NG	neon-green
KEA122-C1NP	neon-red/pink
KEA122-C1NO	neon-orange
KEA122-C1WH	white

Cat. No.	Colour
KEA122-C1YE	yellow
KEA122-C1BU	blue
KEA122-C1GR	green
KEA122-C1RE	red
KEA122-C1VT	violet
KEA122-C1BL	black
KEA122-C1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37



EPi® 128 / 9x9 • height 128 mm fix.

Inner height 125 mm. 9x9 divider holds 81 vials up to 12.5 mm diameter.

Cat. No.	Colour
KEA128-C81NA	natural
KEA128-C81NY	neon-yellow
KEA128-C81NB	neon-blue
KEA128-C81NG	neon-green
KEA128-C81NP	neon-red/pink
KEA128-C81NO	neon-orange
KEA128-C81WH	white

Cat. No.	Colour
KEA128-C81YE	yellow
KEA128-C81BU	blue
KEA128-C81GR	green
KEA128-C81RE	red
KEA128-C81VT	violet
KEA128-C81BL	black
KEA128-C81BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 128 / 8x8 • height 128 mm fix.

Inner height 125 mm. 8x8 holes for 64 vials up to 11 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB128-C064NA	natural
KEB128-C064NY	neon-yellow
KEB128-C064NB	neon-blue
KEB128-C064NG	neon-green
KEB128-C064NP	neon-red/pink
KEB128-C064NO	neon-orange
KEB128-C064WH	white

Cat. No.	Colour
KEB128-C064YE	yellow
KEB128-C064BU	blue
KEB128-C064GR	green
KEB128-C064RE	red
KEB128-C064VT	violet
KEB128-C064BL	black
KEB128-C064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 128 / 7x7 • height 128 mm fix.

Inner height 125 mm. 7x7 divider holds 49 vials up to 16.5 mm diameter.

Cat. No.	Colour
KEA128-C49NA	natural
KEA128-C49NY	neon-yellow
KEA128-C49NB	neon-blue
KEA128-C49NG	neon-green
KEA128-C49NP	neon-red/pink
KEA128-C49NO	neon-orange
KEA128-C49WH	white

Cat. No.	Colour
KEA128-C49YE	yellow
KEA128-C49BU	blue
KEA128-C49GR	green
KEA128-C49RE	red
KEA128-C49VT	violet
KEA128-C49BL	black
KEA128-C49BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 128 / 7x7 • height 128 mm fix.

Inner height 125 mm. 7x7 holes for 49 vials up to 16.5 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB128-C049NA	natural
KEB128-C049NY	neon-yellow
KEB128-C049NB	neon-blue
KEB128-C049NG	neon-green
KEB128-C049NP	neon-red/pink
KEB128-C049NO	neon-orange
KEB128-C049WH	white

Cat. No.	Colour
KEB128-C049YE	yellow
KEB128-C049BU	blue
KEB128-C049GR	green
KEB128-C049RE	red
KEB128-C049VT	violet
KEB128-C049BL	black
KEB128-C049BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

EPi® Box made of Polypropylene

130 x 130 mm



EPi® 128 / 5x5 • height 128 mm fix.

Inner height 125 mm. For 25 centrifuge tubes 15 ml (up to 17 mm diameter).

Cat. No.	Colour
KEA128-C025NA	natural
KEA128-C025NY	neon-yellow
KEA128-C025NB	neon-blue
KEA128-C025NG	neon-green
KEA128-C025NP	neon-red/pink
KEA128-C025NO	neon-orange
KEA128-C025WH	white

Cat. No.	Colour
KEA128-C025YE	yellow
KEA128-C025BU	blue
KEA128-C025GR	green
KEA128-C025RE	red
KEA128-C025VT	violet
KEA128-C025BL	black
KEA128-C025BB	black/black*



EPi® 128 / 10 L • height 128 mm fix.

Inner height 125 mm. For 10 centrifuge tubes 50 ml (up to 30 mm diameter).

Cat. No.	Colour
KEA128-C010NA	natural
KEA128-C010NY	neon-yellow
KEA128-C010NB	neon-blue
KEA128-C010NG	neon-green
KEA128-C010NP	neon-red/pink
KEA128-C010NO	neon-orange
KEA128-C010WH	white

Cat. No.	Colour
KEA128-C010YE	yellow
KEA128-C010BU	blue
KEA128-C010GR	green
KEA128-C010RE	red
KEA128-C010VT	violet
KEA128-C010BL	black
KEA128-C010BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 128 / 1x1 • height 128 mm fix.

Inner height 125 mm. Without grid divider. Useful for EPi® holes inserts or grid dividers (page 27).

Cat. No.	Colour
KEA128-C1NA	natural
KEA128-C1NY	neon-yellow
KEA128-C1NB	neon-blue
KEA128-C1NG	neon-green
KEA128-C1NP	neon-red/pink
KEA128-C1NO	neon-orange
KEA128-C1WH	white

Cat. No.	Colour
KEA128-C1YE	yellow
KEA128-C1BU	blue
KEA128-C1GR	green
KEA128-C1RE	red
KEA128-C1VT	violet
KEA128-C1BL	black
KEA128-C1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 147 / 1x1 • height 147 mm fix.

Inner height 145 mm. Without grid divider. Useful for EPi® holes inserts or grid dividers (page 27).

Cat. No.	Colour
KEA147-C1NA	natural
KEA147-C1NY	neon-yellow
KEA147-C1NB	neon-blue
KEA147-C1NG	neon-green
KEA147-C1NP	neon-red/pink
KEA147-C1NO	neon-orange
KEA147-C1WH	white

Cat. No.	Colour
KEA147-C1YE	yellow
KEA147-C1BU	blue
KEA147-C1GR	green
KEA147-C1RE	red
KEA147-C1VT	violet
KEA147-C1BL	black
KEA147-C1BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 37

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31




Drawer Racks
page 37

Matching Racks

Vertical Racks
Please request!

Horizontal Racks
Please request!

Drawer Racks
Please request!



EPPi® Grid Dividers

K = Polypropylene (PP)
N = Fiberboard

Cells for vials	Heights					
	15 mm	25 mm	30 mm	40 mm	65 mm	80 mm
Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
16x16 up to ca. 7.2 mm Ø	RB10151616N	RB10251616N	RB10301616N	RB10401616N	RB10651616N	RB10801616N
15x15 up to ca. 7.7 mm Ø	RB10151515N	RB10251515N	RB10301515N	RB10401515N	RB10651515N	RB10801515N
14x14 up to ca. 8.3 mm Ø	RB10151414N	RB10251414N	RB10301414N	RB10401414N	RB10651414N	RB10801414N
13x13 up to ca. 9.0 mm Ø	RB10151313N	RB10251313N	RB10301313N	RB10401313N	RB10651313N	RB10801313N
12x12 up to ca. 9.8 mm Ø	RB10151212N	RB10251212N	RB10301212N	RB10401212N	RB10651212N	RB10801212N
11x11 up to ca. 10.8 mm Ø	RB10151111N	RB10251111N	RB10301111N	RB10401111N	RB10651111N	RB10801111N
10x10 up to ca. 12.0 mm Ø	RB10151010N	RB10251010N	RE10301010 K	RB10401010N	RB10651010N	RB10801010N
9x9 up to ca. 13.6 mm Ø	RB10150909N	RB10250909N	RE10300909 K	RB10400909N	RB10650909N	RB10800909N
8x8 up to ca. 15.0 mm Ø	RB10150808N	RB10250808N	RE10300808 K	RB10400808N	RB10650808N	RB10800808N
7x7 up to ca. 17.3 mm Ø	RB10150707N	RB10250707N	RE10300707 K	RB10400707N	RB10650707N	RB10800707N
6x6 up to ca. 20.3 mm Ø	RB10150606N	RB10250606N	RB10300606N	RB10400606N	RB10650606N	RB10800606N
5x5 up to ca. 24.5 mm Ø	RB10150505N	RB10250505N	RE10300505 K	RB10400505N	RB10650505N	RB10800505N
4x4 up to ca. 30.8 mm Ø	RB10150404N	RB10250404N	RE10300404 K	RB10400404N	RB10650404N	RB10800404N
3x3 up to ca. 41.2 mm Ø	RB10150303N	RB10250303N	RE10300303 K	RB10400303N	RB10650303N	RB10800303N





Kryo Box 100 / 10x10 • height 52 mm fix.

100-place cryogenic storage box made of extra strong polycarbonate.
For 10x10 cryo vials 1.0/2.0 ml with female thread. Transparent lid with numeric ID code (1-100). 133x133x52 mm (WxDxH).

Cat. No.	Colour	Height
KPC50-C100BU	blue	52 mm
KPC50-C100GR	green	52 mm
KPC50-C100RE	red	52 mm
KPC50-C100YE	yellow	52 mm



Kryo Box 81 / 9x9 • height 52 mm fix.

81-place cryogenic storage box made of extra strong polycarbonate.
For 9x9 cryo vials 1.0/2.0 ml. Transparent lid with numeric ID code (1-81). 133x133x52 mm (WxDxH).

Cat. No.	Colour	Height
KPC50-C81BU	blue	52 mm
KPC50-C81GR	green	52 mm
KPC50-C81RE	red	52 mm
KPC50-C81YE	yellow	52 mm



Kryo Box 25 / 5x5 • height 52 mm fix.

25-place cryogenic storage box made of extra strong polycarbonate.
For 5x5 cryo vials 1.0/2.0 ml. Transparent lid with numeric ID code (1-25). 76x76x52 mm (WxDxH).

Cat. No.	Colour	Height
KPC50-C25BU	blue	52 mm
KPC50-C25GR	green	52 mm
KPC50-C25RE	red	52 mm
KPC50-C25YE	yellow	52 mm



Kryo Box 81H / 9x9 • height 95 mm fix.

81-place cryogenic storage box made of extra strong polycarbonate.
For 9x9 cryo vials 3.0 up to 5.0 ml. Transparent lid with numeric ID code (1-81). 133x133x95 mm (WxDxH).

Cat. No.	Colour	Height
KPC95-C81BU	blue	95 mm
KPC95-C81GR	green	95 mm
KPC95-C81RE	red	95 mm
KPC95-C81YE	yellow	95 mm

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 64

Horizontal Racks
page 64


Drawer Racks
Please request!

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 36



EPPI® 32 Vertical Racks

Useful for EPPI® 32 and all boxes up to external dimensions 133x133x32 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE32-08AN	8	282
TE32-09AN	9	317
TE32-10AN	10	351
TE32-11AN	11	386
TE32-12AN	12	420
TE32-13AN	13	455
TE32-14AN	14	489
TE32-15AN	15	524
TE32-16AN	16	558
TE32-17AN	17	593
TE32-18AN	18	627
TE32-19AN	19	662
TE32-20AN	20	696
TE32-21AN	21	731

EPPI® 37 Vertical Racks

Useful for EPPI® 37 and all boxes up to external dimensions 133x133x37 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE37-05AN	5	204
TE37-06AN	6	243
TE37-07AN	7	283
TE37-08AN	8	322
TE37-09AN	9	362
TE37-10AN	10	401
TE37-11AN	11	441
TE37-12AN	12	480
TE37-13AN	13	520
TE37-14AN	14	559
TE37-15AN	15	599
TE37-16AN	16	638
TE37-17AN	17	678
TE37-18AN	18	717

EPPI® 45 + 50 Vertical Racks

Useful for EPPI® 45 + 50 and all boxes up to external dimensions 133x133x53 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE50-04AN	4	226
TE50-05AN	5	281
TE50-06AN	6	336
TE50-07AN	7	391
TE50-08AN	8	446
TE50-09AN	9	501
TE50-10AN	10	556
TE50-11AN	11	611
TE50-12AN	12	666
TE50-13AN	13	721
TE50-14AN	14	776

EPPI® 70 + 75 Vertical Racks

Useful for EPPI® 70 + 75 and all boxes up to external dimensions 133x133x76 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE75-03AN	3	242
TE75-04AN	4	320
TE75-05AN	5	399
TE75-06AN	6	477
TE75-07AN	7	556
TE75-08AN	8	634
TE75-09AN	9	713

EPPI® 80 Vertical Racks

Useful for EPPI® 80 and all boxes up to external dimensions 133x133x80 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE80-03AN	3	252
TE80-04AN	4	334
TE80-05AN	5	416
TE80-06AN	6	498
TE80-07AN	7	580
TE80-08AN	8	662
TE80-09AN	9	744

EPPI® 95 + 96 Vertical Racks

Useful for EPPI® 95 + 96 and all boxes up to external dimensions 133x133x97 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE96-03AN	3	303
TE96-04AN	4	402
TE96-05AN	5	501
TE96-06AN	6	600
TE96-07AN	7	699

EPPI® 102 Vertical Racks

Useful for EPPI® 102 and all boxes up to external dimensions 133x133x103 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE102-03AN	3	321
TE102-04AN	4	426
TE102-05AN	5	531
TE102-06AN	6	636
TE102-07AN	7	741

EPPI® 105 Vertical Racks

Useful for EPPI® 105 and all boxes up to external dimensions 133x133x105 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE105-02AN	2	220
TE105-03AN	3	327
TE105-04AN	4	434
TE105-05AN	5	541
TE105-06AN	6	648

EPPI® 122 + 128 Vertical Racks

Useful for EPPI® 122 + 128 and all boxes up to external dimensions 133x133x128 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE128-02AN	2	267
TE128-03AN	3	398
TE128-04AN	4	528
TE128-05AN	5	659



Folding handle positioned crosswise as standard.

Alternatives: Folding handle positioned lengthwise or hand grip. Please enquire if needed.



All vertical racks can be supplied with a safety rod if required.

This increases the depth by 5 mm to 143 mm. Please enquire if needed.


EPPI® Horizontal Racks

Width 139 mm. Stainless steel.



EPPI® Horizontal Racks without intermediate shelf

For EPPI® boxes of different heights. Also useful for boxes with external dimensions of up to 133 x 133 mm (WxD).



Standard Heights


Cat. No.	Bays	Inner Height	Dimensions
	Depth x Height	Bay (mm)	D x H (mm)
SEX169-2-01AN	2 x 1	165	278x169
SEX169-3-01AN	3 x 1	165	414x169
SEX169-4-01AN	4 x 1	165	549x169
SEX169-5-01AN	5 x 1	165	685x169
SEX224-2-01AN	2 x 1	220	278x224
SEX224-3-01AN	3 x 1	220	414x224
SEX224-4-01AN	4 x 1	220	549x224
SEX224-5-01AN	5 x 1	220	685x224
SEX250-2-01AN	2 x 1	246	278x250
SEX250-3-01AN	3 x 1	246	414x250
SEX250-4-01AN	4 x 1	246	549x250
SEX250-5-01AN	5 x 1	246	685x250

Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
	Depth x Height	Bay (mm)	D x H (mm)
SEX279-2-01AN	2 x 1	275	278x279
SEX279-3-01AN	3 x 1	275	414x279
SEX279-4-01AN	4 x 1	275	549x279
SEX279-5-01AN	5 x 1	275	685x279
SEX334-2-01AN	2 x 1	330	278x334
SEX334-3-01AN	3 x 1	330	414x334
SEX334-4-01AN	4 x 1	330	549x334
SEX334-5-01AN	5 x 1	330	685x334
SEX389-2-01AN	2 x 1	385	278x389
SEX389-3-01AN	3 x 1	385	414x389
SEX389-4-01AN	4 x 1	385	549x389
SEX389-5-01AN	5 x 1	385	685x389

EPPI® Horizontal Racks with one intermediate shelf

For EPPI® boxes of different heights. Also useful for boxes with external dimensions of up to 133 x 133 mm (WxD).



Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
	Depth x Height	Bay (mm)	D x H (mm)
SEX169-2-02AN	2x2	82	278x169
SEX169-3-02AN	3x2	82	414x169
SEX169-4-02AN	4x2	82	549x169
SEX169-5-02AN	5x2	82	685x169
SEX224-2-02AN	2x2	109	278x224
SEX224-3-02AN	3x2	109	414x224
SEX224-4-02AN	4x2	109	549x224
SEX224-5-02AN	5x2	109	685x224
SEX250-2-02AN	2x2	122	278x250
SEX250-3-02AN	3x2	122	414x250
SEX250-4-02AN	4x2	122	549x250
SEX250-5-02AN	5x2	122	685x250

Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
	Depth x Height	Bay (mm)	D x H (mm)
SEX279-2-02AN	2x2	137	278x279
SEX279-3-02AN	3x2	137	414x279
SEX279-4-02AN	4x2	137	549x279
SEX279-5-02AN	5x2	137	685x279
SEX334-2-02AN	2x2	164	278x334
SEX334-3-02AN	3x2	164	414x334
SEX334-4-02AN	4x2	164	549x334
SEX334-5-02AN	5x2	164	685x334
SEX389-2-02AN	2x2	192	278x389
SEX389-3-02AN	3x2	192	414x389
SEX389-4-02AN	4x2	192	549x389
SEX389-5-02AN	5x2	192	685x389



Folding handle positioned vertically as standard.
Alternatives: Folding handle positioned horizontally or hand grip. Please enquire if needed.



All horizontal racks can be supplied with safety rods if required.
Add 5 mm to width and height in each case.
Please enquire if needed.



EPPI® 32 Horizontal Racks

Useful for EPPI® 32 and all boxes up to external dimensions 133x133x32 mm (WxDxH).

Cat. No.	Bays Depth x Height	Dimensions D x H (mm)
SE32-37AN	3 x 7	414x245
SE32-38AN	3 x 8	414x280
SE32-39AN	3 x 9	414x314
SE32-45AN	4 x 5	549x176
SE32-46AN	4 x 6	549x211
SE32-47AN	4 x 7	549x245
SE32-48AN	4 x 8	549x280
SE32-49AN	4 x 9	549x314
SE32-55AN	5 x 5	685x176
SE32-56AN	5 x 6	685x211
SE32-57AN	5 x 7	685x245
SE32-58AN	5 x 8	685x280
SE32-59AN	5 x 9	685x314



EPPI® 37 Horizontal Racks

Useful for EPPI® 37 and all boxes up to external dimensions 133x133x37 mm (WxDxH).

Cat. No.	Bays Depth x Height	Dimensions D x H (mm)
SE37-36AN	3 x 6	414x241
SE37-37AN	3 x 7	414x280
SE37-38AN	3 x 8	414x320
SE37-39AN	3 x 9	414x359
SE37-45AN	4 x 5	549x201
SE37-46AN	4 x 6	549x241
SE37-47AN	4 x 7	549x280
SE37-48AN	4 x 8	549x320
SE37-49AN	4 x 9	549x359
SE37-55AN	5 x 5	685x201
SE37-56AN	5 x 6	685x241
SE37-57AN	5 x 7	685x280
SE37-58AN	5 x 8	685x320
SE37-59AN	5 x 9	685x359

EPPI® 45 + 50 Horizontal Racks

Useful for EPPI® 45 + 50 and all boxes up to external dimensions 133x133x53 mm (WxDxH).

Cat. No.	Bays Depth x Height	Dimensions D x H (mm)
SE50-24AN	2 x 4	278x224
SE50-25AN	2 x 5	278x279
SE50-26AN	2 x 6	278x334
SE50-33AN	3 x 3	414x169
SE50-34AN	3 x 4	414x224
SE50-35AN	3 x 5	414x279
SE50-36AN	3 x 6	414x334
SE50-37AN	3 x 7	414x389
SE50-43AN	4 x 3	549x169
SE50-44AN	4 x 4	549x224
SE50-45AN	4 x 5	549x279
SE50-46AN	4 x 6	549x334
SE50-47AN	4 x 7	549x389
SE50-53AN	5 x 3	685x169
SE50-54AN	5 x 4	685x224
SE50-55AN	5 x 5	685x279
SE50-56AN	5 x 6	685x334
SE50-57AN	5 x 7	685x389



Folding handle positioned vertically as standard.
Alternatives: Folding handle positioned horizontally or hand grip. Please enquire if needed.



EPPI® 70 + 75 Horizontal Racks

Useful for EPPI® 70 + 75 and all boxes up to external dimensions 133x133x76 mm (WxDxH).

Cat. No.	Bays Depth x Height	Dimensions D x H (mm)
SE75-32AN	3 x 2	414x161
SE75-33AN	3 x 3	414x239
SE75-34AN	3 x 4	414x318
SE75-42AN	4 x 2	549x161
SE75-43AN	4 x 3	549x239
SE75-44AN	4 x 4	549x318
SE75-52AN	5 x 2	685x161
SE75-53AN	5 x 3	685x239
SE75-54AN	5 x 4	685x318

EPPI® 80 Horizontal Racks

Useful for EPPI® 80 and all boxes up to external dimensions 133x133x80 mm (WxDxH).

Cat. No.	Bays Depth x Height	Dimensions D x H (mm)
SE80-32AN	3 x 2	414x168
SE80-33AN	3 x 3	414x250
SE80-34AN	3 x 4	414x332
SE80-42AN	4 x 2	549x168
SE80-43AN	4 x 3	549x250
SE80-44AN	4 x 4	549x332
SE80-52AN	5 x 2	685x168
SE80-53AN	5 x 3	685x250
SE80-54AN	5 x 4	685x332

EPPI® 95 + 96 Horizontal Racks

Useful for EPPI® 95 + 96 and all boxes up to external dimensions 133x133x97 mm (WxDxH).

Cat. No.	Bays Depth x Height	Dimensions D x H (mm)
SE96-32AN	3 x 2	414x202
SE96-33AN	3 x 3	414x301
SE96-42AN	4 x 2	549x202
SE96-43AN	4 x 3	549x301
SE96-52AN	5 x 2	685x202
SE96-53AN	5 x 3	685x301



EPPI® 102 Horizontal Racks

Useful for EPPI® 102 and all boxes up to external dimensions 133x133x103 mm (WxDxH).

Cat. No.	Bays Depth x Height	Dimensions D x H (mm)
SE102-32AN	3 x 2	414x214
SE102-33AN	3 x 3	414x319
SE102-42AN	4 x 2	549x214
SE102-43AN	4 x 3	549x319
SE102-52AN	5 x 2	685x214
SE102-53AN	5 x 3	685x319

EPPI® 105 Horizontal Racks

Useful for EPPI® 105 and all boxes up to external dimensions 133x133x105 mm (WxDxH).

Cat. No.	Bays Depth x Height	Dimensions D x H (mm)
SE105-32AN	3 x 2	414x218
SE105-33AN	3 x 3	414x325
SE105-42AN	4 x 2	549x218
SE105-43AN	4 x 3	549x325
SE105-52AN	5 x 2	685x218
SE105-53AN	5 x 3	685x325

EPPI® 122 + 128 Horizontal Racks

Useful for EPPI® 122 + 128 and all boxes up to external dimensions 133x133x128 mm (WxDxH).

Cat. No.	Bays Depth x Height	Dimensions D x H (mm)
SE128-32AN	3 x 2	414x265
SE128-33AN	3 x 3	414x395
SE128-42AN	4 x 2	549x265
SE128-43AN	4 x 3	549x395
SE128-52AN	5 x 2	685x265
SE128-53AN	5 x 3	685x395



All horizontal racks can be supplied with safety rods if required. Add 5 mm to width and height in each case. Please enquire if needed.

EPi® Drawer Racks

Width 140 mm. Stainless steel.


With Folding Handle Open Design



Please add to cat. no.:

- L** Base of drawer in middle open.
Facilitates removal and insertion.
- V** Base of drawer closed

With Folding Handle Closed Design




Please add to cat. no.:

- L** Base of drawer in middle open.
Facilitates removal and insertion.
- V** Base of drawer closed

EPi® 45 + 50 Drawer Racks. With Folding Handle

Useful for EPi® 45 + 50 and all boxes with external dimensions of up to 133x133x53 mm (WxDxH). Front closed.

With Safety Stop				
Cat. No.	L	Number	DxH	
	V	Depth	Height	(mm)
SE50-33NAM	■	3	3	434 x 184
SE50-34NAM	■	3	4	434 x 245
SE50-35NAM	■	3	5	434 x 306
SE50-36NAM	■	3	6	434 x 367
SE50-43NAM	■	4	3	567 x 184
SE50-44NAM	■	4	4	567 x 245
SE50-45NAM	■	4	5	567 x 306
SE50-46NAM	■	4	6	567 x 367
SE50-53NAM	■	5	3	700 x 184
SE50-54NAM	■	5	4	700 x 245
SE50-55NAM	■	5	5	700 x 306
SE50-56NAM	■	5	6	700 x 367



Without Safety Stop				
Cat. No.	L	Number	DxH	
	V	Depth	Height	(mm)
SE50-33NAO	■	3	3	416 x 173
SE50-34NAO	■	3	4	416 x 231
SE50-35NAO	■	3	5	416 x 288
SE50-36NAO	■	3	6	416 x 346
SE50-43NAO	■	4	3	549 x 173
SE50-44NAO	■	4	4	549 x 231
SE50-45NAO	■	4	5	549 x 288
SE50-46NAO	■	4	6	549 x 346
SE50-53NAO	■	5	3	682 x 173
SE50-54NAO	■	5	4	682 x 231
SE50-55NAO	■	5	5	682 x 288
SE50-56NAO	■	5	6	682 x 346



EPi® 45 + 50 Drawer Racks. With Folding Handle

Useful for EPi® 45 + 50 and all boxes with external dimensions of up to 133x133x53 mm (WxDxH). Front closed.

With Safety Stop				
Cat. No.	L	Number	DxH	
	V	Depth	Height	(mm)
SE50-33GAM	■	3	3	444 x 184
SE50-34GAM	■	3	4	444 x 245
SE50-35GAM	■	3	5	444 x 306
SE50-36GAM	■	3	6	444 x 367
SE50-43GAM	■	4	3	577 x 184
SE50-44GAM	■	4	4	577 x 245
SE50-45GAM	■	4	5	577 x 306
SE50-46GAM	■	4	6	577 x 367
SE50-53GAM	■	5	3	710 x 184
SE50-54GAM	■	5	4	710 x 245
SE50-55GAM	■	5	5	710 x 306
SE50-56GAM	■	5	6	710 x 367



Without Safety Stop				
Cat. No.	L	Number	DxH	
	V	Depth	Height	(mm)
SE50-33GAO	■	3	3	426 x 173
SE50-34GAO	■	3	4	426 x 231
SE50-35GAO	■	3	5	426 x 288
SE50-36GAO	■	3	6	426 x 346
SE50-43GAO	■	4	3	559 x 173
SE50-44GAO	■	4	4	559 x 231
SE50-45GAO	■	4	5	559 x 288
SE50-46GAO	■	4	6	559 x 346
SE50-53GAO	■	5	3	692 x 173
SE50-54GAO	■	5	4	692 x 231
SE50-55GAO	■	5	5	692 x 288
SE50-56GAO	■	5	6	692 x 346




With Hand Grip Open Design



With Hand Grip Closed Design



Without Handle (Tray) Without Safety Stop



EPi® 45 + 50 Drawer Racks. With Hand Grip

Useful for EPi® 45 + 50 and all boxes with external dimensions of up to 133x133x53 mm (WxDxH). Base of drawers closed.

With Safety Stop				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SE50-33NBMV	3	3	434 x 184	
SE50-34NBMV	3	4	434 x 245	
SE50-35NBMV	3	5	434 x 306	
SE50-36NBMV	3	6	434 x 367	
SE50-43NBMV	4	3	567 x 184	
SE50-44NBMV	4	4	567 x 245	
SE50-45NBMV	4	5	567 x 306	
SE50-46NBMV	4	6	567 x 367	
SE50-53NBMV	5	3	700 x 184	
SE50-54NBMV	5	4	700 x 245	
SE50-55NBMV	5	5	700 x 306	
SE50-56NBMV	5	6	700 x 367	

Without Safety Stop				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SE50-33NB0V	3	3	434 x 173	
SE50-34NB0V	3	4	434 x 231	
SE50-35NB0V	3	5	434 x 288	
SE50-36NB0V	3	6	434 x 346	
SE50-43NB0V	4	3	567 x 173	
SE50-44NB0V	4	4	567 x 231	
SE50-45NB0V	4	5	567 x 288	
SE50-46NB0V	4	6	567 x 346	
SE50-53NB0V	5	3	700 x 173	
SE50-54NB0V	5	4	700 x 231	
SE50-55NB0V	5	5	700 x 288	
SE50-56NB0V	5	6	700 x 346	

EPi® 45 + 50 Drawer Racks. With Hand Grip

Useful for EPi® 45 + 50 and all boxes with external dimensions of up to 133x133x53 mm (WxDxH). Base of drawers closed.

With Safety Stop				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SE50-33GBMV	3	3	444 x 184	
SE50-34GBMV	3	4	444 x 245	
SE50-35GBMV	3	5	444 x 306	
SE50-36GBMV	3	6	444 x 367	
SE50-43GBMV	4	3	577 x 184	
SE50-44GBMV	4	4	577 x 245	
SE50-45GBMV	4	5	577 x 306	
SE50-46GBMV	4	6	577 x 367	
SE50-53GBMV	5	3	710 x 184	
SE50-54GBMV	5	4	710 x 245	
SE50-55GBMV	5	5	710 x 306	
SE50-56GBMV	5	6	710 x 367	

Without Safety Stop				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SE50-33GB0V	3	3	444 x 173	
SE50-34GB0V	3	4	444 x 231	
SE50-35GB0V	3	5	444 x 288	
SE50-36GB0V	3	6	444 x 346	
SE50-43GB0V	4	3	577 x 173	
SE50-44GB0V	4	4	577 x 231	
SE50-45GB0V	4	5	577 x 288	
SE50-46GB0V	4	6	577 x 346	
SE50-53GB0V	5	3	710 x 173	
SE50-54GB0V	5	4	710 x 231	
SE50-55GB0V	5	5	710 x 288	
SE50-56GB0V	5	6	710 x 346	

EPi® 45 + 50 Drawer Racks. Without Handle


Useful for EPi® 45 + 50 and all boxes with external dimensions of up to 133x133x53 mm (WxDxH). Base of drawers closed.

Open Design				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SE50-33NT0V	3	3	416 x 173	
SE50-34NT0V	3	4	416 x 231	
SE50-35NT0V	3	5	416 x 288	
SE50-36NT0V	3	6	416 x 346	
SE50-43NT0V	4	3	549 x 173	
SE50-44NT0V	4	4	549 x 231	
SE50-45NT0V	4	5	549 x 288	
SE50-46NT0V	4	6	549 x 346	
SE50-53NT0V	5	3	682 x 173	
SE50-54NT0V	5	4	682 x 231	
SE50-55NT0V	5	5	682 x 288	
SE50-56NT0V	5	6	682 x 346	

Closed Design				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SE50-33GT0V	3	3	426 x 173	
SE50-34GT0V	3	4	426 x 231	
SE50-35GT0V	3	5	426 x 288	
SE50-36GT0V	3	6	426 x 346	
SE50-43GT0V	4	3	559 x 173	
SE50-44GT0V	4	4	559 x 231	
SE50-45GT0V	4	5	559 x 288	
SE50-46GT0V	4	6	559 x 346	
SE50-53GT0V	5	3	692 x 173	
SE50-54GT0V	5	4	692 x 231	
SE50-55GT0V	5	5	692 x 288	
SE50-56GT0V	5	6	692 x 346	

EPPI® 32 Drawer Racks

With Grip Rail Open Design




Example: Cat. No. **SE32-38NCMV**
For 24 x EPPI® 32 Boxes.

With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x434x325 mm

See illustrations on page 32

EPPI® 37 Drawer Racks

With Grip Rail Open Design




Example: Cat. No. **SE37-46NCML**
For 24 x EPPI® 37 Boxes.

With safety stop. Base of drawers open.
Dimensions: WxDxH = 140x567x274 mm

See illustrations on page 32

EPPI® 32 Drawer Racks

Without Handle (Tray) Open Design




Example: Cat. No. **SE32-28NTOV**
For 16 x EPPI® 32 Boxes.

Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x283x297 mm

EPPI® 37 Drawer Racks


Without Handle (Tray) Open Design



Example: Cat. No. **SE37-35NTOV**
For 15 x EPPI® 37 Boxes.

Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x416x211 mm

With Grip Rail




Without Handle (Tray)



EPPI® 75 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SE75-33NAML
For 9 x EPPI® 75 Boxes.
With safety stop. Base of drawers open.
Dimensions: WxDxH = 140x434x254 mm

See illustrations on page 32

EPPI® 75 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SE75-22NBOV
For 4 x EPPI® 75 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x301x163 mm

EPPI® 75 Drawer Racks


Without Handle (Tray) Open Design¹



Example: Cat. No. SE75-43NTOV
For 12 x EPPI® 75 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x549x244 mm

EPPI® 80 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SE80-42NAOL
For 8 x EPPI® 80 Boxes.
Without safety stop. Base of drawers open.
Dimensions: WxDxH = 140x549x170 mm

See illustrations on page 32

EPPI® 80 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SE80-53NBMV
For 15 x EPPI® 80 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x700x265 mm

EPPI® 80 Drawer Racks

Without Handle (Tray) Open Design¹




Example: Cat. No. SE80-33NTOV
For 9 x EPPI® 80 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x416x254 mm



¹ All drawer racks can also be supplied with a closed design. (See illustrations at page 32+33).
Please enquire if needed.

EPPI® 96 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SE96-43NAOL
For 12 x EPPI® 96 Boxes.
Without safety stop. Base of drawers open.
Dimensions: WxDxH = 140x549x305 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

EPPI® 96 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SE96-33NBMV
For 9 x EPPI® 96 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x434x316 mm

EPPI® 96 Drawer Racks


Without Handle (Tray) Open Design¹



Example: Cat. No. SE96-23NTOV
For 6 x EPPI® 96 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x283x305 mm

EPPI® 102 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SE102-32NAML
For 6 x EPPI® 102 Boxes.
With safety stop. Base of drawers open.
Dimensions: WxDxH = 140x434x223 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

EPPI® 102 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SE102-52NBVO
For 10 x EPPI® 102 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x700x216 mm

EPPI® 102 Drawer Racks.

Without Handle (Tray) Open Design¹




Example: Cat. No. SE102-42NTOV
For 8 x EPPI® 102 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x549x216 mm



- ¹ All drawer racks can also be supplied with a closed design. (See illustrations at page 32+33). Please enquire if needed.
- ² All drawer racks can be produced with additional drawers for shallower boxes (e.g. EPPI® 50) to exploit the interior space of uprights to the maximum.

EPPI® 105 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SE105-43NAML
For 12 x EPPI® 105 Boxes.
With safety stop. Base of drawers open.
Dimensions: WxDxH = 140x567x340 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

EPPI® 105 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SE105-32NBMV
For 6 x EPPI® 105 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x434x227 mm

EPPI® 105 Drawer Racks


Without Handle (Tray) Open Design¹



Example: Cat. No. SE105-53NTOV
For 15 x EPPI® 105 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x682x329 mm

EPPI® 128 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SE128-22NAOL
For 4 x EPPI® 128 Boxes.
Without safety stop. Base of drawers open.
Dimensions: WxDxH = 140x283x267 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

EPPI® 128 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SE128-42NBMV
For 8 x EPPI® 128 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x567x274 mm

EPPI® 128 Drawer Racks.

Without Handle (Tray) Open Design¹



Example: Cat. No. SE128-32NTOV
For 6 x EPPI® 128 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x416x267 mm



- ¹ All drawer racks can also be supplied with a closed design. (See illustrations at page 32+33). Please enquire if needed.
- ² All drawer racks can be produced with additional drawers for shallower boxes (e.g. EPPI® 50) to exploit the interior space of uprights to the maximum.

EPPI® Upright Bins


Width 139 and 278 mm. Stainless steel.



EPPI® Upright Bins

Single Width

With same external dimensions as horizontal racks. Open at top. With hand grip. Ideal for combining with horizontal racks and drawer racks.



Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SEK169-121 BN	1	168	139x278x169
SEK169-131 BN	1	168	139x414x169
SEK169-141 BN	1	168	139x549x169
SEK169-151 BN	1	168	139x685x169
SEK224-121 BN	1	223	139x278x224
SEK224-131 BN	1	223	139x414x224
SEK224-141 BN	1	223	139x549x224
SEK224-151 BN	1	223	139x685x224
SEK250-121 BN	1	249	139x278x250
SEK250-131 BN	1	249	139x414x250
SEK250-141 BN	1	249	139x549x250
SEK250-151 BN	1	249	139x685x250


Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SEK279-121 BN	1	278	139x278x279
SEK279-131 BN	1	278	139x414x279
SEK279-141 BN	1	278	139x549x279
SEK279-151 BN	1	278	139x685x279
SEK334-121 BN	1	333	139x278x334
SEK334-131 BN	1	333	139x414x334
SEK334-141 BN	1	333	139x549x334
SEK334-151 BN	1	333	139x685x334
SEK389-121 BN	1	388	139x278x389
SEK389-131 BN	1	388	139x414x389
SEK389-141 BN	1	388	139x549x389
SEK389-151 BN	1	388	139x685x389

EPPI® Upright Bins

Double Width

With same external dimensions as horizontal racks. Open at top. With hand grip. Ideal for combining with horizontal racks and drawer racks.



Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SEK169-221 BN	1	168	278x278x169
SEK169-231 BN	1	168	278x414x169
SEK169-241 BN	1	168	278x549x169
SEK169-251 BN	1	168	278x685x169
SEK224-221 BN	1	223	278x278x224
SEK224-231 BN	1	223	278x414x224
SEK224-241 BN	1	223	278x549x224
SEK224-251 BN	1	223	278x685x224
SEK250-221 BN	1	249	278x278x250
SEK250-231 BN	1	249	278x414x250
SEK250-241 BN	1	249	278x549x250
SEK250-251 BN	1	249	278x685x250

Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SEK279-221 BN	1	278	278x278x279
SEK279-231 BN	1	278	278x414x279
SEK279-241 BN	1	278	278x549x279
SEK279-251 BN	1	278	278x685x279
SEK334-221 BN	1	333	278x278x334
SEK334-231 BN	1	333	278x414x334
SEK334-241 BN	1	333	278x549x334
SEK334-251 BN	1	333	278x685x334
SEK389-221 BN	1	388	278x278x389
SEK389-231 BN	1	388	278x414x389
SEK389-241 BN	1	388	278x549x389
SEK389-251 BN	1	388	278x685x389

ObiBox made of Polystyrene Slide Boxes



Obi Box 25 Slide Box for 25 Slides

With index card and numbering in lower section.
For slides 76x26 mm. Made of polystyrene.
96x122x35 mm (WxDxH).


Cat. No.	Colour
KOS35-25BU	blue
KOS35-25GR	green
KOS35-25RE	red
KOS35-25YE	yellow

Matching Racks

Vertical Racks
Please request!

Horizontal Racks
Please request!

Drawer Racks
Please request!



Obi Box 50 Slide Box for 50 Slides

With index card and numbering in lower section.
For slides 76x26 mm. Made of polystyrene.
100x230x35 mm (WxDxH).

Cat. No.	Colour
KOS35-50BU	blue
KOS35-50GR	green
KOS35-50RE	red
KOS35-50YE	yellow

Matching Racks

Vertical Racks
Please request!

Horizontal Racks
Please request!

Drawer Racks
Please request!



Obi Box 100 Slide Box for 100 Slides

With index card and numbering in lower section.
For slides 76x26 mm. Made of polystyrene.
185x230x35 mm (WxDxH).

Cat. No.	Colour
KOS35-100BU	blue
KOS35-100GR	green
KOS35-100RE	red
KOS35-100YE	yellow

Matching Racks

Vertical Racks
Please request!

Horizontal Racks
Please request!

Drawer Racks
Please request!



EPPi® Ident Identification Labels down to -90°C

Self-adhesive labels for marking boxes made of plastic and fibreboard as well as stainless steel racks, etc. Paper, white, printable.

Supplied on sheets in A4 format in packs of 50 sheets each.

Cat. No.	W x H	Quantity
LAL2813WH50	28x13 mm	50x154
LAL5230WH50	52x30 mm	50x36
LAL7015WH50	70x15 mm	50x57
LAL7020WH50	70x20 mm	50x42
LAL7035WH50	70x35 mm	50x24
LAL6020WH50	60x20 mm	50x42

Cat. No.	W x H	Quantity
LAL6006WH50	60x 6 mm	50x135
LAL4013WH50	40x13 mm	50x115
LAL4025WH50	40x25 mm	50x60
LAL4812WH50	48x12 mm	50x96
LAL3512WH50	35x12 mm	50x144
LAL010WH50	Ø10 mm	50x100



ALPHA Box 136 x 136 Cryo and Storage Boxes


Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 136 x 136 mm. Height = inner height.

Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.

Colour			Heights						
1 or 2			25 mm	32 mm	50 mm	75 mm	100 mm	110 mm	130 mm
	WH	1 standard	025	032	050	075	100	110	130
		2 special							
	BU	1 standard	025	032	050	075	100	110	130
		2 special							
	YE	1 standard	025	032	050	075	100	110	130
		2 special							
	GR	1 standard	025	032	050	075	100	110	130
		2 special							
	RE	1 standard	025	032	050	075	100	110	130
		2 special							

How to determine the item no.:



Example 50 mm height, white, standard:

KA1050WH1



Grid Dividers for ALPHA Box

			Heights					
			15 mm	25 mm	30 mm	40 mm	65 mm	80 mm
Cells for vials			Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
16x16	up to	ca. 7.5 mm Ø	RA10151616N	RA10251616N	RA10301616N	RA10401616N	RA10651616N	RA10801616N
15x15	up to	ca. 8.0 mm Ø	RA10151515N	RA10251515N	RA10301515N	RA10401515N	RA10651515N	RA10801515N
14x14	up to	ca. 8.7 mm Ø	RA10151414N	RA10251414N	RA10301414N	RA10401414N	RA10651414N	RA10801414N
13x13	up to	ca. 9.4 mm Ø	RA10151313N	RA10251313N	RA10301313N	RA10401313N	RA10651313N	RA10801313N
12x12	up to	ca. 10.2 mm Ø	RA10151212N	RA10251212N	RA10301212N	RA10401212N	RA10651212N	RA10801212N
11x11	up to	ca. 11.2 mm Ø	RA10151111N	RA10251111N	RA10301111N	RA10401111N	RA10651111N	RA10801111N
10x10	up to	ca. 12.5 mm Ø	RA10151010N	RA10251010N	RA10301010N	RA10401010N	RA10651010N	RA10801010N
9x9	up to	ca. 13.9 mm Ø	RA10150909N	RA10250909N	RA10300909N	RA10400909N	RA10650909N	RA10800909N
8x8	up to	ca. 15.7 mm Ø	RA10150808N	RA10250808N	RA10300808N	RA10400808N	RA10650808N	RA10800808N
7x7	up to	ca. 17.9 mm Ø	RA10150707N	RA10250707N	RA10300707N	RA10400707N	RA10650707N	RA10800707N
6x6	up to	ca. 21.0 mm Ø	RA10150606N	RA10250606N	RA10300606N	RA10400606N	RA10650606N	RA10800606N
5x5	up to	ca. 25.5 mm Ø	RA10150505N	RA10250505N	RA10300505N	RA10400505N	RA10650505N	RA10800505N
4x4	up to	ca. 32.0 mm Ø	RA10150404N	RA10250404N	RA10300404N	RA10400404N	RA10650404N	RA10800404N
3x3	up to	ca. 42.9 mm Ø	RA10150303N	RA10250303N	RA10300303N	RA10400303N	RA10650303N	RA10800303N

When used with item number given in the overview (beginning with "RA1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: ALPHA Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 10x10 cells.

KA1050WH1 + 0301010N = Cat. No. KA1050WH10301010N



ALPHA 25 Vertical Racks

Useful for ALPHA 25 and all boxes up to external dimensions 136x136x28 mm (WxDxH). 2 boxes in each bay. Number in ().

Cat. No.	Bays	Height (mm)
TA25-08AN2	4 (8)	238
TA25-10AN2	5 (10)	296
TA25-12AN2	6 (12)	354
TA25-14AN2	7 (14)	412
TA25-16AN2	8 (16)	470
TA25-18AN2	9 (18)	528
TA25-20AN2	10 (20)	586
TA25-22AN2	11 (22)	644
TA25-24AN2	12 (24)	702
TA25-26AN2	13 (26)	760

ALPHA 32 Vertical Racks

Useful for ALPHA 32 and all boxes up to external dimensions 136x136x35 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TA32-06AN	6	228
TA32-07AN	7	265
TA32-08AN	8	302
TA32-09AN	9	339
TA32-10AN	10	376
TA32-11AN	11	413
TA32-12AN	12	450
TA32-13AN	13	487
TA32-14AN	14	524
TA32-15AN	15	561
TA32-16AN	16	598
TA32-17AN	17	635
TA32-18AN	18	672
TA32-19AN	19	709
TA32-20AN	20	746

ALPHA 50 Vertical Racks

Useful for ALPHA 50 and all boxes up to external dimensions 136x136x53 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TA50-04AN	4	226
TA50-05AN	5	281
TA50-06AN	6	336
TA50-07AN	7	391
TA50-08AN	8	446
TA50-09AN	9	501
TA50-10AN	10	556
TA50-11AN	11	611
TA50-12AN	12	666
TA50-13AN	13	721
TA50-14AN	14	776

ALPHA 75 Vertical Racks

Useful for ALPHA 75 and all boxes up to external dimensions 136x136x78 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TA75-03AN	3	246
TA75-04AN	4	326
TA75-05AN	5	406
TA75-06AN	6	486
TA75-07AN	7	566
TA75-08AN	8	646
TA75-09AN	9	726

ALPHA 100 Vertical Racks

Useful for ALPHA 100 and all boxes up to external dimensions 136x136x103 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TA100-03AN	3	321
TA100-04AN	4	426
TA100-05AN	5	531
TA100-06AN	6	636
TA100-07AN	7	741

ALPHA 110 Vertical Racks

Useful for ALPHA 110 and all boxes up to external dimensions 136x136x113 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TA110-03AN	3	351
TA110-04AN	4	466
TA110-05AN	5	581
TA110-06AN	6	696
TA110-07AN	7	811

ALPHA 130 Vertical Racks

Useful for ALPHA 130 and all boxes up to external dimensions 136x136x133 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TA130-03AN	3	411
TA130-04AN	4	546
TA130-05AN	5	681



Folding handle positioned crosswise as standard.

Alternatives: Folding handle positioned lengthwise or hand grip. Please enquire if needed.



All vertical racks can be supplied with a safety rod if required.

This increases the depth by 5 mm to 146 mm. Please enquire if needed.


ALPHA Horizontal Racks

Width 142 mm. Stainless steel.



ALPHA Horizontal Racks without intermediate shelf

For ALPHA boxes of different heights. Also useful for boxes with external dimensions of up to 136 x 136 mm (WxD).



Standard Heights


Cat. No.	Bays	Inner Height	Dimensions
	Depth x Height	Bay (mm)	D x H (mm)
SAX169-2-01AN	2 x 1	165	286x169
SAX169-3-01AN	3 x 1	165	426x169
SAX169-4-01AN	4 x 1	165	565x169
SAX169-5-01AN	5 x 1	165	705x169
SAX224-2-01AN	2 x 1	220	286x224
SAX224-3-01AN	3 x 1	220	426x224
SAX224-4-01AN	4 x 1	220	565x224
SAX224-5-01AN	5 x 1	220	705x224
SAX250-2-01AN	2 x 1	246	286x250
SAX250-3-01AN	3 x 1	246	426x250
SAX250-4-01AN	4 x 1	246	565x250
SAX250-5-01AN	5 x 1	246	705x250

Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
	Depth x Height	Bay (mm)	D x H (mm)
SAX279-2-01AN	2 x 1	275	286x279
SAX279-3-01AN	3 x 1	275	426x279
SAX279-4-01AN	4 x 1	275	565x279
SAX279-5-01AN	5 x 1	275	705x279
SAX334-2-01AN	2 x 1	330	286x334
SAX334-3-01AN	3 x 1	330	426x334
SAX334-4-01AN	4 x 1	330	565x334
SAX334-5-01AN	5 x 1	330	705x334
SAX389-2-01AN	2 x 1	385	286x389
SAX389-3-01AN	3 x 1	385	426x389
SAX389-4-01AN	4 x 1	385	565x389
SAX389-5-01AN	5 x 1	385	705x389

ALPHA Horizontal Racks with one intermediate shelf

For ALPHA boxes of different heights. Also useful for boxes with external dimensions of up to 136 x 136 mm (WxD).



Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
	Depth x Height	Bay (mm)	D x H (mm)
SAX169-2-02AN	2x2	82	286x169
SAX169-3-02AN	3x2	82	426x169
SAX169-4-02AN	4x2	82	565x169
SAX169-5-02AN	5x2	82	705x169
SAX224-2-02AN	2x2	109	286x224
SAX224-3-02AN	3x2	109	426x224
SAX224-4-02AN	4x2	109	565x224
SAX224-5-02AN	5x2	109	705x224
SAX250-2-02AN	2x2	122	286x250
SAX250-3-02AN	3x2	122	426x250
SAX250-4-02AN	4x2	122	565x250
SAX250-5-02AN	5x2	122	705x250

Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
	Depth x Height	Bay (mm)	D x H (mm)
SAX279-2-02AN	2x2	137	286x279
SAX279-3-02AN	3x2	137	426x279
SAX279-4-02AN	4x2	137	565x279
SAX279-5-02AN	5x2	137	705x279
SAX334-2-02AN	2x2	164	286x334
SAX334-3-02AN	3x2	164	426x334
SAX334-4-02AN	4x2	164	565x334
SAX334-5-02AN	5x2	164	705x334
SAX389-2-02AN	2x2	192	286x389
SAX389-3-02AN	3x2	192	426x389
SAX389-4-02AN	4x2	192	565x389
SAX389-5-02AN	5x2	192	705x389



Folding handle positioned vertically as standard.
Alternatives: Folding handle positioned horizontally or hand grip. Please enquire if needed.



All horizontal racks can be supplied with safety rods if required.
Add 5 mm to width and height in each case.
Please enquire if needed.



ALPHA 25 Horizontal Racks

Useful for ALPHA 25 and all boxes up to external dimensions 136x136x28 mm (WxDxH). 2 boxes in each bay. Number in ().

Cat. No.	Bays		Dimensions DxH (mm)
	Depth x Height		
SA25-2x08AN2	2 x 4 (16)		286x236
SA25-2x10AN2	2 x 5 (20)		286x294
SA25-2x12AN2	2 x 6 (24)		286x352
SA25-3x08AN2	3 x 4 (24)		426x236
SA25-3x10AN2	3 x 5 (30)		426x294
SA25-3x12AN2	3 x 6 (36)		426x352
SA25-4x08AN2	4 x 4 (32)		565x236
SA25-4x10AN2	4 x 5 (40)		565x294
SA25-4x12AN2	4 x 6 (48)		565x352
SA25-5x08AN2	5 x 4 (40)		705x236
SA25-5x10AN2	5 x 5 (50)		705x294
SA25-5x12AN2	5 x 6 (60)		705x352



Folding handle positioned vertically as standard.

Alternatives: Folding handle positioned horizontally or hand grip. Please enquire if needed.



ALPHA 75 Horizontal Racks

Useful for ALPHA 75 and all boxes up to external dimensions 136x136x78 mm (WxDxH).

Cat. No.	Bays		Dimensions DxH (mm)
	Depth x Height		
SA75-32AN	3 x 2		426x164
SA75-33AN	3 x 3		426x244
SA75-34AN	3 x 4		426x324
SA75-42AN	4 x 2		565x164
SA75-43AN	4 x 3		565x244
SA75-44AN	4 x 4		565x324
SA75-52AN	5 x 2		705x164
SA75-53AN	5 x 3		705x244
SA75-54AN	5 x 4		705x324



ALPHA 130 Horizontal Racks

Useful for ALPHA 130 and all boxes up to external dimensions 136x136x133 mm (WxDxH).

Cat. No.	Bays		Dimensions DxH (mm)
	Depth x Height		
SA130-32AN	3 x 2		426x279
SA130-33AN	3 x 3		426x409
SA130-42AN	4 x 2		565x279
SA130-43AN	4 x 3		565x409
SA130-52AN	5 x 2		705x279
SA130-53AN	5 x 3		705x409



ALPHA 32 Horizontal Racks

Useful for ALPHA 32 and all boxes up to external dimensions 136x136x35 mm (WxDxH).

Cat. No.	Bays		Dimensions DxH (mm)
	Depth x Height		
SA32-35AN	3 x 5		426x189
SA32-36AN	3 x 6		426x226
SA32-37AN	3 x 7		426x263
SA32-38AN	3 x 8		426x300
SA32-39AN	3 x 9		426x337
SA32-45AN	4 x 5		565x189
SA32-46AN	4 x 6		565x226
SA32-47AN	4 x 7		565x263
SA32-48AN	4 x 8		565x300
SA32-49AN	4 x 9		565x337
SA32-55AN	5 x 5		705x189
SA32-56AN	5 x 6		705x226
SA32-57AN	5 x 7		705x263
SA32-58AN	5 x 8		705x300
SA32-59AN	5 x 9		705x337

ALPHA 100 Horizontal Racks

Useful for ALPHA 100 and all boxes up to external dimensions 136x136x103 mm (WxDxH).

Cat. No.	Bays		Dimensions DxH (mm)
	Depth x Height		
SA100-32AN	3 x 2		426x214
SA100-33AN	3 x 3		426x319
SA100-42AN	4 x 2		565x214
SA100-43AN	4 x 3		565x319
SA100-52AN	5 x 2		705x214
SA100-53AN	5 x 3		705x319



All horizontal racks can be supplied with safety rods if required. Add 5 mm to width and height in each case. Please enquire if needed.



ALPHA 50 Horizontal Racks

Useful for ALPHA 50 and all boxes up to external dimensions 136x136x53 mm (WxDxH).

Cat. No.	Bays		Dimensions DxH (mm)
	Depth x Height		
SA50-24AN	2 x 4		286x224
SA50-25AN	2 x 5		286x279
SA50-26AN	2 x 6		286x334
SA50-33AN	3 x 3		426x169
SA50-34AN	3 x 4		426x224
SA50-35AN	3 x 5		426x279
SA50-36AN	3 x 6		426x334
SA50-37AN	3 x 7		426x389
SA50-43AN	4 x 3		565x169
SA50-44AN	4 x 4		565x224
SA50-45AN	4 x 5		565x279
SA50-46AN	4 x 6		565x334
SA50-47AN	4 x 7		565x389
SA50-53AN	5 x 3		705x169
SA50-54AN	5 x 4		705x224
SA50-55AN	5 x 5		705x279
SA50-56AN	5 x 6		705x334
SA50-57AN	5 x 7		705x389

ALPHA 110 Horizontal Racks


Useful for ALPHA 110 and all boxes up to external dimensions 136x136x113 mm (WxDxH).

Cat. No.	Bays		Dimensions DxH (mm)
	Depth x Height		
SA110-32AN	3 x 2		426x234
SA110-33AN	3 x 3		426x349
SA110-42AN	4 x 2		565x234
SA110-43AN	4 x 3		565x349
SA110-52AN	5 x 2		705x234
SA110-53AN	5 x 3		705x349

ALPHA Drawer Racks

Width 145 mm. Stainless steel.


With Folding Handle Open Design



Please add to cat. no.:

- L** Base of drawer in middle open.
Facilitates removal and insertion.
- V** Base of drawer closed

With Folding Handle Closed Design





Please add to cat. no.:

- L** Base of drawer in middle open.
Facilitates removal and insertion.
- V** Base of drawer closed

ALPHA 50 Drawer Racks. With Folding Handle

Useful for ALPHA 50 and all boxes with external dimensions of up to 136x136x53 mm (WxDxH). Front closed.



With Safety Stop					Without Safety Stop				
Cat. No.	L	Number	DxH	V	Cat. No.	L	Number	DxH	V
				Depth Height (mm)					Depth Height (mm)
SA50-33NAM	■	3	3	444 x 184	SA50-33NAO	■	3	3	429 x 173
SA50-34NAM	■	3	4	444 x 245	SA50-34NAO	■	3	4	429 x 231
SA50-35NAM	■	3	5	444 x 306	SA50-35NAO	■	3	5	429 x 288
SA50-36NAM	■	3	6	444 x 367	SA50-36NAO	■	3	6	429 x 346
SA50-43NAM	■	4	3	580 x 184	SA50-43NAO	■	4	3	565 x 173
SA50-44NAM	■	4	4	580 x 245	SA50-44NAO	■	4	4	565 x 231
SA50-45NAM	■	4	5	580 x 306	SA50-45NAO	■	4	5	565 x 288
SA50-46NAM	■	4	6	580 x 367	SA50-46NAO	■	4	6	565 x 346
SA50-53NAM	■	5	3	716 x 184	SA50-53NAO	■	5	3	701 x 173
SA50-54NAM	■	5	4	716 x 245	SA50-54NAO	■	5	4	701 x 231
SA50-55NAM	■	5	5	716 x 306	SA50-55NAO	■	5	5	701 x 288
SA50-56NAM	■	5	6	716 x 367	SA50-56NAO	■	5	6	701 x 346



ALPHA 50 Drawer Racks. With Folding Handle

Useful for ALPHA 50 and all boxes with external dimensions of up to 136x136x53 mm (WxDxH). Front closed.

With Safety Stop					Without Safety Stop				
Cat. No.	L	Number	DxH	V	Cat. No.	L	Number	DxH	V
				Depth Height (mm)					Depth Height (mm)
SA50-33GAM	■	3	3	454 x 184	SA50-33GAO	■	3	3	439 x 173
SA50-34GAM	■	3	4	454 x 245	SA50-34GAO	■	3	4	439 x 231
SA50-35GAM	■	3	5	454 x 306	SA50-35GAO	■	3	5	439 x 288
SA50-36GAM	■	3	6	454 x 367	SA50-36GAO	■	3	6	439 x 346
SA50-43GAM	■	4	3	590 x 184	SA50-43GAO	■	4	3	575 x 173
SA50-44GAM	■	4	4	590 x 245	SA50-44GAO	■	4	4	575 x 231
SA50-45GAM	■	4	5	590 x 306	SA50-45GAO	■	4	5	575 x 288
SA50-46GAM	■	4	6	590 x 367	SA50-46GAO	■	4	6	575 x 346
SA50-53GAM	■	5	3	726 x 184	SA50-53GAO	■	5	3	711 x 173
SA50-54GAM	■	5	4	726 x 245	SA50-54GAO	■	5	4	711 x 231
SA50-55GAM	■	5	5	726 x 306	SA50-55GAO	■	5	5	711 x 288
SA50-56GAM	■	5	6	726 x 367	SA50-56GAO	■	5	6	711 x 346




ALPHA Drawer Racks

Width 145 mm. Stainless steel.


With Hand Grip Open Design



With Hand Grip Closed Design



Without Handle (Tray) Without Safety Stop



ALPHA 50 Drawer Racks. With Hand Grip

Useful for ALPHA 50 and all boxes with external dimensions of up to 136x136x53 mm (WxDxH). Base of drawers closed.

With Safety Stop				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SA50-33NBMV	3	3	444 x 184	
SA50-34NBMV	3	4	444 x 245	
SA50-35NBMV	3	5	444 x 306	
SA50-36NBMV	3	6	444 x 367	
SA50-43NBMV	4	3	580 x 184	
SA50-44NBMV	4	4	580 x 245	
SA50-45NBMV	4	5	580 x 306	
SA50-46NBMV	4	6	580 x 367	
SA50-53NBMV	5	3	716 x 184	
SA50-54NBMV	5	4	716 x 245	
SA50-55NBMV	5	5	716 x 306	
SA50-56NBMV	5	6	716 x 367	

Without Safety Stop				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SA50-33NBOV	3	3	444 x 173	
SA50-34NBOV	3	4	444 x 231	
SA50-35NBOV	3	5	444 x 288	
SA50-36NBOV	3	6	444 x 346	
SA50-43NBOV	4	3	580 x 173	
SA50-44NBOV	4	4	580 x 231	
SA50-45NBOV	4	5	580 x 288	
SA50-46NBOV	4	6	580 x 346	
SA50-53NBOV	5	3	716 x 173	
SA50-54NBOV	5	4	716 x 231	
SA50-55NBOV	5	5	716 x 288	
SA50-56NBOV	5	6	716 x 346	

ALPHA 50 Drawer Racks. With Hand Grip

Useful for ALPHA 50 and all boxes with external dimensions of up to 136x136x53 mm (WxDxH). Base of drawers closed.

With Safety Stop				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SA50-33GBMV	3	3	454 x 184	
SA50-34GBMV	3	4	454 x 245	
SA50-35GBMV	3	5	454 x 306	
SA50-36GBMV	3	6	454 x 367	
SA50-43GBMV	4	3	590 x 184	
SA50-44GBMV	4	4	590 x 245	
SA50-45GBMV	4	5	590 x 306	
SA50-46GBMV	4	6	590 x 367	
SA50-53GBMV	5	3	726 x 184	
SA50-54GBMV	5	4	726 x 245	
SA50-55GBMV	5	5	726 x 306	
SA50-56GBMV	5	6	726 x 367	

Without Safety Stop				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SA50-33GBOV	3	3	454 x 173	
SA50-34GBOV	3	4	454 x 231	
SA50-35GBOV	3	5	454 x 288	
SA50-36GBOV	3	6	454 x 346	
SA50-43GBOV	4	3	590 x 173	
SA50-44GBOV	4	4	590 x 231	
SA50-45GBOV	4	5	590 x 288	
SA50-46GBOV	4	6	590 x 346	
SA50-53GBOV	5	3	726 x 173	
SA50-54GBOV	5	4	726 x 231	
SA50-55GBOV	5	5	726 x 288	
SA50-56GBOV	5	6	726 x 346	

ALPHA 50 Drawer Racks. Without Handle


Useful for ALPHA 50 and all boxes with external dimensions of up to 136x136x53 mm (WxDxH). Base of drawers closed.

Open Design				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SA50-33NTOV	3	3	429 x 173	
SA50-34NTOV	3	4	429 x 231	
SA50-35NTOV	3	5	429 x 288	
SA50-36NTOV	3	6	429 x 346	
SA50-43NTOV	4	3	565 x 173	
SA50-44NTOV	4	4	565 x 231	
SA50-45NTOV	4	5	565 x 288	
SA50-46NTOV	4	6	565 x 346	
SA50-53NTOV	5	3	701 x 173	
SA50-54NTOV	5	4	701 x 231	
SA50-55NTOV	5	5	701 x 288	
SA50-56NTOV	5	6	701 x 346	

Closed Design				
Cat. No.	Number		DxH	
	Depth	Height	(mm)	
SA50-33GTOV	3	3	439 x 173	
SA50-34GTOV	3	4	439 x 231	
SA50-35GTOV	3	5	439 x 288	
SA50-36GTOV	3	6	439 x 346	
SA50-43GTOV	4	3	575 x 173	
SA50-44GTOV	4	4	575 x 231	
SA50-45GTOV	4	5	575 x 288	
SA50-46GTOV	4	6	575 x 346	
SA50-53GTOV	5	3	711 x 173	
SA50-54GTOV	5	4	711 x 231	
SA50-55GTOV	5	5	711 x 288	
SA50-56GTOV	5	6	711 x 346	

ALPHA 25 Drawer Racks

With Grip Rail Open Design




Example: Cat. No. SA25-47NCMV
For 28 x ALPHA 25 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x580x253 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 44

ALPHA 32 Drawer Racks

With Grip Rail Open Design




Example: Cat. No. SA32-37NCML
For 21 x ALPHA 32 Boxes.
With safety stop. Base of drawers open.
Dimensions: WxDxH = 145x444x302 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 44

ALPHA 25 Drawer Racks


Without Handle (Tray) Open Design



Example: Cat. No. SA25-56NTOV
For 30 x ALPHA 25 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x701x196 mm


ALPHA 32 Drawer Racks

Without Handle (Tray) Open Design



Example: Cat. No. SA32-25NTOV
For 10 x ALPHA 32 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x293x198 mm

With Grip Rail




Without Handle (Tray)



ALPHA 75 Drawer Racks

With Folding Handle Open Design¹



Example: Cat. No. **SA75-42NAML**
For 8 x ALPHA 75 Boxes.


With safety stop. Base of drawers open.
Dimensions: WxDxH = 145x580x173 mm

- L** Base of drawer in middle open. Facilitates removal and insertion.
- V** Base of drawer closed.

See illustrations on page 44

ALPHA 75 Drawer Racks

**With Hand Grip
Open Design¹**




Example: Cat. No. **SA75-53NBOV**
For 15 x ALPHA 75 Boxes.

Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x716x248 mm

ALPHA 75 Drawer Racks

Without Handle (Tray)




Example: Cat. No. **SA75-33NTOV**
For 9 x ALPHA 75 Boxes.

Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x429x248 mm

ALPHA 100 Drawer Racks

With Folding Handle Open Design¹



Example: Cat. No. **SA100-22NAOL**
For 4 x ALPHA 100 Boxes.


Without safety stop. Base of drawers open.
Dimensions: WxDxH = 145x293x216 mm

- L** Base of drawer in middle open. Facilitates removal and insertion.
- V** Base of drawer closed.

See illustrations on page 44

ALPHA 100 Drawer Racks

**With Hand Grip
Open Design¹**




Example: Cat. No. **SA100-32NAMV**
For 6 x ALPHA 100 Boxes.

With safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x444x223 mm

ALPHA 100 Drawer Racks

Without Handle (Tray)

Open Design¹



Example: Cat. No. **SA100-42NTOV**
For 8 x ALPHA 100 Boxes.

Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x565x216 mm




¹ All drawer racks can also be supplied with a closed design. (See illustrations at page 44+45). Please enquire if needed.

2 All drawer racks can be produced with additional drawers for shallower boxes (e.g. ALPHA 50) to exploit the interior space of uprights to the maximum.

ALPHA 110 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SA110-52NAML
For 10 x ALPHA 110 Boxes.
With safety stop. Base of drawers open.
Dimensions: WxDxH = 145x716x243 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 44

ALPHA 110 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SA110-33NBMV
For 9 x ALPHA 110 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x444x364 mm

ALPHA 110 Drawer Racks


Without Handle (Tray) Open Design¹



Example: Cat. No. SA110-43NTOV
For 12 x ALPHA 110 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x565x353 mm

ALPHA 130 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SA130-22NAOL
For 4 x ALPHA 130 Boxes.
Without safety stop. Base of drawers open.
Dimensions: WxDxH = 145x293x276 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 44

ALPHA 130 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SA130-42NBMV
For 8 x ALPHA 130 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x580x283 mm

ALPHA 130 Drawer Racks

Without Handle (Tray) Open Design¹



Example: Cat. No. SA130-52NTOV
For 10 x ALPHA 130 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 145x701x276 mm



- ¹ All drawer racks can also be supplied with a closed design. (See illustrations at page 44+45). Please enquire if needed.
- ² All drawer racks can be produced with additional drawers for shallower boxes (e.g. ALPHA 50) to exploit the interior space of uprights to the maximum.

ALPHA Upright Bins


Width 142 und 284 mm. Stainless steel.



ALPHA Upright Bins

Single Width

With same external dimensions as horizontal racks. Open at top. With hand grip. Ideal for combining with horizontal racks and drawer racks.



Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SAK169-121BN	1	168	142x286x169
SAK169-131BN	1	168	142x426x169
SAK169-141BN	1	168	142x565x169
SAK169-151BN	1	168	142x705x169
SAK224-121BN	1	223	142x286x224
SAK224-131BN	1	223	142x426x224
SAK224-141BN	1	223	142x565x224
SAK224-151BN	1	223	142x705x224
SAK250-121BN	1	249	142x286x250
SAK250-131BN	1	249	142x426x250
SAK250-141BN	1	249	142x565x250
SAK250-151BN	1	249	142x705x250


Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SAK279-121BN	1	278	142x286x279
SAK279-131BN	1	278	142x426x279
SAK279-141BN	1	278	142x565x279
SAK279-151BN	1	278	142x705x279
SAK334-121BN	1	333	142x286x334
SAK334-131BN	1	333	142x426x334
SAK334-141BN	1	333	142x565x334
SAK334-151BN	1	333	142x705x334
SAK389-121BN	1	388	142x286x389
SAK389-131BN	1	388	142x426x389
SAK389-141BN	1	388	142x565x389
SAK389-151BN	1	388	142x705x389

ALPHA Upright Bins

Double Width

With same external dimensions as horizontal racks. Open at top. With hand grip. Ideal for combining with horizontal racks and drawer racks.



Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SAK169-221BN	1	168	284x286x169
SAK169-231BN	1	168	284x426x169
SAK169-241BN	1	168	284x565x169
SAK169-251BN	1	168	284x705x169
SAK224-221BN	1	223	284x286x224
SAK224-231BN	1	223	284x426x224
SAK224-241BN	1	223	284x565x224
SAK224-251BN	1	223	284x705x224
SAK250-221BN	1	249	284x286x250
SAK250-231BN	1	249	284x426x250
SAK250-241BN	1	249	284x565x250
SAK250-251BN	1	249	284x705x250

Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SAK279-221BN	1	278	284x286x279
SAK279-231BN	1	278	284x426x279
SAK279-241BN	1	278	284x565x279
SAK279-251BN	1	278	284x705x279
SAK334-221BN	1	333	284x286x334
SAK334-231BN	1	333	284x426x334
SAK334-241BN	1	333	284x565x334
SAK334-251BN	1	333	284x705x334
SAK389-221BN	1	388	284x286x389
SAK389-231BN	1	388	284x426x389
SAK389-241BN	1	388	284x565x389
SAK389-251BN	1	388	284x705x389



ALPHA Box long2 136 x 268 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 136 x 268 mm. Height = inner height.


Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.


Matching Racks

Please request!

			Heights						
Colour	1 or 2		25 mm	32 mm	50 mm	75 mm	100 mm	110 mm	130 mm
 WH	1 standard 2 special		025	032	050	075	100	110	130
 BU	1 standard 2 special		025	032	050	075	100	110	130
 YE	1 standard 2 special		025	032	050	075	100	110	130
 GR	1 standard 2 special		025	032	050	075	100	110	130
 RE	1 standard 2 special		025	032	050	075	100	110	130



How to determine
the item no.:



Example 50 mm height,
white, standard:

KA 2 0 50 WH 1

Grid Dividers for ALPHA Box long2

			Heights					
			15 mm	25 mm	30 mm	40 mm	65 mm	80 mm
Cells for vials			Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
16x32 up to ca. 7.5 mm Ø			RA20151632N	RA20251632N	RA20301632N	RA20401632N	RA20651632N	RA20801632N
15x30 up to ca. 8.0 mm Ø			RA20151530N	RA20251530N	RA20301530N	RA20401530N	RA20651530N	RA20801530N
14x28 up to ca. 8.7 mm Ø			RA20151428N	RA20251428N	RA20301428N	RA20401428N	RA20651428N	RA20801428N
13x26 up to ca. 9.4 mm Ø			RA20151326N	RA20251326N	RA20301326N	RA20401326N	RA20651326N	RA20801326N
12x24 up to ca. 10.2 mm Ø			RA20151224N	RA20251224N	RA20301224N	RA20401224N	RA20651224N	RA20801224N
11x22 up to ca. 11.2 mm Ø			RA20151122N	RA20251122N	RA20301122N	RA20401122N	RA20651122N	RA20801122N
10x20 up to ca. 12.5 mm Ø			RA20151020N	RA20251020N	RA20301020N	RA20401020N	RA20651020N	RA20801020N
9x18 up to ca. 13.9 mm Ø			RA20150918N	RA20250918N	RA20300918N	RA20400918N	RA20650918N	RA20800918N
8x16 up to ca. 15.7 mm Ø			RA20150816N	RA20250816N	RA20300816N	RA20400816N	RA20650816N	RA20800816N
7x14 up to ca. 17.9 mm Ø			RA20150714N	RA20250714N	RA20300714N	RA20400714N	RA20650714N	RA20800714N
6x12 up to ca. 21.0 mm Ø			RA20150612N	RA20250612N	RA20300612N	RA20400612N	RA20650612N	RA20800612N
5x10 up to ca. 25.5 mm Ø			RA20150510N	RA20250510N	RA20300510N	RA20400510N	RA20650510N	RA20800510N
4x8 up to ca. 32.0 mm Ø			RA20150408N	RA20250408N	RA20300408N	RA20400408N	RA20650408N	RA20800408N
3x6 up to ca. 42.9 mm Ø			RA20150306N	RA20250306N	RA20300306N	RA20400306N	RA20650306N	RA20800306N

When used with item number given in the overview (beginning with "RA2..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: ALPHA Box long2 of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 10x20 cells.

KA2050WH1 + 0301020N = Cat. No. KA2050WH10301020N



ALPHA Box short 136 x 94 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 136 x 94 mm. Height = inner height.

Other heights can be produced on request.


Grid dividers can be ordered separately or can be supplied fitted.

Matching Racks

Please request!

			Heights						
Colour	1 or 2		25 mm	32 mm	50 mm	75 mm	100 mm	110 mm	130 mm
 WH	1 standard 2 special		025	032	050	075	100	110	130
 BU	1 standard 2 special		025	032	050	075	100	110	130
 YE	1 standard 2 special		025	032	050	075	100	110	130
 GR	1 standard 2 special		025	032	050	075	100	110	130
 RE	1 standard 2 special		025	032	050	075	100	110	130

How to determine
the item no.:



Example 50 mm height,
white, standard:

KAK050WH1

Grid Dividers for ALPHA Box short

			Heights					
			15 mm	25 mm	30 mm	40 mm	65 mm	80 mm
Cells for vials			Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
16x11 up to ca. 7.5 mm Ø			RAK0151611N	RAK0251611N	RAK0301611N	RAK0401611N	RAK0651611N	RAK0801611N
15x10 up to ca. 8.0 mm Ø			RAK0151510N	RAK0251510N	RAK0301510N	RAK0401510N	RAK0651510N	RAK0801510N
14x9 up to ca. 8.7 mm Ø			RAK0151409N	RAK0251409N	RAK0301409N	RAK0401409N	RAK0651409N	RAK0801409N
13x8 up to ca. 9.4 mm Ø			RAK0151308N	RAK0251308N	RAK0301308N	RAK0401308N	RAK0651308N	RAK0801308N
12x8 up to ca. 10.2 mm Ø			RAK0151208N	RAK0251208N	RAK0301208N	RAK0401208N	RAK0651208N	RAK0801208N
11x7 up to ca. 11.2 mm Ø			RAK0151107N	RAK0251107N	RAK0301107N	RAK0401107N	RAK0651107N	RAK0801107N
10x6 up to ca. 12.5 mm Ø			RAK0151006N	RAK0251006N	RAK0301006N	RAK0401006N	RAK0651006N	RAK0801006N
9x6 up to ca. 13.9 mm Ø			RAK0150906N	RAK0250906N	RAK0300906N	RAK0400906N	RAK0650906N	RAK0800906N
8x5 up to ca. 15.7 mm Ø			RAK0150805N	RAK0250805N	RAK0300805N	RAK0400805N	RAK0650805N	RAK0800805N
7x4 up to ca. 17.9 mm Ø			RAK0150704N	RAK0250704N	RAK0300704N	RAK0400704N	RAK0650704N	RAK0800704N
6x3 up to ca. 21.0 mm Ø			RAK0150603N	RAK0250603N	RAK0300603N	RAK0400603N	RAK0650603N	RAK0800603N
5x3 up to ca. 25.5 mm Ø			RAK0150503N	RAK0250503N	RAK0300503N	RAK0400503N	RAK0650503N	RAK0800503N
4x2 up to ca. 32.0 mm Ø			RAK0150402N	RAK0250402N	RAK0300402N	RAK0400402N	RAK0650402N	RAK0800402N
3x2 up to ca. 42.9 mm Ø			RAK0150302N	RAK0250302N	RAK0300302N	RAK0400302N	RAK0650302N	RAK0800302N

When used with item number given in the overview (beginning with "RAK..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: ALPHA Box short of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 10x6 cells.

KAK050WH1 + 0301006N = Cat. No. KAK050WH10301006N



BRAVO Box 133 x 133 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 133 x 133 mm. Height = inner height.


Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.

			Heights						
Colour	1 or 2		25 mm	32 mm	50 mm*	75 mm	100 mm	110 mm	130 mm
white WH	1 standard 2 special		025	032	050	075	100	110	130
blue BU	1 standard 2 special		025	032	050	075	100	110	130
yellow YE	1 standard 2 special		025	032	050	075	100	110	130
green GR	1 standard 2 special		025	032	050	075	100	110	130
red RE	1 standard 2 special		025	032	050	075	100	110	130

* Produced traditionally to the American design.
(Actual internal height: 45 mm)

How to determine
the item no.:



Example 50 mm height,
white, standard:

KB 1 050 WH 1



Grid Dividers for BRAVO Box

			Heights					
			15 mm	25 mm	30 mm	40 mm	65 mm	80 mm
Cells for vials			Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
16x16 up to ca. 7.2 mm Ø			RB10151616N	RB10251616N	RB10301616N	RB10401616N	RB10651616N	RB10801616N
15x15 up to ca. 7.7 mm Ø			RB10151515N	RB10251515N	RB10301515N	RB10401515N	RB10651515N	RB10801515N
14x14 up to ca. 8.3 mm Ø			RB10151414N	RB10251414N	RB10301414N	RB10401414N	RB10651414N	RB10801414N
13x13 up to ca. 9.0 mm Ø			RB10151313N	RB10251313N	RB10301313N	RB10401313N	RB10651313N	RB10801313N
12x12 up to ca. 9.8 mm Ø			RB10151212N	RB10251212N	RB10301212N	RB10401212N	RB10651212N	RB10801212N
11x11 up to ca. 10.8 mm Ø			RB10151111N	RB10251111N	RB10301111N	RB10401111N	RB10651111N	RB10801111N
10x10 up to ca. 12.0 mm Ø			RB10151010N	RB10251010N	RB10301010N	RB10401010N	RB10651010N	RB10801010N
9x9 up to ca. 13.6 mm Ø			RB10150909N	RB10250909N	RB10300909N	RB10400909N	RB10650909N	RB10800909N
8x8 up to ca. 15.0 mm Ø			RB10150808N	RB10250808N	RB10300808N	RB10400808N	RB10650808N	RB10800808N
7x7 up to ca. 17.3 mm Ø			RB10150707N	RB10250707N	RB10300707N	RB10400707N	RB10650707N	RB10800707N
6x6 up to ca. 20.3 mm Ø			RB10150606N	RB10250606N	RB10300606N	RB10400606N	RB10650606N	RB10800606N
5x5 up to ca. 24.5 mm Ø			RB10150505N	RB10250505N	RB10300505N	RB10400505N	RB10650505N	RB10800505N
4x4 up to ca. 30.8 mm Ø			RB10150404N	RB10250404N	RB10300404N	RB10400404N	RB10650404N	RB10800404N
3x3 up to ca. 41.2 mm Ø			RB10150303N	RB10250303N	RB10300303N	RB10400303N	RB10650303N	RB10800303N

When used with item number given in the overview (beginning with "RB1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: BRAVO Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 9x9 cells.

KB1050WH1 + 0300909N = Cat. No. KB1050WH10300909N



BRAVO 25 Vertical Racks

Useful for BRAVO 25 and all boxes up to external dimensions 133x133x28 mm (WxDxH). 2 boxes in each bay. Number in ().

Cat. No.	Bays	Height (mm)
TB25-08AN2	4 (8)	238
TB25-10AN2	5 (10)	296
TB25-12AN2	6 (12)	354
TB25-14AN2	7 (14)	412
TB25-16AN2	8 (16)	470
TB25-18AN2	9 (18)	528
TB25-20AN2	10 (20)	586
TB25-22AN2	11 (22)	644
TB25-24AN2	12 (24)	702
TB25-26AN2	13 (26)	760

BRAVO 32 Vertical Racks

Useful for BRAVO 32 and all boxes up to external dimensions 133x133x35 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TB32-06AN	6	228
TB32-07AN	7	265
TB32-08AN	8	302
TB32-09AN	9	339
TB32-10AN	10	376
TB32-11AN	11	413
TB32-12AN	12	450
TB32-13AN	13	487
TB32-14AN	14	524
TB32-15AN	15	561
TB32-16AN	16	598
TB32-17AN	17	635
TB32-18AN	18	672
TB32-19AN	19	709
TB32-20AN	20	746

BRAVO 50 Vertical Racks

Useful for BRAVO 50 and all boxes up to external dimensions 133x133x53 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE50-04AN	4	226
TE50-05AN	5	281
TE50-06AN	6	336
TE50-07AN	7	391
TE50-08AN	8	446
TE50-09AN	9	501
TE50-10AN	10	556
TE50-11AN	11	611
TE50-12AN	12	666
TE50-13AN	13	721
TE50-14AN	14	776

BRAVO 75 Vertical Racks

Useful for BRAVO 75 and all boxes up to external dimensions 133x133x78 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TB75-03AN	3	246
TB75-04AN	4	326
TB75-05AN	5	406
TB75-06AN	6	486
TB75-07AN	7	566
TB75-08AN	8	646
TB75-09AN	9	726

BRAVO 100 Vertical Racks

Useful for BRAVO 100 and all boxes up to external dimensions 133x133x103 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TE102-03AN	3	321
TE102-04AN	4	426
TE102-05AN	5	531
TE102-06AN	6	636
TE102-07AN	7	741

BRAVO 110 Vertical Racks

Useful for BRAVO 110 and all boxes up to external dimensions 133x133x113 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TB110-03AN	3	351
TB110-04AN	4	466
TB110-05AN	5	581
TB110-06AN	6	696
TB110-07AN	7	811

BRAVO 130 Vertical Racks

Useful for BRAVO 130 and all boxes up to external dimensions 133x133x133 mm (WxDxH).

Cat. No.	Bays	Height (mm)
TB130-03AN	3	411
TB130-04AN	4	546
TB130-05AN	5	681



Folding handle positioned crosswise as standard.

Alternatives: Folding handle positioned lengthwise or hand grip. Please enquire if needed.



All vertical racks can be supplied with a safety rod if required.

This increases the depth by 5 mm to 143 mm. Please enquire if needed.



BRAVO 25 Horizontal Racks

Useful for BRAVO 25 and all boxes up to external dimensions 133x133x28 mm (WxDxH). 2 boxes in each bay. Number in ().

Cat. No.	Bays	Dimensions Depth x Height	Dimensions D x H (mm)
SB25-2x08AN2	2 x 4 (16)	278x236	
SB25-2x10AN2	2 x 5 (20)	278x294	
SB25-2x12AN2	2 x 6 (24)	278x352	
SB25-3x08AN2	3 x 4 (24)	414x236	
SB25-3x10AN2	3 x 5 (30)	414x294	
SB25-3x12AN2	3 x 6 (36)	414x352	
SB25-4x08AN2	4 x 4 (32)	549x236	
SB25-4x10AN2	4 x 5 (40)	549x294	
SB25-4x12AN2	4 x 6 (48)	549x352	
SB25-5x08AN2	5 x 4 (40)	685x236	
SB25-5x10AN2	5 x 5 (50)	685x294	
SB25-5x12AN2	5 x 6 (60)	685x352	



Folding handle positioned vertically as standard.

Alternatives: Folding handle positioned horizontally or hand grip. Please enquire if needed.



BRAVO 32 Horizontal Racks

Useful for BRAVO 32 and all boxes up to external dimensions 133x133x35 mm (WxDxH).

Cat. No.	Bays	Dimensions Depth x Height	Dimensions D x H (mm)
SB32-35AN	3 x 5	414x189	
SB32-36AN	3 x 6	414x226	
SB32-37AN	3 x 7	414x263	
SB32-38AN	3 x 8	414x300	
SB32-39AN	3 x 9	414x337	
SB32-45AN	4 x 5	549x189	
SB32-46AN	4 x 6	549x226	
SB32-47AN	4 x 7	549x263	
SB32-48AN	4 x 8	549x300	
SB32-49AN	4 x 9	549x337	
SB32-55AN	5 x 5	685x189	
SB32-56AN	5 x 6	685x226	
SB32-57AN	5 x 7	685x263	
SB32-58AN	5 x 8	685x300	
SB32-59AN	5 x 9	685x337	

BRAVO 75 Horizontal Racks

Useful for BRAVO 75 and all boxes up to external dimensions 133x133x78 mm (WxDxH).

Cat. No.	Bays	Dimensions Depth x Height	Dimensions D x H (mm)
SB75-32AN	3 x 2	414x164	
SB75-33AN	3 x 3	414x244	
SB75-34AN	3 x 4	414x324	
SB75-42AN	4 x 2	549x164	
SB75-43AN	4 x 3	549x244	
SB75-44AN	4 x 4	549x324	
SB75-52AN	5 x 2	685x164	
SB75-53AN	5 x 3	685x244	
SB75-54AN	5 x 4	685x324	



BRAVO 100 Horizontal Racks

Useful for BRAVO 100 and all boxes up to external dimensions 133x133x103 mm (WxDxH).

Cat. No.	Bays	Dimensions Depth x Height	Dimensions D x H (mm)
SE102-32AN	3 x 2	414x214	
SE102-33AN	3 x 3	414x319	
SE102-42AN	4 x 2	549x214	
SE102-43AN	4 x 3	549x319	
SE102-52AN	5 x 2	685x214	
SE102-53AN	5 x 3	685x319	



All horizontal racks can be supplied with safety rods if required. Add 5 mm to width and height in each case. Please enquire if needed.

BRAVO 130 Horizontal Racks

Useful for BRAVO 130 and all boxes up to external dimensions 133x133x133 mm (WxDxH).

Cat. No.	Bays	Dimensions Depth x Height	Dimensions D x H (mm)
SB130-32AN	3 x 2	414x279	
SB130-33AN	3 x 3	414x409	
SB130-42AN	4 x 2	549x279	
SB130-43AN	4 x 3	549x409	
SB130-52AN	5 x 2	685x279	
SB130-53AN	5 x 3	685x409	



BRAVO 50 Horizontal Racks

Useful for BRAVO 50 and all boxes up to external dimensions 133x133x53 mm (WxDxH).

Cat. No.	Bays	Dimensions Depth x Height	Dimensions D x H (mm)
SE50-24AN	2 x 4	278x224	
SE50-25AN	2 x 5	278x279	
SE50-26AN	2 x 6	278x334	
SE50-33AN	3 x 3	414x169	
SE50-34AN	3 x 4	414x224	
SE50-35AN	3 x 5	414x279	
SE50-36AN	3 x 6	414x334	
SE50-37AN	3 x 7	414x389	
SE50-43AN	4 x 3	549x169	
SE50-44AN	4 x 4	549x224	
SE50-45AN	4 x 5	549x279	
SE50-46AN	4 x 6	549x334	
SE50-47AN	4 x 7	549x389	
SE50-53AN	5 x 3	685x169	
SE50-54AN	5 x 4	685x224	
SE50-55AN	5 x 5	685x279	
SE50-56AN	5 x 6	685x334	
SE50-57AN	5 x 7	685x389	

BRAVO 110 Horizontal Racks

Useful for BRAVO 110 and all boxes up to external dimensions 133x133x113 mm (WxDxH).


Cat. No.	Bays	Dimensions Depth x Height	Dimensions D x H (mm)
SB110-32AN	3 x 2	414x234	
SB110-33AN	3 x 3	414x349	
SB110-42AN	4 x 2	549x234	
SB110-43AN	4 x 3	549x349	
SB110-52AN	5 x 2	685x234	
SB110-53AN	5 x 3	685x349	



Horizontal Racks without and with one intermediate shelf see page 30.

BRAVO 25 Drawer Racks

With Grip Rail Open Design




Example: Cat. No. **SB25-38NCMV**
For 24 x BRAVO 25 Boxes.

With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x434x325 mm

BRAVO 32 Drawer Racks

With Grip Rail Open Design




Example: Cat. No. **SB32-46NCML**
For 24 x BRAVO 32 Boxes.

With safety stop. Base of drawers open.
Dimensions: WxDxH = 140x567x274 mm

BRAVO 25 Drawer Racks

Without Handle (Tray) Open Design




Example: Cat. No. **SB25-28NTOV**
For 16 x BRAVO 25 Boxes.

Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x283x297 mm

BRAVO 32 Drawer Racks


Without Handle (Tray) Open Design



Example: Cat. No. **SB32-35NTOV**
For 15 x BRAVO 32 Boxes.

Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x416x211 mm

With Grip Rail





Without Handle (Tray)



BRAVO 50 Drawer Racks




Please use the EPPi® 50 Drawer Racks on pages 32+33 for the BRAVO 50 Box.



¹ All drawer racks can also be supplied with a closed design. (See illustrations on page 32+33). Please enquire if needed.


BRAVO 75 Drawer Racks

With Folding Handle Open Design¹




BRAVO 75 Drawer Racks

With Hand Grip Open Design¹




BRAVO 75 Drawer Racks

Without Handle (Tray) Open Design¹



BRAVO 110 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SB110-43NAML
For 12 x BRAVO 110 Boxes.
With safety stop. Base of drawers open.
Dimensions: WxDxH = 140x567x340 mm

L Base of drawer in middle open.
Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

BRAVO 110 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SB110-32NBMV
For 6 x BRAVO 110 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x434x227 mm

BRAVO 110 Drawer Racks


Without Handle (Tray) Open Design¹



Example: Cat. No. SB110-53NTOV
For 15 x BRAVO 110 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x682x329 mm

BRAVO 130 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SB130-22NAOL
For 4 x BRAVO 130 Boxes.
Without safety stop. Base of drawers open.
Dimensions: WxDxH = 140x283x276 mm

L Base of drawer in middle open.
Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

BRAVO 130 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SB130-42NBMV
For 8 x BRAVO 130 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x567x283 mm

BRAVO 130 Drawer Racks

Without Handle (Tray) Open Design¹



Example: Cat. No. SB130-32NTOV
For 6 x BRAVO 130 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x416x276 mm



- All drawer racks can also be supplied with a closed design. (See illustrations on page 32+33). Please enquire if needed.
- All drawer racks can be produced with additional drawers for shallower boxes (e.g. BRAVO 50) to exploit the interior space of uprights to the maximum.

BRAVO Box long2 133 x 257 Cryo and Storage Boxes


Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 133 x 257 mm. Height = inner height.

Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.


Matching Racks

Please request!




			Heights						
Colour		1 or 2	25 mm	32 mm	50 mm*	75 mm	100 mm	110 mm	130 mm
	WH	1 standard 2 special	025	032	050	075	100	110	130
	BU	1 standard 2 special	025	032	050	075	100	110	130
	YE	1 standard 2 special	025	032	050	075	100	110	130
	GR	1 standard 2 special	025	032	050	075	100	110	130
	RE	1 standard 2 special	025	032	050	075	100	110	130

* Produced traditionally to the American design.
(Actual internal height: 45 mm)



How to determine
the item no.:



Example 50 mm height,
white, standard:

KB 2 0 50 WH 1

Grid Dividers for BRAVO Box long2


			Heights					
			15 mm	25 mm	30 mm	40 mm	65 mm	80 mm
Cells for vials			Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
16x32 up to ca. 7.2 mm Ø			RB20151632N	RB20251632N	RB20301632N	RB20401632N	RB20651632N	RB20801632N
15x30 up to ca. 7.7 mm Ø			RB20151530N	RB20251530N	RB20301530N	RB20401530N	RB20651530N	RB20801530N
14x28 up to ca. 8.3 mm Ø			RB20151428N	RB20251428N	RB20301428N	RB20401428N	RB20651428N	RB20801428N
13x26 up to ca. 9.0 mm Ø			RB20151326N	RB20251326N	RB20301326N	RB20401326N	RB20651326N	RB20801326N
12x24 up to ca. 9.8 mm Ø			RB20151224N	RB20251224N	RB20301224N	RB20401224N	RB20651224N	RB20801224N
11x22 up to ca. 10.8 mm Ø			RB20151122N	RB20251122N	RB20301122N	RB20401122N	RB20651122N	RB20801122N
10x20 up to ca. 12.0 mm Ø			RB20151020N	RB20251020N	RB20301020N	RB20401020N	RB20651020N	RB20801020N
9x18 up to ca. 13.6 mm Ø			RB20150918N	RB20250918N	RB20300918N	RB20400918N	RB20650918N	RB20800918N
8x16 up to ca. 15.0 mm Ø			RB20150816N	RB20250816N	RB20300816N	RB20400816N	RB20650816N	RB20800816N
7x14 up to ca. 17.3 mm Ø			RB20150714N	RB20250714N	RB20300714N	RB20400714N	RB20650714N	RB20800714N
6x12 up to ca. 20.3 mm Ø			RB20150612N	RB20250612N	RB20300612N	RB20400612N	RB20650612N	RB20800612N
5x10 up to ca. 24.5 mm Ø			RB20150510N	RB20250510N	RB20300510N	RB20400510N	RB20650510N	RB20800510N
4x8 up to ca. 30.8 mm Ø			RB20150408N	RB20250408N	RB20300408N	RB20400408N	RB20650408N	RB20800408N
3x6 up to ca. 41.2 mm Ø			RB20150306N	RB20250306N	RB20300306N	RB20400306N	RB20650306N	RB20800306N

When used with item number given in the overview (beginning with "RB2..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: BRAVO Box long2 of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 9x18 cells.

KB2050WH1 + 0300918N = Cat. No. KB2050WH10300918N



BRAVO Box short 133 x 92 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 133 x 92 mm. Height = inner height.

Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.


Matching Racks

Please request!

			Heights						
Colour	1 or 2		25 mm	32 mm	50 mm*	75 mm	100 mm	110 mm	130 mm
 WH	1 standard 2 special		025	032	050	075	100	110	130
 BU	1 standard 2 special		025	032	050	075	100	110	130
 YE	1 standard 2 special		025	032	050	075	100	110	130
 GR	1 standard 2 special		025	032	050	075	100	110	130
 RE	1 standard 2 special		025	032	050	075	100	110	130


* Produced traditionally to the American design.
(Actual internal height: 45 mm)

How to determine
the item no.:



Example 50 mm height,
white, standard:

KBK050WH1



Grid Dividers for BRAVO Box short

			Heights					
			15 mm	25 mm	30 mm	40 mm	65 mm	80 mm
Cells for vials			Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
16x12 up to ca. 7.2 mm Ø			RBK0151612N	RBK0251612N	RBK0301612N	RBK0401612N	RBK0651612N	RBK0801612N
15x11 up to ca. 7.7 mm Ø			RBK0151511N	RBK0251511N	RBK0301511N	RBK0401511N	RBK0651511N	RBK0801511N
14x10 up to ca. 8.3 mm Ø			RBK0151410N	RBK0251410N	RBK0301410N	RBK0401410N	RBK0651410N	RBK0801410N
13x9 up to ca. 9.0 mm Ø			RBK0151309N	RBK0251309N	RBK0301309N	RBK0401309N	RBK0651309N	RBK0801309N
12x8 up to ca. 9.8 mm Ø			RBK0151208N	RBK0251208N	RBK0301208N	RBK0401208N	RBK0651208N	RBK0801208N
11x7 up to ca. 10.8 mm Ø			RBK0151107N	RBK0251107N	RBK0301107N	RBK0401107N	RBK0651107N	RBK0801107N
10x6 up to ca. 12.0 mm Ø			RBK0151006N	RBK0251006N	RBK0301006N	RBK0401006N	RBK0651006N	RBK0801006N
9x6 up to ca. 13.6 mm Ø			RBK0150906N	RBK0250906N	RBK0300906N	RBK0400906N	RBK0650906N	RBK0800906N
8x5 up to ca. 15.0 mm Ø			RBK0150805N	RBK0250805N	RBK0300805N	RBK0400805N	RBK0650805N	RBK0800805N
7x4 up to ca. 17.3 mm Ø			RBK0150704N	RBK0250704N	RBK0300704N	RBK0400704N	RBK0650704N	RBK0800704N
6x4 up to ca. 20.3 mm Ø			RBK0150604N	RBK0250604N	RBK0300604N	RBK0400604N	RBK0650604N	RBK0800604N
5x3 up to ca. 24.5 mm Ø			RBK0150503N	RBK0250503N	RBK0300503N	RBK0400503N	RBK0650503N	RBK0800503N
4x2 up to ca. 30.8 mm Ø			RBK0150402N	RBK0250402N	RBK0300402N	RBK0400402N	RBK0650402N	RBK0800402N
3x2 up to ca. 41.2 mm Ø			RBK0150302N	RBK0250302N	RBK0300302N	RBK0400302N	RBK0650302N	RBK0800302N

When used with item number given in the overview (beginning with "RBK..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: BRAVO Box short of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 9x6 cells.

KBK050WH1 + 0300906N = Cat. No. KBK050WH10300906N




CellBox Cryo and Storage Boxes for 15 ml and 50 ml tubes

Useful for all tubes up to approx. 22 mm (15 ml) or 34 mm (50 ml) diameter and 125 mm maximum length. Standard box water-repellent fiberboard boxes suitable for ultra-low temperatures. Quality >special< with additional water-proof coating made of plastic. With and without inserted grid divider. Dimensions = External dimensions incl. lid.

		15 ml	
		CellBox Mini 5 x 5 Cells (122x122x128 mm)	CellBox Maxi 6 x 6 Cells (148x148x128 mm)
Colour	1 or 2	Cat. No.	Cat. No.
 WH	1 standard 2 special	KCMi1WH1-5x5 KCMi1WH2-5x5	KCMa1WH1-6x6 KCMa1WH2-6x6
 BU	1 standard 2 special	KCMi1BU1-5x5 KCMi1BU2-5x5	KCMa1BU1-6x6 KCMa1BU2-6x6
 YE	1 standard 2 special	KCMi1YE1-5x5 KCMi1YE2-5x5	KCMa1YE1-6x6 KCMa1YE2-6x6
 GR	1 standard 2 special	KCMi1GR1-5x5 KCMi1GR2-5x5	KCMa1GR1-6x6 KCMa1GR2-6x6
 RE	1 standard 2 special	KCMi1RE1-5x5 KCMi1RE2-5x5	KCMa1RE1-6x6 KCMa1RE2-6x6



		15 ml	
		CellBox Mini long 5 x 10 Cells (122x237x128 mm)	CellBox Maxi long 6 x 12 Cells (148x287x128 mm)
Colour	1 or 2	Cat. No.	Cat. No.
 WH	1 standard 2 special	KCMi2WH1-5x10 KCMi2WH2-5x10	KCMa2WH1-6x12 KCMa2WH2-6x12
 BU	1 standard 2 special	KCMi2BU1-5x10 KCMi2BU2-5x10	KCMa2BU1-6x12 KCMa2BU2-6x12
 YE	1 standard 2 special	KCMi2YE1-5x10 KCMi2YE2-5x10	KCMa2YE1-6x12 KCMa2YE2-6x12
 GR	1 standard 2 special	KCMi2GR1-5x10 KCMi2GR2-5x10	KCMa2GR1-6x12 KCMa2GR2-6x12
 RE	1 standard 2 special	KCMi2RE1-5x10 KCMi2RE2-5x10	KCMa2RE1-6x12 KCMa2RE2-6x12



		15 ml	
		CellBox Mini without Cells (122x122x128 mm)	CellBox Maxi without Cells (148x148x128 mm)
Colour	1 or 2	Cat. No.	Cat. No.
 WH	1 standard 2 special	KCMi1WH1 KCMi1WH2	KCMa1WH1 KCMa1WH2
 BU	1 standard 2 special	KCMi1BU1 KCMi1BU2	KCMa1BU1 KCMa1BU2
 YE	1 standard 2 special	KCMi1YE1 KCMi1YE2	KCMa1YE1 KCMa1YE2
 GR	1 standard 2 special	KCMi1GR1 KCMi1GR2	KCMa1GR1 KCMa1GR2
 RE	1 standard 2 special	KCMi1RE1 KCMi1RE2	KCMa1RE1 KCMa1RE2

		50 ml	
		CellBox Mini 3 x 3 Cells (122x122x128 mm)	CellBox Maxi 4 x 4 Cells (148x148x128 mm)
Colour	1 or 2	Cat. No.	Cat. No.
 WH	1 standard 2 special	KCMi1WH1-3x3 KCMi1WH2-3x3	KCMa1WH1-4x4 KCMa1WH2-4x4
 BU	1 standard 2 special	KCMi1BU1-3x3 KCMi1BU2-3x3	KCMa1BU1-4x4 KCMa1BU2-4x4
 YE	1 standard 2 special	KCMi1YE1-3x3 KCMi1YE2-3x3	KCMa1YE1-4x4 KCMa1YE2-4x4
 GR	1 standard 2 special	KCMi1GR1-3x3 KCMi1GR2-3x3	KCMa1GR1-4x4 KCMa1GR2-4x4
 RE	1 standard 2 special	KCMi1RE1-3x3 KCMi1RE2-3x3	KCMa1RE1-4x4 KCMa1RE2-4x4

		50 ml	
		CellBox Mini long 3 x 6 Cells (122x237x128 mm)	CellBox Maxi long 4 x 8 Cells (148x287x128 mm)
Colour	1 or 2	Cat. No.	Cat. No.
 WH	1 standard 2 special	KCMi2WH1-3x6 KCMi2WH2-3x6	KCMa2WH1-4x8 KCMa2WH2-4x8
 BU	1 standard 2 special	KCMi2BU1-3x6 KCMi2BU2-3x6	KCMa2BU1-4x8 KCMa2BU2-4x8
 YE	1 standard 2 special	KCMi2YE1-3x6 KCMi2YE2-3x6	KCMa2YE1-4x8 KCMa2YE2-4x8
 GR	1 standard 2 special	KCMi2GR1-3x6 KCMi2GR2-3x6	KCMa2GR1-4x8 KCMa2GR2-4x8
 RE	1 standard 2 special	KCMi2RE1-3x6 KCMi2RE2-3x6	KCMa2RE1-4x8 KCMa2RE2-4x8

		50 ml	
		CellBox Mini long without Cells (122x237x128 mm)	CellBox Maxi long without Cells (148x287x128 mm)
Colour	1 or 2	Cat. No.	Cat. No.
 WH	1 standard 2 special	KCMi2WH1 KCMi2WH2	KCMa2WH1 KCMa2WH2
 BU	1 standard 2 special	KCMi2BU1 KCMi2BU2	KCMa2BU1 KCMa2BU2
 YE	1 standard 2 special	KCMi2YE1 KCMi2YE2	KCMa2YE1 KCMa2YE2
 GR	1 standard 2 special	KCMi2GR1 KCMi2GR2	KCMa2GR1 KCMa2GR2
 RE	1 standard 2 special	KCMi2RE1 KCMi2RE2	KCMa2RE1 KCMa2RE2

Grid Dividers for CellBox (only for empty boxes or repeat orders)

Cells	Cat. No.	Cells	Cat. No.	Cells	Cat. No.	Cells	Cat. No.
Mini 15 ml 5x5	RC0900505N	Mini 50 ml 3x3	RC0900303N	Maxi 15 ml 6x6	RC0900606N	Maxi 50 ml 4x4	RC0900404N
Mini long 15 ml 5x10	RC0900510N	Mini long 50 ml 3x6	RC0900306N	Maxi long 15 ml 6x12	RC0900612N	Maxi long 50 ml 4x8	RC0900408N
Mini long 10 ml 8x12	RC0900812N						



CellBox Mini Vertical Racks

Useful for CellBox Mini and all boxes up to external dimensions 122x122x128 mm (WxDxH). *

Cat. No.	Bays	WxDxH (mm)
TCMi1-02AN	2	127x127x266
TCMi1-03AN	3	127x127x396
TCMi1-04AN	4	127x127x526
TCMi1-05AN	5	127x127x656
TCMi1-06AN	6	127x127x786

CellBox Mini long Vertical Racks

Useful for CellBox Mini long and all boxes up to external dimensions 122x237x128 mm (WxDxH). *

Cat. No.	Bays	WxDxH (mm)
TCMi2-02AN	2	127x245x266
TCMi2-03AN	3	127x245x396
TCMi2-04AN	4	127x245x526
TCMi2-05AN	5	127x245x656
TCMi2-06AN	6	127x245x786



Vertical racks with folding handle positioned crosswise as standard.

Alternatives: Folding handle positioned lengthwise or hand grip. Please enquire if needed.



CellBox Mini Horizontal Racks

Useful for CellBox Mini and all boxes up to external dimensions 122x122x128 mm (WxDxH). *

Cat. No.	Bays Depth x Height	WxDxH (mm)
SCMi1-22AN	2 x 2	128x258x264
SCMi1-23AN	2 x 3	128x258x394
SCMi1-32AN	3 x 2	128x384x264
SCMi1-33AN	3 x 3	128x384x394
SCMi1-42AN	4 x 2	128x509x264
SCMi1-43AN	4 x 3	128x509x394
SCMi1-52AN	5 x 2	128x635x264
SCMi1-53AN	5 x 3	128x635x394

CellBox Mini long Horizontal Racks

Useful for CellBox Mini long and all boxes up to external dimensions 122x237x128 mm (WxDxH). *

Cat. No.	Bays Depth x Height	WxDxH (mm)
SCMi2-22AN	2 x 2	128x494x264
SCMi2-23AN	2 x 3	128x494x394
SCMi2-32AN	3 x 2	128x738x264
SCMi2-33AN	3 x 3	128x738x394



Horizontal racks with folding handle positioned vertically as standard.

Alternatives: Folding handle positioned horizontally or hand grip. Please enquire if needed.



* Dimensions = External dimensions incl. lid.

CellBox Maxi Vertical Racks

Useful for CellBox Maxi and all boxes up to external dimensions 148x148x128 mm (WxDxH). *

Cat. No.	Bays	WxDxH (mm)
TCMa1-02AN	2	153x153x266
TCMa1-03AN	3	153x153x396
TCMa1-04AN	4	153x153x526
TCMa1-05AN	5	153x153x656
TCMa1-06AN	6	153x153x786

CellBox Maxi long Vertical Racks

Useful for CellBox Maxi long and all boxes up to external dimensions 148x287x128 mm (WxDxH). *

Cat. No.	Bays	WxDxH (mm)
TCMa2-02AN	2	153x296x266
TCMa2-03AN	3	153x296x396
TCMa2-04AN	4	153x296x526
TCMa2-05AN	5	153x296x656
TCMa2-06AN	6	153x296x786

CellBox Maxi Horizontal Racks

Useful for CellBox Maxi and all boxes up to external dimensions 148x148x128 mm (WxDxH). *


Cat. No.	Bays Depth x Height	WxDxH (mm)
SCMa1-22AN	2 x 2	154x310x264
SCMa1-23AN	2 x 3	154x310x394
SCMa1-32AN	3 x 2	154x462x264
SCMa1-33AN	3 x 3	154x462x394
SCMa1-42AN	4 x 2	154x613x264
SCMa1-43AN	4 x 3	154x613x394
SCMa1-52AN	5 x 2	154x765x264
SCMa1-53AN	5 x 3	154x765x394

CellBox Maxi long Horizontal Racks

Useful for CellBox Maxi long and all boxes up to external dimensions 148x287x128 mm (WxDxH). *

Cat. No.	Bays Depth x Height	WxDxH (mm)
SCMa2-22AN	2 x 2	154x596x264
SCMa2-23AN	2 x 3	154x596x394

Drawer Racks for CellBox on request.



DELTA Box 133 x 133 Cryo and Storage Boxes


Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With extended lid. External dimensions 133 x 133 mm. Height = inner height.

Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.


Colour			1 or 2	25 mm	32 mm	50 mm	75 mm	100 mm	110 mm	130 mm
				Heights						
white	WH	1	standard	see BRAVO Box	see BRAVO Box	050	075	100	110	130
		2	special							
blue	BU	1	standard			050	075	100	110	130
		2	special							
yellow	YE	1	standard			050	075	100	110	130
		2	special							
green	GR	1	standard			050	075	100	110	130
		2	special							
red	RE	1	standard			050	075	100	110	130
		2	special							

How to determine the item no.:



Example 50 mm height, white, standard:

KD1050WH1



Grid Dividers for DELTA Box

			15 mm	25 mm	30 mm	40 mm	65 mm	80 mm
			Heights					
Cells for vials			Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
16x16	up to	ca. 7.2 mm Ø	RB10151616N	RB10251616N	RB10301616N	RB10401616N	RB10651616N	RB10801616N
15x15	up to	ca. 7.7 mm Ø	RB10151515N	RB10251515N	RB10301515N	RB10401515N	RB10651515N	RB10801515N
14x14	up to	ca. 8.3 mm Ø	RB10151414N	RB10251414N	RB10301414N	RB10401414N	RB10651414N	RB10801414N
13x13	up to	ca. 9.0 mm Ø	RB10151313N	RB10251313N	RB10301313N	RB10401313N	RB10651313N	RB10801313N
12x12	up to	ca. 9.8 mm Ø	RB10151212N	RB10251212N	RB10301212N	RB10401212N	RB10651212N	RB10801212N
11x11	up to	ca. 10.8 mm Ø	RB10151111N	RB10251111N	RB10301111N	RB10401111N	RB10651111N	RB10801111N
10x10	up to	ca. 12.0 mm Ø	RB10151010N	RB10251010N	RB10301010N	RB10401010N	RB10651010N	RB10801010N
9x9	up to	ca. 13.6 mm Ø	RB10150909N	RB10250909N	RB10300909N	RB10400909N	RB10650909N	RB10800909N
8x8	up to	ca. 15.0 mm Ø	RB10150808N	RB10250808N	RB10300808N	RB10400808N	RB10650808N	RB10800808N
7x7	up to	ca. 17.3 mm Ø	RB10150707N	RB10250707N	RB10300707N	RB10400707N	RB10650707N	RB10800707N
6x6	up to	ca. 20.3 mm Ø	RB10150606N	RB10250606N	RB10300606N	RB10400606N	RB10650606N	RB10800606N
5x5	up to	ca. 24.5 mm Ø	RB10150505N	RB10250505N	RB10300505N	RB10400505N	RB10650505N	RB10800505N
4x4	up to	ca. 30.8 mm Ø	RB10150404N	RB10250404N	RB10300404N	RB10400404N	RB10650404N	RB10800404N
3x3	up to	ca. 41.2 mm Ø	RB10150303N	RB10250303N	RB10300303N	RB10400303N	RB10650303N	RB10800303N

When used with item number given in the overview (beginning with "RB1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: DELTA Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 9x9 cells.

KD1050WH1 + 0300909N = Cat. No. KD1050WH10300909N



MIKE Box 142 x 142 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 142 x 142 mm. Height = inner height.

Other heights can be produced on request.


Grid dividers can be ordered separately or can be supplied fitted.

Matching Racks

Please request!

			Heights						
Colour	1 or 2		25 mm	32 mm	50 mm	75 mm	100 mm	110 mm	130 mm
 WH	1 standard 2 special		025	032	050	075	100	110	130
 BU	1 standard 2 special		025	032	050	075	100	110	130
 YE	1 standard 2 special		025	032	050	075	100	110	130
 GR	1 standard 2 special		025	032	050	075	100	110	130
 RE	1 standard 2 special		025	032	050	075	100	110	130

How to determine the item no.:



Example 50 mm height, white, standard:

KM1050WH1



Grid Dividers for MIKE Box

			Heights					
			15 mm	25 mm	30 mm	40 mm	65 mm	80 mm
Cells for vials			Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
16x16 up to ca. 7.7 mm Ø			RM10151616N	RM10251616N	RM10301616N	RM10401616N	RM10651616N	RM10801616N
15x15 up to ca. 8.3 mm Ø			RM10151515N	RM10251515N	RM10301515N	RM10401515N	RM10651515N	RM10801515N
14x14 up to ca. 9.0 mm Ø			RM10151414N	RM10251414N	RM10301414N	RM10401414N	RM10651414N	RM10801414N
13x13 up to ca. 9.7 mm Ø			RM10151313N	RM10251313N	RM10301313N	RM10401313N	RM10651313N	RM10801313N
12x12 up to ca. 10.6 mm Ø			RM10151212N	RM10251212N	RM10301212N	RM10401212N	RM10651212N	RM10801212N
11x11 up to ca. 11.6 mm Ø			RM10151111N	RM10251111N	RM10301111N	RM10401111N	RM10651111N	RM10801111N
10x10 up to ca. 12.8 mm Ø			RM10151010N	RM10251010N	RM10301010N	RM10401010N	RM10651010N	RM10801010N
9x9 up to ca. 14.3 mm Ø			RM10150909N	RM10250909N	RM10300909N	RM10400909N	RM10650909N	RM10800909N
8x8 up to ca. 16.2 mm Ø			RM10150808N	RM10250808N	RM10300808N	RM10400808N	RM10650808N	RM10800808N
7x7 up to ca. 18.6 mm Ø			RM10150707N	RM10250707N	RM10300707N	RM10400707N	RM10650707N	RM10800707N
6x6 up to ca. 21.8 mm Ø			RM10150606N	RM10250606N	RM10300606N	RM10400606N	RM10650606N	RM10800606N
5x5 up to ca. 26.3 mm Ø			RM10150505N	RM10250505N	RM10300505N	RM10400505N	RM10650505N	RM10800505N
4x4 up to ca. 33.0 mm Ø			RM10150404N	RM10250404N	RM10300404N	RM10400404N	RM10650404N	RM10800404N
3x3 up to ca. 44.2 mm Ø			RM10150303N	RM10250303N	RM10300303N	RM10400303N	RM10650303N	RM10800303N

When used with item number given in the overview (beginning with "RM1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: MIKE Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 10x10 cells.

KM1050WH1 + 0301010N = Cat. No. KM1050WH10301010N



NANU Box 76 x 76 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With extended lid. External dimensions 76 x 76 mm. Height = inner height.


Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.


Matching Racks

See below

Colour			1 or 2	25 mm	32 mm	50 mm	Heights 75 mm	100 mm	110 mm	130 mm
	WH	1	standard	025	032	050	075	100	110	130
		2	special							
	BU	1	standard	025	032	050	075	100	110	130
		2	special							
	YE	1	standard	025	032	050	075	100	110	130
		2	special							
	GR	1	standard	025	032	050	075	100	110	130
		2	special							
	RE	1	standard	025	032	050	075	100	110	130
		2	special							



How to determine
the item no.:



Example 50 mm height,
white, standard:

KN 1 050 WH 1

Grid Dividers for NANU Box

			15 mm	25 mm	30 mm	40 mm	65 mm	80 mm
Cells for vials			Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.	Cat. No.
7x7	up to	ca. 8.0 mm Ø	RN10150707N	RN10250707N	RN10300707N	RN10400707N	RN10650707N	RN10800707N
6x6	up to	ca. 10.0 mm Ø	RN10150606N	RN10250606N	RN10300606N	RN10400606N	RN10650606N	RN10800606N
5x5	up to	ca. 13.5 mm Ø	RN10150505N	RN10250505N	RN10300505N	RN10400505N	RN10650505N	RN10800505N
4x4	up to	ca. 15.0 mm Ø	RN10150404N	RN10250404N	RN10300404N	RN10400404N	RN10650404N	RN10800404N
3x3	up to	ca. 21.0 mm Ø	RN10150303N	RN10250303N	RN10300303N	RN10400303N	RN10650303N	RN10800303N

When used with item number given in the overview (beginning with "RN1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: NANU Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 5x5 cells.

KN1050WH1 + 0300505N = Cat. No. KN1050WH10300505N

NANU 50 Vertical Racks • 82x82 mm (WxD)

NANU 50 Horizontal Racks • Width 83 mm

>> Useful for NANU 50 and all boxes up to external dimensions 76x76x53 mm (WxDxH). <<

Cat. No.	Bays	Height (mm)	Cat. No.	Bays	Height (mm)
TN50-04AN	4	226	TN50-10AN	10	556
TN50-05AN	5	281	TN50-11AN	11	611
TN50-06AN	6	336	TN50-12AN	12	666
TN50-07AN	7	391	TN50-13AN	13	721
TN50-08AN	8	446	TN50-14AN	14	776
TN50-09AN	9	501			



All vertical racks can be supplied with a safety rod if required.
This increases the depth by 5 mm to 87 mm.

Cat. No.	Bays D x H	Dimensions DxH (mm)	Cat. No.	Bays D x H	Dimensions DxH (mm)
SN50-44AN	4 x 4	329x224	SN50-65AN	6 x 5	490x279
SN50-45AN	4 x 5	329x279	SN50-66AN	6 x 6	490x334
SN50-46AN	4 x 6	329x334	SN50-74AN	7 x 4	571x224
SN50-54AN	5 x 4	410x224	SN50-75AN	7 x 5	571x279
SN50-55AN	5 x 5	410x279	SN50-76AN	7 x 6	571x334
SN50-56AN	5 x 6	410x334	SN50-85AN	8 x 5	651x279
SN50-64AN	6 x 4	490x224	SN50-86AN	8 x 6	651x334

Additional vertical and horizontal racks as well as drawer racks for NANU box on request.

Mikrotiter Vertical Racks


For MT-Plates 86 x 128 mm (WxD). Stainless steel.

For Microtiter Plates or DeepWells of all familiar brands.




Up to 18 mm plate height.

Single



Double




Example:
Cat. No. TM18-22AN
Vertical rack with
22 bays
(Individual bays)



Example:
Cat. No. TM18-19x2AN
Vertical rack with
19 bays (2 plates next to
each other in each bay)


Up to 39 mm plate height.

Single



Cat. No. TM39- . . AN


Double



Cat. No. TM39- AN


Up to 53 mm plate height.

Single



Cat. No. TM53- . . AN

Double




Cat. No. TM53- AN

Cat. No. TM53- . . AN

Cat. No. TM53- AN

Personal Rack



If the examples shown
here do not meet your
storage needs,
please let us know
your requirements:

- Maximum height of the rack?
- Height of the plate stack or required internal height of the bays?
- How many racks are needed?


You will receive an individual quotation from us.



All vertical racks can be supplied with a safety rod if required.
Depth plus 5 mm.
Please enquire if needed.


Up to 45 mm plate height.

Single



Cat. No. TM45- . . AN


Double



Cat. No. TM45- AN


Up to 58 mm plate height.

Single



Cat. No. TM58- . . AN

Double




Cat. No. TM58- AN

For Microtiter Plates or DeepWells of all familiar brands.



Mikrotiter Horizontal Racks with one intermediate shelf

For MT plates/DeepWells of different heights or large stacks of plates.



Folding handle positioned vertically as standard.

Alternatives:
Folding handle positioned horizontally or hand grip.

Please enquire if needed.

Standard Heights


Cat. No.	Bays	Inner Height	Dimensions
Depth x Height	Bay (mm)	DxH (mm)	
SMX169-3-02AN	3x2	82	276x169
SMX169-4-02AN	4x2	82	365x169
SMX169-5-02AN	5x2	82	455x169
SMX169-6-02AN	6x2	82	544x169
SMX169-7-02AN	7x2	82	634x169
SMX169-8-02AN	8x2	82	723x169
SMX224-3-02AN	3x2	109	276x224
SMX224-4-02AN	4x2	109	365x224
SMX224-5-02AN	5x2	109	455x224
SMX224-6-02AN	6x2	109	544x224
SMX224-7-02AN	7x2	109	634x224
SMX224-8-02AN	8x2	109	723x224

Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
Depth x Height	Bay (mm)	DxH (mm)	
SMX250-3-02AN	3x2	122	276x250
SMX250-4-02AN	4x2	122	365x250
SMX250-5-02AN	5x2	122	455x250
SMX250-6-02AN	6x2	122	544x250
SMX250-7-02AN	7x2	122	634x250
SMX250-8-02AN	8x2	122	723x250
SMX279-3-02AN	3x2	137	276x279
SMX279-4-02AN	4x2	137	365x279
SMX279-5-02AN	5x2	137	455x279
SMX279-6-02AN	6x2	137	544x279
SMX279-7-02AN	7x2	137	634x279
SMX279-8-02AN	8x2	137	723x279


Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
Depth x Height	Bay (mm)	DxH (mm)	
SMX334-3-02AN	3x2	164	276x334
SMX334-4-02AN	4x2	164	365x334
SMX334-5-02AN	5x2	164	455x334
SMX334-6-02AN	6x2	164	544x334
SMX334-7-02AN	7x2	164	634x334
SMX334-8-02AN	8x2	164	723x334
SMX389-3-02AN	3x2	192	276x389
SMX389-4-02AN	4x2	192	365x389
SMX389-5-02AN	5x2	192	455x389
SMX389-6-02AN	6x2	192	544x389
SMX389-7-02AN	7x2	192	634x389
SMX389-8-02AN	8x2	192	723x389



Mikrotiter Horizontal Racks with two intermediate shelves

For MT plates/DeepWells of different heights or large stacks of plates.



All horizontal racks can be supplied with safety rods if required.

Add 5 mm to height.

Please enquire if needed.

Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
Depth x Height	Bay (mm)	DxH (mm)	
SMX169-3-03AN	3x3	54	276x169
SMX169-4-03AN	4x3	54	365x169
SMX169-5-03AN	5x3	54	455x169
SMX169-6-03AN	6x3	54	544x169
SMX169-7-03AN	7x3	54	634x169
SMX169-8-03AN	8x3	54	723x169
SMX224-3-03AN	3x3	73	276x224
SMX224-4-03AN	4x3	73	365x224
SMX224-5-03AN	5x3	73	455x224
SMX224-6-03AN	6x3	73	544x224
SMX224-7-03AN	7x3	73	634x224
SMX224-8-03AN	8x3	73	723x224

Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
Depth x Height	Bay (mm)	DxH (mm)	
SMX250-3-03AN	3x3	81	276x250
SMX250-4-03AN	4x3	81	365x250
SMX250-5-03AN	5x3	81	455x250
SMX250-6-03AN	6x3	81	544x250
SMX250-7-03AN	7x3	81	634x250
SMX250-8-03AN	8x3	81	723x250
SMX279-3-03AN	3x3	91	276x279
SMX279-4-03AN	4x3	91	365x279
SMX279-5-03AN	5x3	91	455x279
SMX279-6-03AN	6x3	91	544x279
SMX279-7-03AN	7x3	91	634x279
SMX279-8-03AN	8x3	91	723x279


Standard Heights

Cat. No.	Bays	Inner Height	Dimensions
Depth x Height	Bay (mm)	DxH (mm)	
SMX334-3-03AN	3x3	109	276x334
SMX334-4-03AN	4x3	109	365x334
SMX334-5-03AN	5x3	109	455x334
SMX334-6-03AN	6x3	109	544x334
SMX334-7-03AN	7x3	109	634x334
SMX334-8-03AN	8x3	109	723x334
SMX389-3-03AN	3x3	128	276x389
SMX389-4-03AN	4x3	128	365x389
SMX389-5-03AN	5x3	128	455x389
SMX389-6-03AN	6x3	128	544x389
SMX389-7-03AN	7x3	128	634x389
SMX389-8-03AN	8x3	128	723x389

Mikrotiter Horizontal Racks

For MT-Plates 86 x 128 mm (WxD). Stainless steel.

For Microtiter Plates or DeepWells of all familiar brands.



Up to 39 mm height.

Cat. No. SM39-6x6 AN

334	8x08	7x08	6x08	5x08	4x08	3x08
291	8x07	7x07	6x07	5x07	4x07	3x07
250	8x06	7x06	6x06	5x06	4x06	3x06
209	8x05	7x05	6x05	5x05	4x05	3x05
169	8x04	7x04	6x04	5x04	4x04	3x04

723 634 544 455 365 276 139



Example: Cat. No. **SM39 - 6 x 0 6 AN**

Example: Cat. No. **SM33 - 3x33AN**
MT horizontal rack with 36 bays. For plates or stacks of up to 39 mm in height.
Dimensions: WxDxH = 139x544x250 mm

Up to 45 mm height.

Cat. No. SM45- x  AN

334	8x07	7x07	6x07	5x07	4x07	3x07
286	8x06	7x06	6x06	5x06	4x06	3x06
239	8x05	7x05	6x05	5x05	4x05	3x05
192	8x04	7x04	6x04	5x04	4x04	3x04

723 634 544 455 365 276 139

Up to 53 mm height.

Cat. No. SM53-AN

Height (mm)	8x06	7x06	6x06	5x06	4x06	3x06
334	8x06	7x06	6x06	5x06	4x06	3x06
279	8x05	7x05	6x05	5x05	4x05	3x05
224	8x04	7x04	6x04	5x04	4x04	3x04
169	8x03	7x03	6x03	5x03	4x03	3x03

Dimensions (mm): 723, 634, 544, 455, 365, 276, 139

Up to 58 mm height.

Cat. No. **SM58-** x **AN**

364	8x06	7x06	6x06	5x06	4x06	3x06
304	8x05	7x05	6x05	5x05	4x05	3x05
244	8x04	7x04	6x04	5x04	4x04	3x04
	723	634	544	455	365	276

1390



Folding handle positioned vertically as standard. **Alternatives:** Folding handle positioned horizontally or hand grip. Please enquire if needed.



All horizontal racks can be supplied with safety rods if required.
Add 5 mm to height.
Please enquire if needed.


Mikrotiter Drawer Racks

For MT-Plates 86 x 128 mm (WxD). Stainless steel.

**For Microtiter Plates or DeepWells of all familiar brands.
Base of drawers closed. Draw front closed to maximum.**

Up to 45 mm plate height. With grip rail.

Open Design




M = With Stop **O** = Without Stop



Example: Cat. No. **SM45 - 5 x 0 5 PCOV**
MT drawer rack for 25 microtiter plates without safety stop.
Dimensions: WxDxH = 138x446x249 mm

Closed Design




M = With Stop **O** = Without Stop



Example: Cat. No. **SM45 - 4 x 0 6 HCM V**
MT drawer rack for 24 microtiter plates with safety stop.
Dimensions: WxDxH = 138x390x320 mm


Up to 53 mm plate height. With grip rail.

Open Design



M = With Stop **O** = Without Stop


Closed Design



M = With Stop **O** = Without Stop


Up to 58 mm plate height. With grip rail.

Open Design





M = With Stop O = Without Stop

Closed Design



M = With Stop O = Without Stop



Nothing is impossible!

Well-organized freezers and refrigerators are an essential prerequisite for systematic, efficient and ultimately successful work – not only in scientific laboratories.


We provide the means to fulfil this prerequisite and devote ourselves to the task from a sense of conviction and with passion.


In fact, it is a passion for practical details and our liking for intelligent solutions.

This catalogue contains many suggestions for the interior fittings of all brands of refrigerator and freezer cabinets as well as freezer chests.

If you are looking for a solution to a problem or need assistance in setting up your refrigeration equipment, please contact us.

We will always be able to offer you a solution and will be glad to help you.






National Lab GmbH
Refrigeration Technology
 Grambeker Weg 157
 D-23879 Moelln
 Tel: +49 45 42/84 91-70 · Fax: -71
 E-mail: LabSales@NationalLab.com
 www.Kryobox.de

Inventory Systems

for Freezers and Refrigerators

2011

Supplement



Online Information
www.Kryobox.de

EPi® Box made of Polypropylene

130 x 130 mm



EPi® 37 / 10x10 • height 37 mm fix.

Inner height 35 mm. 10x10 holes for 100 vials 0.5 ml Safe-Lock® up to 8.5 mm diameter. Without ID code.

Cat. No.	Colour
KEA37-C0100NA	natural
KEA37-C0100NY	neon-yellow
KEA37-C0100NB	neon-blue
KEA37-C0100NG	neon-green
KEA37-C0100NP	neon-red/pink
KEA37-C0100NO	neon-orange
KEA37-C0100WH	white

Cat. No.	Colour
KEA37-C0100YE	yellow
KEA37-C0100BU	blue
KEA37-C0100GR	green
KEA37-C0100RE	red
KEA37-C0100VT	violet
KEA37-C0100BL	black
KEA37-C0100BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 37 / 10x10 • height 37 mm fix.

Inner height 35 mm. 10x10 holes for 100 vials 0.5 ml Safe-Lock® up to 8.5 mm diameter. With alphanumeric ID code.

Cat. No.	Colour
KEB37-C0100NA	natural
KEB37-C0100NY	neon-yellow
KEB37-C0100NB	neon-blue
KEB37-C0100NG	neon-green
KEB37-C0100NP	neon-red/pink
KEB37-C0100NO	neon-orange
KEB37-C0100WH	white

Cat. No.	Colour
KEB37-C0100YE	yellow
KEB37-C0100BU	blue
KEB37-C0100GR	green
KEB37-C0100RE	red
KEB37-C0100VT	violet
KEB37-C0100BL	black
KEB37-C0100BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples




EPi® 45 / 8x8 • height 45-53 mm variable.

8x8 holes for 64 reaction vessels 1.5 ml and 2.0 ml Safe-Lock® up to 11 mm diameter. Without ID code.

Cat. No.	Colour
KEA45-V064NA	natural
KEA45-V064NY	neon-yellow
KEA45-V064NB	neon-blue
KEA45-V064NG	neon-green
KEA45-V064NP	neon-red/pink
KEA45-V064NO	neon-orange
KEA45-V064WH	white

Cat. No.	Colour
KEA45-V064YE	yellow
KEA45-V064BU	blue
KEA45-V064GR	green
KEA45-V064RE	red
KEA45-V064VT	violet
KEA45-V064BL	black
KEA45-V064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples



EPi® 50 / 8x8 • height 52 mm fix.

Inner height 50 mm. 8x8 holes for 64 reaction vessels 1.5 ml and 2.0 ml Safe-Lock® up to 11 mm diameter. Without ID code.

Cat. No.	Colour
KEA50-C064NA	natural
KEA50-C064NY	neon-yellow
KEA50-C064NB	neon-blue
KEA50-C064NG	neon-green
KEA50-C064NP	neon-red/pink
KEA50-C064NO	neon-orange
KEA50-C064WH	white

Cat. No.	Colour
KEA50-C064YE	yellow
KEA50-C064BU	blue
KEA50-C064GR	green
KEA50-C064RE	red
KEA50-C064VT	violet
KEA50-C064BL	black
KEA50-C064BB	black/black*

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 34

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 34

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

Matching Racks

Vertical Racks
page 29

Horizontal Racks
page 30-31

Drawer Racks
page 32-33

VICTOR Box + STAR Box

made of Polypropylene, with hinged lid



VICTOR Box 57 / 10x10 • height 57 mm fix.

With hinged lid. Inner height 50 mm. 10x10 divider holds 100 vials up to 12,5 mm diameter. With alphanumeric ID code on the lid.

Dimensions WxDxH = 141x151x57 mm

Cat. No.	Colour
KVH57-S100NA	natural
KVH57-S100NY	neon-yellow in preparation
KVH57-S100NB	neon-blue in preparation
KVH57-S100NG	neon-green in preparation
KVH57-S100NP	neon-red/pink in preparation
KVH57-S100NO	neon-orange in preparation
KVH57-S100BB	black/black* in preparation

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
Please request!

Horizontal Racks
Please request!

Drawer Racks
Please request!



VICTOR Box 59 / 5x10 • height 59 mm fix.

With hinged lid. Inner height 50 mm. 5x10 divider holds 50 vials up to 12,0 mm diameter. With alphanumeric ID code on 2 sides.

Dimensions WxDxH = 141x85x59 mm

Cat. No.	Colour
KVA59-S50NA	natural
KVA59-S50NY	neon-yellow in preparation
KVA59-S50NB	neon-blue in preparation
KVA59-S50NG	neon-green in preparation
KVA59-S50NP	neon-red/pink in preparation
KVA59-S50NO	neon-orange in preparation
KVA59-S50BB	black/black* in preparation

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
Please request!

Horizontal Racks
Please request!

Drawer Racks
Please request!



STAR Box 57 / 9x9 • height 57 mm fix.

With hinged lid. Inner height 50 mm. 9x9 divider holds 81 vials up to 12,5 mm diameter. With alphanumeric ID code on the lid.

Dimensions WxDxH = 128x138x57 mm

Cat. No.	Colour
KSH57-S81NA	natural
KSH57-S81NY	neon-yellow in preparation
KSH57-S81NB	neon-blue in preparation
KSH57-S81NG	neon-green in preparation
KSH57-S81NP	neon-red/pink in preparation
KSH57-S81NO	neon-orange in preparation
KSH57-S81BB	black/black* in preparation

*black/black: Lower half and lid black for light-sensitive samples

Matching Racks

Vertical Racks
Please request!

Horizontal Racks
Please request!

Drawer Racks
Please request!



Straight-Side Wide-Mouth Jars

Large wide-necked opening for easy filling. Made of transparent polycarbonate with white or blue PP lid. Window-clear. High impact resistance. Non-toxic. Stackable. Popular screw-top jar also for freezing down to -90°C. Fully unscrew the lid prior to autoclaving.

Cat. No.	Content	Pack	Ø Lid	Height
LPB21160125	125 ml	4 pcs.	70 mm	74 mm
LPB21160250	250 ml	4 pcs.	70 mm	119 mm
LPB21160500	500 ml	4 pcs.	120 mm	86 mm
LPB21161000	1000 ml	4 pcs.	120 mm	142 mm

Isotherm Quick Freeze + Isotherm Mini Cooler



Isotherm Quick Freeze -70°C

For 12 tubes from 0.5 to 2.0 ml. With alphanumeric ID code on the transparent lid and the lower section.

For rapid and repeated refrigeration in the molecular biology laboratory. Can be used repeatedly for up to 45 minutes. Made of resistant polycarbonate filled with non-toxic insulating fluid. No direct contact between the tube and the insulating fluid. With non-slip feet.

Volume	Time taken to attain			
	-20°C (isopropanol)		0°C (water)	
1.5 ml*	with isopropanol	2.0 min.	with water	3.2 min.
0.5 ml*	with isopropanol/adaptor	2.8 min.	with water	4.0 min.

* Time taken with micro centrifuge tubes (ca.)

Cat. No.	Colour	Cells	Vials (ml)	WxDxH (mm)
LCCQ51140012	white	4x3	0.5 to 2.0	151x108x125



Isotherm Mini Cooler 0.5 ml / 2.0 ml

For tubes from 0.5 to 2.0 ml. With alphanumeric ID code on the lid and the lower section.

Temperature-conditioning of reagents (enzymes, DNA, RNA, cell suspensions, etc.) on the laboratory bench. Made of resistant polycarbonate filled with non-toxic insulating fluid. No direct contact between the tube and the insulating fluid. With non-slip feet.

Cat. No.	Temp.	Colour	Cells	Lid	Storage time	WxDxH (mm)
LCCM51160012	0°C	red	4x3	without gel	<1°C/3.5 hrs.	151x108x125
LCCM51150012	-20°C	yellow	4x3	without gel	<-15°C/1 hrs.	151x108x125
LCCM51160032	0°C	red	8x4	with gel	<1°C/5 hrs.	243x157x100**
LCCM51160032C	0°C	red	8x4	without gel	<1°C/3.5 hrs.	243x157x100**
LCCM51150020	-20°C	yellow	5x4	without gel	<-15°C/1 hrs.	220x163x113**
LCCM51150032	-20°C	yellow	8x4	with gel	<-15°C/2 hrs.	243x157x100**
LCCM51150032C	-20°C	yellow	8x4	without gel	<-15°C/1 hrs.	243x157x100**

** Height without fastening clip. With fastening clip in place, the height increases by approx. 50 mm.



Isotherm Mini Cooler 15 ml + Vacutainer¹⁾

For 12 tubes, e. g. 15 ml vials alternatively Vacutainer¹⁾. With alphanumeric ID code on the transparent lid and the lower section.

Made of resistant polycarbonate filled with non-toxic insulating fluid. No direct contact between the tube and the insulating fluid. With non-slip feet.

Cat. No.	Temp.	Colour	Cells	Lid	Storage time	WxDxH (mm)
LCCM51161600	0°C	red	4x3	without gel	<1°C/5 hrs.	197x140x190 for vials Ø16–17 mm (e. g. 15 ml centrifuge tubes)
LCCM51161300	0°C	red	4x3	without gel	<1°C/5 hrs.	197x140x190 for vials Ø12–13 mm (e. g. Vacutainer ¹⁾)

¹⁾ Vacutainer is a registered trademark of Becton Dickinson.

Isotherm Mini Cooler + Cryo-Freezer



Isotherm Mini Cooler -20°C PCR

Specially developed to protect enzymes and other biological specimens while working at the laboratory bench. Removable aluminium block with 96 holes in universal 96-well standard format. Mini Cooler made of resistant polycarbonate filled with non-toxic insulating fluid. With non-slip feet.

With alphanumeric ID code on the aluminium block.

0.2-ml PCR tubes individually or in strips as well as 96-well plates are maintained at between -15°C and -20°C for up to 2 hours.

Cat. No.	Temp.	Cells	Lid	Storage time	WxDxH (mm)
LCCM96PCR	-20°C	12x8	without gel	<-20°C/2 hrs.	220x163x150**

** Height without fastening clip. With fastening clip in place, the height increases by approx. 50 mm.



Mr. Frosty polycarbonate cryo-freezer

Lid and tube holders made of HDPE with foam liner. Enables critical, repeatable cooling rate of approx. 1°C/minute for successful cryo-preservation and optimum cell yield.

Easy to use with 100% isopropanol and a freezing device.

For up to 18 cryo-tubes from 1.0 to 2.0 ml. No direct contact between the tubes and the alcohol. Instructions for use printed on device. Stackable.

Cat. No.	Vials	Dimensions (mm)
LCC5100	18	H x Ø = 86 x 117, Lid-Ø 120

Straight-Side Wide-Mouth Jars for specimen storage



Straight-Side Wide-Mouth Jars

Large wide-necked opening for easy filling. Made of transparent polycarbonate with white or blue PP lid. Window-clear. High impact resistance. Non-toxic. Stackable. Popular screw-top jar also for freezing down to -90°C. Fully unscrew the lid prior to autoclaving.

Cat. No.	Content	Pack	Ø Lid	Height
LPB21160125	125 ml	4 pcs.	70 mm	74 mm
LPB21160250	250 ml	4 pcs.	70 mm	119 mm
LPB21160500	500 ml	4 pcs.	120 mm	86 mm
LPB21161000	1000 ml	4 pcs.	120 mm	142 mm

Thorbi insulated container + Blood Bag and Plasma Storage Racks



Thorbi insulated container without lid

Suitable for ice, dry ice and liquid nitrogen. One-piece foamed PVC, thick-walled with very good insulating properties. Unbreakable during normal use. Chemically resistant to acids, alkalis, salts, alcohol, acetone, etc. Temperature-resistant from -196°C to +93°C.

Content: 1 litre
WxDxH: 200x175x91 mm

Cat. No.	Colour
LCB110BU	blue
LCB110GR	green
LCB110RE	red
LCB110YE	yellow
LCB110BL	black

Content: 4 litres
WxDxH: 350x240x122 mm

Cat. No.	Colour
LCB140BU	blue
LCB140GR	green
LCB140RE	red
LCB140YE	yellow
LCB140BL	black

Content: 9 litres
WxDxH: 490x370x140 mm

Cat. No.	Colour
LCB190BU	blue
LCB190GR	green
LCB190RE	red
LCB190YE	yellow
LCB190BL	black



Thorbi insulated container with lid

Suitable for ice, dry ice and liquid nitrogen. One-piece foamed PVC, thick-walled with very good insulating properties. Unbreakable during normal use. Chemically resistant to acids, alkalis, salts, alcohol, acetone, etc. Temperature-resistant from -196°C to +93°C.

Content: 2,5 litres
WxDxH: 330x280x122 mm

Cat. No.	Colour
LCB225BU	blue
LCB225GR	green
LCB225RE	red
LCB225YE	yellow
LCB225BL	black

Content: 4,5 litres
WxDxH: 330x280x185 mm

Cat. No.	Colour
LCB245BU	blue
LCB245GR	green
LCB245RE	red
LCB245YE	yellow
LCB245BL	black



Blood Bag and Plasma Storage Racks

Blood bag and plasma storage racks made of tough 3-mm acrylic glass, transparent. Handles at the two ends. With 2 movable partitions in the standard unit. Labelling tabs and cards optional. Card size WxH = 25x50 mm. Modified designs possible.

Cat. No.	Inside WxDxH (mm)	Outside WxDxH (mm)
LPB900-450104	134x444x100 ¹⁾	140x450x104 ¹⁾
LPB900-500104	134x494x100 ¹⁾	140x500x104 ¹⁾
LPB900-550104	134x544x100 ¹⁾	140x550x104 ¹⁾
LPB900-600104	134x594x100 ¹⁾	140x600x104 ¹⁾
LPB900-650104	134x644x100 ¹⁾	140x650x104 ¹⁾

¹⁾ Different dimensions possible.



Cat. No.	Accessories
LPH800-100130G	Partition H = 100, W = 130
LPH901	Clip-on labelling tab
LPH910PW	Labelling cards, plastic, white, 50 pieces
LPH910PY	Labelling cards, plastic, yellow, 50 pieces
LPH910KW	Labelling cards, cardboard, white, 35 pieces/sheets

