

CERTIFICATE OF CALIBRATION

09-400-5061	Impföse/Nadel PS, 1µl, fest, Ø4x200 mm,
09-401-5001	Impföse/Nadel PS, 1µl, flexibel, Ø4x200 mm, 20
09-420-5061	Impföse/Nadel PS, 10 µl, fest, Ø4x200 mm,
09-421-5065	Impföse/Nadel PS, 10 µl, flexibel, Ø4x200 mm,

INTENDED USE

Calibrated Disposable Inoculating Loops are suitable for use in semi-quantitative inoculation procedures, such as Urine counts, Sampling, Serial dilution and Bacteriological Inoculation.

CALIBRATION STANDARD

The delivery volume of each batch of calibrated Disposable Inoculating Loops are tested per the "EVANS BLUE DYE COLORIMETRIC PROCEDURE" as described in the "Clinical Microbiology Procedures Handbook" (American Society for Microbiology, Washington, DC.).

Inaccuracy % : Acceptable Tolerance +/- 20%

Nerbe plus GmbH Co. KG

14.11.2022


nerbe plus GmbH & Co. KG
Porschestraße 25 · 21423 Winsen/Luhe
Germany
Sylvia Walter
Quality Management