



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	174CS05
GENERAL DESCRIPTION	SPREADERS BUSTA 5PCS 1000PKG
DETAILED DESCRIPTION	1 box containing 200 paper-film pouches each containing 5 L-shape spreaders

PRODUCT COMPONENTS

SPREADERS	SPREDERS A L	DIMENSIONAL SPECIFICATIONS
		TOTAL LENGTH: 145.0 mm \pm 1.0 mm DIAMETER: 4.0 mm \pm 0.2 mm LENGTH "L": 38.0 mm \pm 0.5 mm : 7.0 mm \pm 0.5 mm

PACKAGING SPECIFICATION

PACKAGING	PAPER-FILM POUCH
BOX QTY	1000

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS