PCR Tubes

0.2 ml | 0.5 ml flat and domed cap













- Effortless and secure working
 Ergonomic cap design allows
 opening and closing with one hand
- Elimination of contamination and loss of sample
 High quality raw material meets standards of purity and is highly resistant to extreme mechanical strain and has excellent chemical resistance

Centrifuge specifications: up to 25.000 x g - fixed-angle rotor up to 70.000 x g - swing out rotor

- Maintain sample integrity sterile version available (by beta irradiation)
- Optical measurement and unsurpassed clarity
- Fully autoclavable tube

PCR Tube		Tube Height	Cat. No.	Pack Type	Sterilized	Sales Unit
0.2 ml, flat cap, clear PCR Tube	T	21.5 mm	3-010-50-0 3-010-80-0 3-011-50-0	Bag Bag Bag		1 x 1000 10 x 1000 1 x 1000
0.2 ml, domed cap, clear PCR Tube	IV	23.3 mm	3-020-50-0 3-020-80-0 3-021-50-0	Bag Bag Bag		1 x 1000 10 x 1000 1 x 1000
0.5 ml, flat cap, clear PCR Tube	il	32.1 mm	3-110-50-0 3-110-80-0 3-111-50-0	Bag Bag Bag		1 x 1000 25 x 1000 1 x 1000
0.5 ml, domed cap, clear PCR Tube		34.3 mm	3-120-50-0 3-120-80-0 3-121-50-0	Bag Bag Bag		1 x 1000 25 x 1000 1 x 1000