



## **TECHNICAL DATASHEET**

**PART NUMBER:** 

131054C



## **ClearLine® Cell Scrapers**

**Description:** Sterile cell scrapers, two-position blade

## Features:

- ✓ Two-position blade : cell scraper or cell lifter
- ✓ Handle in ABS, blade in TPE
- ✓ Individually wrapped
- ✓ Lot n° on each bag, certificates are available on demand
- √ Sterilised by gamma irradiation
- ✓ Endotoxins certified < 0,5 EU /ml
  </p>

Cat No	Total length (mm)	Blade width (mm)	Units /bag	Units /box
131054C	250	20	1	100