



TECHNICAL DATASHEET

PART NUMBER:

131056C



ClearLine® Cell Scrapers

Description: Sterile cell scrapers, two-position blade

Features:

- ✓ Two-position blade : cell scraper or cell lifter
- ✓ Handle in ABS, blade in TPE
- ✓ Individually wrapped
- ✓ Lot n° on each bag, certificates are available on demand
- √ Sterilised by gamma irradiation
- ✓ Endotoxins certified < 0,5 EU /ml
 </p>

Cat No	Total length (mm)	Blade width (mm)	Units /bag	Units /box
131056C	390	30	1	100