

Product Name: Capinsert™

Catalogue No.: T312_Series

This document replaces any previous version

1. Product Description:

- **Color Capinsert™:** Disposable, non-sterile for color coding of Cryovial® and sample tube cap
- **Compatible with:**
 - T301, T308, T309, T310 and T311 Cryovial® Series (Except T310-10A)
 - T500, T501 and T502 Screw Cap Series

2. Packaging:

- **Case:** 1 bag of 500 units

3. Product Specifications:

- **Material:** Polypropylene
- **Writing surface** on top
- **Temperature range:** -196°C to +121°C
- **Autoclavable** at 121°C, for up to 30 minutes

4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **CONEG / RoHs:** Plastics and colorants are in conformity with CONEG /RoHS standards for heavy meta
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **REACH (SVHC):** Plastic is in conformity to REACH standards
- **LATEX:** Material is Latex Free
- **BSE / TSE:** Material is BSE / TSE Free

5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimensional measurements
- Insert verification

6. Traceability:

- **Lot No. Composition:** 8 or 9 digits

- **The lot number can be found in one or all of these locations:**
 1. On exterior case label
 2. On label inserted inside the master case
 3. On the inner bag

7. Storage Conditions:

- Store at room temperature in normal warehouse conditions
- Avoid temperature variations and humidity
- Protect from any possible contamination
- Protect from any damage to the packaging

8. Recommended Use:

- None

OWNERSHIP OF MATERIALS:

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*