

# Product Specifications

**CL02.0504**

## **Ethanol 70 vol. % in Water solution**

Etanol 70 vol. % solución en Agua  
Alcole etílico 70 vol. % in soluzione Acqua  
Ethanol 70 vol. % in Wasser Lösung  
Ethanol 70 vol. % in Water oplossing  
Ethanol 70 vol. % dans l'Eau solution

**\*For laboratory use, Belg. Pharm. V**  
**700ml C<sub>2</sub>H<sub>5</sub>OH + 300ml H<sub>2</sub>O + 3% (C<sub>2</sub>H<sub>5</sub>)<sub>2</sub>O**

Mol weight: 46.07 g/mol  
Density: 0.85 g/ml  
HS Nr: 22072000  
CAS Nr: 64-17-5  
UN: 1993  
ADR: 3,II



### **Specifications:**

### **NOVOLAB nv**

Diebeke 7 | t 0032(0)54 42.15.80  
Industriezone | f 0032(0)54 42.17.25  
9500 Geraardsbergen | e info@novolab.be  
Belgium | w [www.novolab.be](http://www.novolab.be)

Antigifcentrum (BE) - Centre Antipoisons (BE)  
070 245 245

Nationaal Vergiftigingen Informatie Centrum(NL):  
088 755 8000