

Labcon, North America 3700 Lakeville Highway Suite 200 Petaluma, CA 94954-5671 U.S.A. 707.766.2100 Fax 707.766.2199 www.labcon.com

Product Description

Catalog Number:

3184-345-001-9

Product Description:

50mL Metal-Free Centrifuge Tube Polypropylene. Sterile. Graduations. Flat natural-colored polyethylene cap

Product Dimensions:

Length of Tube: Outer Diameter: Inner Diameter: 4.51" (114.55mm) 1.16" (29.46mm) 1.08" (27.43mm)

Packaging:50 tubes per bag (with attached caps)10 bags per case (500 tubes total)

Autoclavable. Freezable. Tubes can be centrifuged to 20,000 RCF.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNAse and DNAse), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Sterility:

Certified as sterile with a sterility assurance level (SAL) of 10⁻⁶