




designed for scientists



R 3004 Blade stirrer

/// Data Sheet

Generates a radial flow. For drawing the material to be mixed from the top and bottom. For aerating of low viscosity media and for homogenizing solids in liquids. Used at medium to high speeds.



designed for scientists

Technical Data

Stirrer diameter [mm]	30
Shaft diameter [mm]	8
Shaft length [mm]	359
Material in contact with medium	AISI 316L / 1.4571
Weight [kg]	0.156
Permissible ambient temperature [°C]	5 - 40