

**Reduce your plastic consumption and your carbon footprint by 80% thanks to HybridRack!**

HybridRack consists of a rack and a cover made of 100% recyclable cardboard

The base of the rack is reinforced to guarantee the same stability and pipetting comfort as traditional plastic racks



**GoGREEN**

Per box, with HybridRack, the reduction in plastic amounts to 3 kg for 50 racks, which also represents a saving of 1,801 litres of water\*.

In fact, between 77% and 86% less plastic is used in the manufacture of a HybridRack compared with a conventional rack.


Still fitted with the same coloured polypropylene trays, the quality of the HybridRacks tips remains unchanged!

Available with sterile filter or without sterile filter




\* Complete box of 40 racks of 300 µl tips

Plastic weights  
of Hybridrack and conventional racks



### What's more:

- HybridRacks are lighter, helping to reduce the carbon footprint associated with transport
- The manufacture of this product benefits from a third of renewable energy and a reduction in the use of water and plastic
- The packaging is made from 100% recycled material



### **FILTER TIPS**

The filter does not trap the sample so, in the event of over-pipetting, some of the liquid can be recovered

**Low-binding** tips

### **NON-FILTER TIPS**

Packaged in **sterile rack**



Sterile rack:  
**Standard** tip  
Reusable racks with MultiRack refills



## Tips, 10 $\mu$ l

### Long reach



- Available in filtered and unfiltered versions
- Long and graduated at 2, 5, 10  $\mu$ l
- Ideal for reaching the bottom of 1.5 and 2 ml microtubes without contaminating the micropipette
- Part number: **087240T** and **072120**



## Tips in 20, 100 and 200 $\mu$ l

### Z-tips



- Available in 20 and 100  $\mu$ l with filter and 200  $\mu$ l without filter
- Narrow neck requiring low tip insertion and ejection force
- High pipetting accuracy
- Scaled to 10, 50 and 100  $\mu$ l
- Part number: **087250T**, **087260T**, **072130**



## Tips in 200 and 300 $\mu$ l



### NX tips

- Available in 200  $\mu$ l filtered or 300  $\mu$ l unfiltered versions
- Fine point for precise dispensing, extended length to reach the bottom of tubes
- Scaled to 10, 50, 100  $\mu$ l
- Part number: **087270T** and **072140**



### 300 $\mu$ l XT tips



- Available in filter version
- Flexible neck fits most 300  $\mu$ l and 200  $\mu$ l micropipettes and ejects easily
- Fine tip for precise dispensing
- Scaled to 10, 50, 100 and 200  $\mu$ l
- Part number: **087280T** and **072150**



## Tips, 1250 $\mu$ l

- Available in filtered and unfiltered versions
- Long tip with exclusive design: ideal for 10 and 50 ml tubes, assay blocks, etc
- Scaled to 200, 500 and 1000  $\mu$ l
- Part number: **087290T** and **072160**

### 1250 $\mu$ l VS Std 1000 $\mu$ l



## NOVOLAB nv

Diebeke 7 | t 0032(0)54 42.15.80  
 Industriezone | f 0032(0)54 42.17.25  
 9500 Geraardsbergen | e info@novolab.be  
 Belgium | w [www.novolab.be](http://www.novolab.be)