Document Number: **BD-T0048**

Document Date: January 2011

Page 1/2

Cell Scrapers & Lifters

I. APPLICATIONS

- 1. Cell culture
- 2. Cutting and scraping of surfaces

II. FEATURES

- 1. Materials: polyethylene
- 2. Gamma radiation sterilized
- 3. Non-pyrogenic

III. SCRAPER/LIFTER DESIGN

- 1. Scrapers available in 18cm and 25cm, lifter is 18cm
- 2. 1.8cm blade length for scraper, 2.0cm blade length for lifter
- 3. Pivoting blade on scraper for precise cutting and scraping
- 4. Individually wrapped

IV. PACKAGING

1. Bulk, 100 pieces/case or 200 pieces/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
70-1180	Cell scraper. Handle length: 18cm, blade	Individually wrapped, 100 pieces/case	Yes
	width: 16.68	individually wrapped, 100 pieces/case	
70-1250	Cell scraper. Handle length: 24.3cm,	Individually wrapped 200 pieces/case	Yes
	blade width: 16.68	Individually wrapped, 200 pieces/case	
70-2180	Cell lifter. Handle length: 18cm, blade	Individually wrapped, 100 pieces/cook	Yes
	width: 18.74	Individually wrapped, 100 pieces/case	

VI. STERILITY

- 1. Products are Gamma Radiation Sterilized according to ISO13485 standard.
- 2. For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2015 standard.

VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PHOTOS

Tel: +86-531-67802668 Fax: +86-531-67803768 www.BiologixGroup.com

Document Number: **BD-T0048**

Document Date: January 2011

Page 2**/2**






Scrapers

Lifters

X. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	Α	В	С	D	E
70-1250	Cell Scraper	240.4	16.68	9	89.0	O.D. 7.7
70-1180	Cell Scraper	180.0	16.68	9	85.0	O.D. 7.4
70-2180	Cell Lifter	180.0	18.74	7.4	-	-



Tel: +86-531-67802668 Fax: +86-531-67803768 www.BiologixGroup.com