

VP-8/5C, block

DESCRIPTION

8 × 5 ml conical tubes

SPECIFICATIONS

Max. RPM	1200
Max. temperature	100°C



CAT. NUMBER

BS-010176-SK

VP-8/5C

ACCESSORIES



TS-100C

Without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples down ...

[read more](#)



TS-100C Smart

Software included, without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C Smart** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples ...

[read more](#)