

Technical Datasheet

Date: 08.04.21

Page 1 of 1

EN ISO 13485 certified		
Article - No:	09-812-8071	Revision: 1
Remarks		single use only

Specification

Product	np GREEN SWAB flocked swab, very fine nylon head, ABS stick, breakpoint at 80mm, sterile R, CE2797
Classification	CE2797
Dimensions	Length: 150 mm (+/- 3mm)
	Width: 2,5 mm (+/- 0,2mm)
	Diameter head: 3 mm (+/- 1mm)
Material Characteristics	head: nylon
	stick: Acrylonitrile butadiene styren (ABS)
Product Characteristics	biologically inert
	physiologically harmless
	free of detectable phthalate acc. to REACH regulation
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per pallet: 129.600
	pcs per inner cs: 100
	pcs per carton: 1.200
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	